



Dec 01, 2025

Cryopreservation of Cells

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g7n2qqllwz/v1>

Mohit Sharma¹

¹Jiwaji University



Mohit Sharma

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7n2qqllwz/v1>

Protocol Citation: Mohit Sharma 2025. Cryopreservation of Cells. [protocols.io](https://dx.doi.org/10.17504/protocols.io.q26g7n2qqllwz/v1)

<https://dx.doi.org/10.17504/protocols.io.q26g7n2qqllwz/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 01, 2025



Last Modified: December 01, 2025

Protocol Integer ID: 233855

Keywords: Cell Culture, Cryopreservation , Animal Cell culture, cells trypsinization, adherent cells from the vessel, cryopreservation of cell, cryopreservation, animal cell, systematic approach for the preservation, cell, maintaining cellular integrity, preservation, subsequent thawing, functionality during storage, cellular integrity, storage

Abstract

This protocol outlines a systematic approach for the preservation of animal cells (especially Vero E6), aimed at maintaining cellular integrity, viability, and functionality during storage and subsequent thawing. The primary goal is to optimize long-term storage methods to support cell-based research and clinical applications

Guidelines

Store the frozen cells below -70°C ; frozen cells begin to deteriorate above -50°C .

Materials

- Culture vessels containing cultured cells in log-phase of growth
- Complete growth medium, pre-warmed to 37°C
- Balanced salt solution such as Phosphate Buffered Saline (PBS), containing no calcium, magnesium, or phenol red
- Dissociation reagent such as trypsin
- Cryoprotective agent such as DMSO
- Disposable, sterile 15-mL or 50-mL conical tubes
- Equipment to determine viable and total cell counts such as hemacytometer or Coulter Counter
- Sterile cryogenic storage vials (cryovials)
- Controlled rate freezing apparatus or isopropanol chamber
- Liquid nitrogen storage container
- **Freezing medium-** 20% FBS + 10% DMSO in cell culture media

Safety warnings

- ! DMSO at concentrations higher than 10% (v/v) is reported to induce membrane pore formation and apoptosis, thus Time is thus a critical factor when you freeze cells to minimize their exposure to DMSO.

Before start

All solutions and equipment that come in contact with the cells must be sterile. Always use proper sterile technique and work in a laminar flow hood.



- 1 Prepare freezing medium and store at 2° to 8°C until use
- 2 Remove the media from the Flask
- 3 Treat the cells with Trypsin/Versene solution to detach the cells from the substratum.
- 4 Determine the total number of cells and percent viability using a hemacytometer or cell counter
- 5 Resuspend the cell pellet in cold freezing medium at the recommended viable cell density for the specific cell type.
- 6 Dispense aliquots of the cell suspension into cryogenic storage vials.
- 7 Freeze the cells in a controlled rate freezing apparatus, decreasing the temperature approximately 1°C per minute. Alternatively, place the cryovials containing the cells in an isopropanol chamber and store them at -80°C overnight.
- 8 Transfer frozen cells to liquid nitrogen and store them