


Aug 03, 2020

Cryopreservation of organoid cultures

DOI

dx.doi.org/10.17504/protocols.io.bh4ij8ue



Emily Souster¹, Hazel Rogers¹, Laura Letchford¹, Sara Vieira¹, Maria Garcia-Casado¹, Mya Fekry-Troll¹, Charlotte Beaver¹, Rachel Nelson¹, Hayley Francies¹, Mathew Garnett¹

¹Wellcome Sanger Institute

Cellular Generation and ...

Organoid and Assembloid



Emily Souster

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bh4ij8ue>

Protocol Citation: Emily Souster, Hazel Rogers, Laura Letchford, Sara Vieira, Maria Garcia-Casado, Mya Fekry-Troll, Charlotte Beaver, Rachel Nelson, Hayley Francies, Mathew Garnett 2020. Cryopreservation of organoid cultures. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bh4ij8ue>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 30, 2020

Last Modified: August 03, 2020



Protocol Integer ID: 38762

Keywords: cryopreservation of organoid culture, organoid culture, organoid derivation team within the cellular generation, organoid derivation team, organoid, organoid model, cancer organoid, expanding organoid model, procedure for cryopreservation, cryopreservation, cellular generation

Abstract

This protocol defines the procedure for cryopreservation of organoid cultures. It has been developed by the organoid derivation team within the Cellular Generation and Phenotyping Group at the Wellcome Sanger Institute. The team has extensive experience passaging and expanding organoid models. The method described has mainly been used for cancer organoids derived from colon, pancreas and oesophageal tumours.

Guidelines

- We use 5 ml Eppendorf tubes to help with sterility. However, if you do not have access to these tubes any alternative sterile tubes of appropriate volume can be used.
- It is useful to keep a bottle of cold PBS in the fridge as this can be used to resolve pelleting issues. Resuspending a 'hazy' pellet in cold PBS can help to re-melt the BME2, resulting in a more distinct cell pellet after centrifugation.
- Appropriate freezing containers will ensure that the liquid freezes at a controlled rate of around  -1 °C per minute at  -80 °C .



Materials

MATERIALS

✕ DPBS, no calcium, no magnesium **Thermo Fisher Catalog #14190136**

✕ Eppendorf tube- 5ml **Eppendorf Catalog #0030122321**

✕ Recovery Cell Freezing Media **Gibco - Thermo Fisher Scientific Catalog #12648-010**

✕ Nunc™ Biobanking and Cell Culture Cryogenic Tubes **Thermo Fisher Scientific Catalog #375418**



✕ Nunc™ Cell Scrapers **Thermo Fisher Scientific Catalog #179693**

Equipment

- Sterile cell culture hood
- Centrifuge
- P1000, P200, tips
- Cell counting equipment
- CoolCell or appropriate freezing container
- -80°C freezer
- Liquid Nitrogen storage

Troubleshooting

Before start

- Before cryopreservation of the line, ensure the organoids are of good quality and avoid banking organoids that have become very large or dark.
- Place a labelled CoolCell or appropriate freezing container in the fridge and make sure it is at  4 °C before commencing work.
- Thaw an appropriate amount of freezing media (1ml per cryovial) and keep at  4 °C before use.


Note

The freezing media can be kept at  4 °C for a maximum of 7 days.

Process diagram

1

Banking of organoid cell cultures:



Protocol

- Using a P1000 pipette, or a cell-scraper, detach BME2 drops from as many wells as required (generally 1 well into 1 cryovial- we recommend aiming for a density of 1×10^6 cells per cryovial).

Note


If you plan to bank a specific number of cells, be sure to take a cell count at Step 7.

The count can be taken from the cell suspension, or alternatively you could detach a representative well at this stage, from which to take the count.

- Transfer the organoid suspension to a 5ml Eppendorf or 15ml Falcon tube.
- Wash the wells with cold PBS (~ 4 °C) to collect any remaining organoids and add to the collection tube.



Make sure organoids are dissociated from BME2 drops by aspirating up and down multiple times with a stripette.

- 5 Centrifuge at  800 x g for 2 minutes (with the brake set on 4).


Safety information

Centrifuge buckets must be sealed using safety caps, which must only be opened inside a microbiological safety cabinet.

- 6 Aspirate off the supernatant, leaving the organoid cell pellet at the bottom of the tube. Organoid cell pellets are loose pellets so take care when aspirating.

Note

If the pellet is not clear, aspirate as close to the BME2 'haze' as possible and re-suspend in cold PBS, spin down again (800 x g for 2 minutes, brake set on 4) and then aspirate the supernatant.



BME2 layer above organoid pellet after centrifugation.

- 7 If you **do not** want to bank a specific number of cells, proceed to Step 8.

If banking a specific number of cells, either take a cell count from your whole organoid suspension, or from a 'representative well' by resuspending the pellet in an appropriate amount of organoid media. Take a cell count from this suspension.

From this cell count you can calculate how much organoid suspension to use to freeze the desired cell density.

Note


Organoids are typically counted as aggregates because they are difficult to break down to single cells.

We take two aggregate cell counts which must be within 10% of each other, otherwise a third count is taken and an average is calculated between the three counts.

We recommend freezing a cell density of 1×10^6 cells per cryovial.



8 Resuspend the organoid pellet in as much Recovery Cell Freezing Media as required (1ml per cryovial) and transfer into the prelabelled cryovials.

9 Place the cryovials into a CoolCell or appropriate freezing container and place into  -80 °C freezer.

Note

If there is an odd number of cryovials, use a 'balance' cryovial for better cryopreservation.

10 After 24 hours (up to 72 hours) transfer cryovials to liquid nitrogen storage.