



Apr 25, 2025

Version 2

🌐 Crystallisation of Enterovirus D68 3C protease apo structure and fragment screen in space group $P2_12_12_1$ (PDB code 8CNX) V.2

DOI


<https://dx.doi.org/10.17504/protocols.io.n92ldm6jnl5b/v2>Ryan Lithgo^{1,2,3}, Peter Marples^{1,2,3}, Lizbé Koekemoer^{4,5,3}, Daren Fearon^{1,3}¹XChem, Diamond Light Source; ²Research complex at Harwell; ³ASAP Discovery Consortium;⁴Centre of Medicines Discovery; ⁵University of Oxford

Ryan Lithgo: The principle crystallographer on the Enterovirus 3C protease project.;

ASAP Discovery

**Lizbé Koekemoer**

University of Oxford



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)OPEN  ACCESS**DOI:** <https://dx.doi.org/10.17504/protocols.io.n92ldm6jnl5b/v2>**External link:** <https://asapdiscovery.org/outputs/target-enabling-packages/#ASAP-EVD68-3CLPRO>

Protocol Citation: Ryan Lithgo, Peter Marples, Lizbé Koekemoer, Daren Fearon 2025. Crystallisation of Enterovirus D68 3C protease apo structure and fragment screen in space group P212121 (PDB code 8CNX). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.n92ldm6jnl5b/v2> Version created by **Mary-Ann Xavier**

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 08, 2025

Last Modified: April 25, 2025

Protocol Integer ID: 119833

Keywords: crystallisation, 3C protease, 3CNY, XChem, ASAP, AViDD, Diamond Light Source, Enterovirus D68, crystallisation of enterovirus d68 3c protease, enterovirus d68 3c protease crystal, enterovirus d68 3c protease, enterovirus d68, crystallization protocol, fragment screen in space group p212121, space group p212121, crystal structure, rcsb pdb with code 8cnx, crystal form, pdb code 8cnx, crystal form suitable for fragment screen, crystallisation, rcsb pdb, protease, protein expression, purification protocol, protein

Funders Acknowledgements:

National Institutes of Health/National Institute Of Allergy and Infectious Diseases (NIH/NIAID)

Grant ID: Grant ID: U19AI171399

Abstract

The crystallization protocol and buffer conditions used to obtain Enterovirus D68 3C protease crystals. The crystal structure was deposited to RCSB PDB with code **8CNX**.

In this new version we added the Addgene ID and the protein expression and purification protocol. This construct gave us a crystal form suitable for fragment screen and initial compound soaking.

Materials

 EV-D68 3C protease; aliases: Enterovirus A71 3C protease, non cleavable C-term 6xHis tag **addgene Catalog #**(Plasmid #204817)

SwissCI 3 lens crystallization plates <https://swissci.com/product/3-lens-crystallisation-plate/>

1 M Tris-HCl pH 8.14

1 M Ammonium acetate

50% w/v PEG 3350

Protocol materials

 EV-D68 3C protease; aliases: Enterovirus A71 3C protease, non cleavable C-term 6xHis tag **addgene Catalog #**(Plasmid #204817)



Protein expression and purification

2d

1 Small scale protein purification and expression protocols used for this assay



EV-D68 3C protease; aliases: Enterovirus A71 3C protease, non cleavable C-term
6xHis tag **addgene Catalog #**(Plasmid #204817)

Protocol

NAME

Enterovirus D68 3C protease small scale expression and
purification protocol

CREATED BY

korvus.wang

[Preview](#)

Equipment needed

2 Formulatrix Rock Imager (or incubator of choice)

SPT mosquito

P100 multi-channel pipette

SwissCI 3 lens plate

Crystallisation experiment

1d

3 Protein and buffer requirements:

32 μ L

[M] 1 millimolar (mM)



Sample



2.88 mL Crystallization screen

16 μ L

Sample

seeds, dilution 1:100 000

4 Crystallisation screen composition:

0.1 M Tris-HCl (pH 8.14)



0.2 M Ammonium acetate
25% w/v PEG 3350

Stock solutions used:

1 M Tris-HCl adjusted to pH 8.14 with HCl

1 M Ammonium acetate

50% w/v PEG 3350

Note

The crystallisation screen can be stored in a duran bottle or aliquoted into 96 deep well block for easy dispensing into SwissCI 3 lens plates.

For long term storage keep the Crystallisation screen in the fridge at 4°C.

5 Prepare seed stock:

1: 100 000 dilution  Sample seeds

6 Dispense  30 µL Crystallisation screen into SwissCI 3 lens plate reservoir wells using a 100 µl multi-channel pipette.

Dispense  100 µL  1 millimolar (mM)  Sample to each lens using the SPT mosquito.

Dispense  200 µL Crystallisation screen to each lens using the SPT mosquito.

Dispense  50 µL Seeds to each lens using the SPT mosquito.

Drop ratio: 1:2:05 ratio (100 nl  Sample : 200 nl reservoir solution: 50 nl seeds)

Final drop volume: 350 nl

7 Incubate at  20 °C for  24:00:00 h in Formulatrix Rock Imager.

1d

Imaging Schedule: The first images are taken after 12 h and the imaging schedule follows a Fibonacci sequence of days for further collections.

8 Crystal form after ~12 h.

Expected result

The crystals reach their maximum size after 24 h and most of the precipitant has gone.

Crystals typically form either as single crystals or in small clusters containing 4-6 crystals in a flower-like arrangement. Upon mounting a gentle nudge dissociates the cluster.

Morphology: typically thin rectangles with pointed ends.


Size: ~100 μm in length and ~50 μm in width, depth of the crystals is ~10 μm

Appearance: glass shard.

Average resolution: 1.5 $\text{\AA}$

Space group: $P2_12_12_1$

Unit cell: 42.82 $\text{\AA}$, 62.53 $\text{\AA}$, 147.36 $\text{\AA}$
90.00°, 90.00°, 90.00°



An example of a drop containing EV D68 3C protease crystals.

Soaking tolerance

2h 30m



9 Final condition identified: 20% DMSO, 2hr 30mins using the protocol bellow.

2h 30m

Protocol



NAME

XChem solvent test

CREATED BY

Peter Marples

Preview

Fragment and compound soaking

2h 30m

10 20% fragment/compound, 2hr 30mins using the protocol bellow.

2h 30m

Protocol



NAME

XChem crystallographic fragment screening

CREATED BY

Peter Marples

Preview

Protocol references

N/A