

Oct 16, 2025

Version 2

Crystallisation of SARS-CoV-2 C terminal nucleocapsid protein for fragment screening and compound soaking V.2

Forked from a private protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.8epv5r8ydg1b/v2>

Peter Marples^{1,2,3}, Lizbé Koekemoer^{4,3}, Daren Fearon^{1,2,3}

¹Diamond Light Source; ²Research Complex at Harwel; ³ASAP Discovery Consortium;

⁴Centre of Medicines Discovery, University of Oxford

ASAP Discovery



Peter Marples

Diamond Light Source



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5r8ydg1b/v2>

External link: <https://asapdiscovery.org/outputs/target-enabling-packages/#ASAP-SARS-COV-2-NPROTEIN>

Protocol Citation: Peter Marples, Lizbé Koekemoer, Daren Fearon 2025. Crystallisation of SARS-CoV-2 C terminal nucleocapsid protein for fragment screening and compound soaking. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.8epv5r8ydg1b/v2> Version created by **Mary-Ann Xavier**

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 20, 2025

Last Modified: October 16, 2025

Protocol Integer ID: 225096

Keywords: crystallisation, XChem, ASAP, AVIDD, CMD, Diamond Light Source, i04-1, SARS CoV-2 Nucleocapsid, Nucleocapsid, N-protein, Research complex at Harwell, crystallisation of sar, crystals suitable for xchem fragment screening, reproducible sar, xchem fragment screening, crystallization protocol, crystallisation, sar, crystal, fragment screening, fragment soaking condition

Funders Acknowledgements:

National Institutes of Health/National Institute Of Allergy and Infectious Diseases (NIH/NIAID)

Grant ID: Grant ID: U19AI171399

Disclaimer

The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.

Acknowledgements:

Diamond Light Source Ltd, Harwell Science and Innovation Campus, Didcot OX11 0QX, UK
Research Complex at Harwell, Harwell Science and Innovation Campus, Didcot OX11 0FA, UK
Oxford Lab Technologies crystal shifter <https://doi.org/10.1107/S2059798320014114>

Abstract

The crystallization protocol and buffer conditions used to obtain reproducible SARS CoV-2 Nucleocapsid crystals suitable for **XChem** fragment screening. In this new version we added the Addgene id, together with the expression and purification protocol. We also added the solvent tolerance test and the compound/fragment soaking conditions.



Materials

SwissCI 3 lens crystallization plates <https://swissci.com/product/3-lens-crystallisation-plate/> Codes: *Midi*: UVXPO-3LENS 3W96T-PS 3W96T-UVP

[M] 1 Molarity (M) HEPES adjusted to pH 7.8 with NaOH, Molecular Dimensions, Catalog # MD2-011-PH 7.8
70 % isopropanol, Sigma Aldrich, Catalog # **563935**
50% w/v PEG 4000, Molecular Dimensions, Catalog # MD2-250-11

Purified SARS CoV-2 Nucleocapsid protein ([M] 20 mg/mL) in [M] 10 millimolar (mM) HEPES, pH 7.5 ,
[M] 0.5 Molarity (M) NaCl, 5% glycerol, [M] 0.5 millimolar (mM) TCEP

Construct used protein residues 250-364

Protocol materials

⊗ SARS-CoV-2 isolate : Wuhan-Hu-1 Nucleocapsid Protein C terminal **addgene** Catalog #228641

Safety warnings

⚠ Follow all handling warning for the chemicals used in the crystallisation screen composition.



C-terminal protein expression and purification

- 1 For protein heterolog expression we used the plasmid



SARS-CoV-2 isolate : Wuhan-Hu-1 Nucleocapsid Protein C
terminal **addgene** Catalog #228641

and the following protocol.

Equipment needed

- 2 **Formulatrix Rock Imager** (or incubator of choice)
SPT mosquito
P100 8 multi-channel pipette
SwissCI 3 lens plate

Crystallisation experiment

1d

- 3 **Prepare seed stock:**

Protocol



NAME

Diamond XChem Seeding Protocol

CREATED BY

Peter Marples

Preview

1: 100 dilution Sample seeds

- 4 **Protein and buffer requirements:**

57.6 μ L

20 mg/mL



Sample



2.88 mL Crystallization screen

5.76 μ L

Sample

seeds, dilution 1:100

- 5 **Crystallisation screen composition:**

0.1 Molarity (M)

HEPES



7.8

10 % isopropanol

23% w/v PEG 4000

**Stock solutions used:**

[M] 1 Molarity (M) HEPES adjusted to  7.8 with NaOH

70 % isopropanol

50% w/v PEG 4000

Note

The crystallisation screen can be stored in a duran bottle or aliquoted into 96 deep well block for easy dispensing into SwissCI 3 lens plates.

For long term storage keep the Crystallisation screen in the fridge at 4°C.

- 6 Dispense  30 µL Crystallisation screen into SwissCI 3 lens plate reservoir wells using a 100 µl multi-channel pipette.

Dispense  200 µL [M] 20 mg/mL  Sample to each lens using the SPT mosquito.

Dispense  180 µL Crystallisation screen to each lens using the SPT mosquito.

Dispense  20 µL Seeds to each lens using the SPT mosquito.

Drop ratio: 5:9:1 ratio (100 nL  Sample : 180 nL reservoir solution: 20 nL seeds)

Final drop volume: 300 nL

- 7 Incubate at  20 °C for  24:00:00 h in Formulatrix Rock Imager.

1d

Imaging Schedule: The first images are taken after 12 h and the imaging schedule follows a Fibonacci sequence of days for further collections.

- 8 Crystal form after ~12 h.

Expected result

The crystals reach their maximum size after 24 h.

Crystals typically form as single crystals as six sided shards

Morphology: six sided shard


Size: ~100 μm in length and ~30 μm in width, depth of the crystals is ~30 μm

Appearance: glass shard.

Average resolution: 2.0 $\text{\AA}$

Space group: $I4_1$

Unit cell: 89, 89, 40
90.00°, 90.00°, 90.00°



An example of a drop containing SARS N-protein crystals.



Data collection at Synchrotron

- 9 Diamond Light Source
Unattended Data Collection (UDC)
Data Collection Temperature: 100K
Detector: DECTRIS EIGER2 X 9M
Beamline: I04-1
Wavelength: 0.9212 Å
Resolution (Å): 1.68
Beam Size (μm): 60 X 50
Number of images: 3600
Oscillation: 0.10°
Exposure (s): 0.0020
Transmission (%): 100
Flux (ph/s): 9.50e+11

Soaking tolerance test

- 10 Crystals can tolerate up to 20% DMSO, for up to  03:00:00

Protocol



NAME

XChem solvent test

CREATED BY

Peter Marples

[Preview](#)

Compound soaking conditions

- 11 Crystals can tolerate up to 20% DMSO, for up to  02:00:00



Protocol



NAME

XChem crystallographic fragment screening

CREATED BY

Peter Marples

Preview

Protocol references

N/A