

May 13, 2025

Version 2

🌐 Crystallisation of SARS-CoV-2 MPro in P2₁2₁2₁ for compound soaking V.2

🔗 Version 1 is forked from [Crystallization of SARS-CoV-2 Mpro](#)

DOI

dx.doi.org/10.17504/protocols.io.261ge5emyg47/v2

Blake Balcomb^{1,2,3}, Peter Marples^{1,2,3}, Lizbé Koekemoer^{4,5,3}, Charlie Tomlinson^{1,2,3}, Daren Fearon^{1,2,3}

¹Diamond Light Source; ²Research Complex at Harwell; ³ASAP Discovery Consortium;

⁴Centre of Medicines Discovery; ⁵University of Oxford

Blake Balcomb: The principle crystallographer on the SARS-CoV-2 Mpro project.;

ASAP Discovery



Peter Marples

Diamond Light Source



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261ge5emyg47/v2>

External link: <https://asapdiscovery.org/outputs/target-enabling-packages/#ASAP-COV-MPRO>

Protocol Citation: Blake Balcomb, Peter Marples, Lizbé Koekemoer, Charlie Tomlinson, Daren Fearon 2025. Crystallisation of SARS-CoV-2 MPro in P212121 for compound soaking. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.261ge5emyg47/v2> Version created by **Mary-Ann Xavier**

Manuscript citation:

Melissa L. Bobby *et al.*,
Open science discovery of potent noncovalent SARS-CoV-2 main protease inhibitors. *Science* **382**, eabo7201 (2023). DOI: [10.1126/science.abo7201](https://doi.org/10.1126/science.abo7201)

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 24, 2025

Last Modified: May 13, 2025

Protocol Integer ID: 123294

Keywords: crystallisation, XChem, ASAP, AVIDD, Diamond Light Source, i04-1, Mpro, SARS-CoV-2, crystallisation of sar, p212121 for compound, reproducible sar, crystallization protocol, xchem fragment screening, severe acute respiratory syndrome coronavirus, crystallisation, solvent test, strategies for pandemic preparedness, pandemic preparedness, sar, cov mpro

Funders Acknowledgements:

National Institutes of Health/National Institute Of Allergy and Infectious Diseases (NIH/NIAID)

Grant ID: Grant ID: U19AI171399

Disclaimer

The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.

Acknowledgements:

Diamond Light Source Ltd, Harwell Science and Innovation Campus, Didcot OX11 0QX, UK
Research Complex at Harwell, Harwell Science and Innovation Campus, Didcot OX11 0FA, UK
Oxford Lab Technologies crystal shifter <https://doi.org/10.1107/S2059798320014114>

Abstract

The COVID-19 pandemic has highlighted the need to identify novel therapeutic interventions and strategies for pandemic preparedness against Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2). This protocol outlines the crystallization protocol and buffer conditions used to obtain reproducible SARS-CoV-2 Mpro crystals suitable for **XChem** fragment screening. In this new version we added the solvent test, soaking conditions, seeding reference paper and addgene ids used in this experiment.

Materials

SwissCI 3 lens crystallization plates <https://swissci.com/product/3-lens-crystallisation-plate/> **Codes:**

Midi: UVXPO-3LENS 3W96T-PS 3W96T-UVP

MES adjusted to with HCl, Molecular Dimensions, Catalog # MD2-013-PH 6.7
 50% w/v PEG 4000, Molecular Dimensions, Catalog # MD2-250-11
 99.9% DMSO, Molecular Dimensions, Catalog # MD2-250-159

Purified SARS-CoV-2 Mpro protein () in HEPES, ,
 NaCl, 5% glycerol, TCEP

Protocol materials

 SARS-CoV-2 isolate: Wuhan-Hu-1 Mpro (3CL protease) **addgene Catalog #228666**

 SARS-CoV-2 Mpro protease (C145S mutant - catalytically inactive); aliases: 3CLpro, 3C-like **addgene Catalog #204769**

Troubleshooting

Safety warnings

 Follow all handling warning for the chemicals used in the crystallisation screen composition.



SARS-CoV-2 Mpro expression and purification

- 1 The protein used for crystallisation was expressed and purified using the following protocol and plasmid



SARS-CoV-2 isolate: Wuhan-Hu-1 Mpro (3CL protease) **addgene Catalog #228666**

Protocol

NAME

SARS-CoV-2 Mpro small scale expression and purification protocol

CREATED BY

Korvus Wang

[Preview](#)

Equipment needed

- 2 **Formulatrix Rock Imager** (or incubator of choice)
SPT mosquito

Equipment

Mosquito HV

NAME

High Volume 16-Channel Robotic Liquid Handler

TYPE

SPT LabTech

BRAND

3097-01057

SKU

<https://www.sptlabtech.com/products/liquid-handling/mosquito-hv/>^{LINK}

P100 8 multi-channel pipette



SwissCI 3 lens plate

Crystallisation experiment

1d

3 Preparation of seed stock using imature



SARS-CoV-2 Mpro protease (C145S mutant - catalytically inactive); aliases:
3CLpro, 3C-like **addgene Catalog #204769**

Mpro and mature protein: Protocol for obtaining SARS-CoV-2 Mpro P2₁2₁ crystal form can be found here: DOI: <https://doi.org/10.1016/j.jmb.2021.167118>.

4 Prepare **seed** stock:

17m 40s

Protocol



NAME

Diamond XChem Seeding Protocol

CREATED BY

Peter Marples

[Preview](#)

1: 250 dilution Sample seeds

5 Protein and buffer requirements:



43.2 µL

[M] 5 mg/mL



Sample



1.92 mL Crystallisation screen



14.4 µL

Sample

P2₁2₁ crystal form seeds (see step 3), dilution 1:250

6 Crystallisation screen composition:



0.1 Molarity (M)

MES



6.7

11% w/v PEG 4000

5% DMSO

Stock solutions used:



1 Molarity (M)

MES adjusted to



6.7

with HCl



50% w/v PEG 4000

99.9% DMSO

Note

The crystallisation screen can be stored in a duran bottle or aliquoted into 96 deep well block for easy dispensing into SwissCI 3 lens plates.

For long term storage keep the crystallisation screen in the fridge at 4°C.

- 7 Dispense  20 µL Crystallisation screen into SwissCI 3 lens plate reservoir wells using a 100 µl multi-channel pipette.

10m

Dispense  150 nL  Sample to each lens using the SPT mosquito.

Dispense  150 nL Crystallisation screen to each lens using the SPT mosquito.

Dispense  50 nL Seeds P2₁2₁2₁ stock to each lens using the SPT mosquito.

Drop ratio: 3:3:1 ratio (150 nL  Sample : 150 nL reservoir solution: 50 nl seeds)

Final drop volume: 350 nl

- 8 Incubate at  20 °C for  24:00:00 h in Formulatrix Rock Imager.

1d

Imaging Schedule: The first images are taken after 12 h and the imaging schedule follows a Fibonacci sequence of days for further collections.

- 9 Crystal form after ~24 h.

Expected result

The crystals reach their maximum size after 48 h.

Crystals typically form either as stacked plates or in small clusters containing 4-6 crystals.

Morphology: plates.


Size: ~100 μm in length and ~20 μm in width, depth of the crystals is ~10 μm

Appearance: glass shard.

Average resolution: 1.8 $\text{\AA}$

Space group: $P2_12_12_1$

Unit cell: 67 $\text{\AA}$, 99 $\text{\AA}$, 102 $\text{\AA}$
90.00°, 90.00°, 90.00°



An example of a drop containing SARS-CoV-2 Mpro protease crystals.



Data collection at Synchrotron

- 10 Diamond Light Source
Unattended Data Collection (UDC)
Data Collection Temperature: 100K
Detector: DECTRIS EIGER2 X 9M
Beamline: I04-1
Wavelength: 0.9212 Å
Resolution (Å): 1.78
Beam Size (µm): 60 X 50
Number of images: 3600
Oscillation: 0.10°
Exposure (s): 0.0020
Transmission (%): 100
Flux (ph/s): 9.50e+11

Soaking Tolerance

- 11 Crystals can tolerate up to 20% DMSO, 1-2h

Protocol



NAME

XChem solvent test

CREATED BY

Peter Marples

[Preview](#)

Compound Soaking

- 12 Crystals can tolerate up to 20% DMSO, 1-2h



Protocol



NAME

XChem crystallographic fragment screening

CREATED BY

Peter Marples

Preview

Protocol references

For SARS-CoV-2 Main protease in orthorhombic space group

<https://www.rcsb.org/structure/7MBG>

<https://www.sciencedirect.com/science/article/pii/S0022283621003429#s0045>

Group deposition SARS-CoV-2 main protease in complex with inhibitors from the COVID Moonshot

<https://www.rcsb.org/structure/7GMB>

<https://www.science.org/doi/10.1126/science.abo7201>