

Jun 12, 2019

CTAB DNA Extraction for high quality/molecular weight DNA

DOI

<https://dx.doi.org/10.17504/protocols.io.3rsgm6e>



Yaowu Yuan¹

¹University of Connecticut

High molecular weight D...

Mimulus



Andrea Sweigart

University of Georgia

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3rsgm6e>

External link:

[http://mimubase.org/FTP/Protocols/DNA_extraction/CTAB%20DNA%20Extraction%20\(High%20Molecular%20Weight\).pdf](http://mimubase.org/FTP/Protocols/DNA_extraction/CTAB%20DNA%20Extraction%20(High%20Molecular%20Weight).pdf)

Protocol Citation: Yaowu Yuan 2019. CTAB DNA Extraction for high quality/molecular weight DNA. protocols.io
<https://dx.doi.org/10.17504/protocols.io.3rsgm6e>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 05, 2019

Last Modified: June 12, 2019

Protocol Integer ID: 24082

Keywords: DNA extraction, Mimulus, ctab dna extraction for high quality, ctab dna extraction, dna

Attachments



[CTAB DNA Extraction ...](#)

28KB

Guidelines

For Safety Warnings and Hazard Information please refer to the SDS (Satety Data Sheet).



Materials

MATERIALS

⊗ liquid nitrogen

⊗ 70% Ethanol

⊗ CTAB DNA Extraction buffer

⊗ Chloroform: IsoAmyl Alcohol (24:1)

⊗ 7.5M Ammonium acetate

⊗ 100% Ethanol

⊗ dH₂O

CTAB DNA Extraction Buffer (Recipe to make 100 mL)

10 mL 1 M Tris Buffer

8.3 g NaCl (1.4 M)

0.744 g EDTA

2 g CTAB

2 g PVP

0.088 g Asorbic acid



- 1 Grind plant tissue in a mortar cooled with liquid nitrogen.
 - 2 Add  750 μ L CTAB DNA Extraction buffer .
 - 3 Wait until it warms up and becomes a green paste, then transfer to an eppie tube.
 - 4 Incubate the CTAB/plant extract mixture for  00:15:00 at  55 °C in the heat block and invert to mix throughout the 15 minutes.
 - 5 Add  500 μ L Chloroform : IsoAmyl Alcohol (24:1) in the hood and mix the solution by inverting the tubes (**do not vortex**).
 - 6 Centrifuge at  13000 rpm for  00:10:00 .
 - 7 Transfer the upper aqueous phase **only** to a new eppie tube (~  500 μ L).
 - 8 Add RNase A (10 μ g/ml).
- Note**
-  5 μ L of 1mg/ml stock if you have  500 μ L of sample.
- 9 Incubate at  37 °C for  00:30:00 .
 - 10 Add  50 μ L 7.5M Ammonium acetate followed by  500 μ L ice cold 100% ethanol and invert to mix.



11 Put tubes in -20 °C freezer for 1 hour (or longer) to precipitate the DNA.

12 Centrifuge at 13000 rpm for 00:15:00 .

Note

You should see a pellet at the bottom (align the tubes so that you know where the pellet is in case you can't see it very well).

13 Remove the supernatant and wash the DNA pellet as follows. (1/2)

13.1 Add 500 µL ice cold 70% ethanol . (1/2)

13.2 Centrifuge at 13000 rpm for 00:05:00 . (1/2)

14 Remove the supernatant and wash the DNA pellet as follows. (2/2)

14.1 Add 500 µL ice cold 70% ethanol . (2/2)

14.2 Centrifuge at 13000 rpm for 00:05:00 . (2/2)

15 Remove all the supernatant and allow the DNA pellet to dry in the hood (approx. 00:20:00).

Note

Do not over dry the pellet since it will be hard to re-dissolve.



16 Resuspend the DNA in  50 μ L dH₂O .

17 NanoDrop the sample to estimate the concentration.

Note

Alternatively the DNA can be run on a gel to estimate the concentration or size of the DNA. Running a 0.4% gel overnight with the lambda DNA mono-cut ladder can give you an estimate of the size (it still doesn't separate the large bands very well).