

Oct 12, 2022

CTAB DNA Extraction Protocol for Mollusks

DOI

dx.doi.org/10.17504/protocols.io.e6nvwjpr2lmk/v1

Lavanya M Vythalingam¹, Tim Regan¹, Tim Bean¹, Carolina Peñaloza¹

¹The Roslin Institute and Royal (Dick) School of Veterinary Studies, University of Edinburgh



Lavanya M Vythalingam

Heriot-Watt University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwjpr2lmk/v1>

Protocol Citation: Lavanya M Vythalingam, Tim Regan, Tim Bean, Carolina Peñaloza 2022. CTAB DNA Extraction Protocol for Mollusks. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.e6nvwjpr2lmk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 06, 2022

Last Modified: October 12, 2022

Protocol Integer ID: 69614

Keywords: comprehensive ctab extraction protocol, dna extraction of mollusk, dna extraction protocol, dna extraction, mollusk, extraction, dna

Abstract

This is a comprehensive CTAB extraction protocol intended for DNA extraction of mollusks. It includes detailed methodology and several reagent preparation.

Protocol materials

⊗ CTAB extraction solution **Teknova Catalog #C2190**

⊗ Proteinase K

⊗ Chloroform/Isoamyl Alcohol (24:1) **Acros Organics Catalog #327155000**

⊗ Chloroform/Isoamyl Alcohol (24:1) **Acros Organics Catalog #327155000**

⊗ RNase A

⊗ Ethanol

Troubleshooting



Method

- 1
Slice tissue sample using a clean, sterile scalpel. Put sample into a 1.5 ml tube with  400 μL  CTAB extraction solution **Teknova Catalog #C2190** .
- 2
Add  10 μL  Proteinase K  20 mg/mL and invert tube to mix.
- 3
Incubate at  50 $^{\circ}\text{C}$ for 60 minutes until lysed (or longer depending on tissue thickness). Invert tube occasionally to mix sample.
- 4
Add  400 μL  Chloroform/Isoamyl Alcohol (24:1) **Acros Organics Catalog #327155000** . Gently shake tube to emulsify sample.
- 5
Leave at room temperature for 2 minutes and shake tube to mix.
- 6
Spin tube at 13 000 rpm for 2 minutes and transfer upper layer to a 1.5 ml tube.
- 7
Add  5 μL of  RNase A  10 mg/mL . Leave for 30 minutes at  37 $^{\circ}\text{C}$.
- 8
Add 1 equal volume of  Chloroform/Isoamyl Alcohol (24:1) **Acros Organics Catalog #327155000** .
- 9
Spin tube at 13 000 rpm for 2 minutes and transfer upper layer to a new 1.5 ml tube containing  900 μL **CTAB dilution solution**. Gently invert tube to mix. Leave in fridge overnight (**optional**).
- 10



Spin tube at 13 000 rpm for 10 minutes.

11

Discard supernatant with a micropipette and add  1000 μL **0.4 M NaCl in TE**. Invert tube to wash.

12

Spin tube at 13 000 rpm for 5 minutes and remove supernatant.

13

Add  300 μL **1.42 M NaCl in TE** to dissociate the CTAB from the DNA. Invert tube until pellet becomes transparent.

Note

May require gentle shaking at  37 °C

14

Add  600 μL  Ethanol (previously left at  -20 °C). Invert manually until precipitation appears.

Note

Leave precipitating overnight for a higher DNA concentration

15

Centrifuge at full speed for 8 minutes, then remove supernatant. Dry at room temperature or in heat block at  37 °C

16

Resuspend in  50 μL **EB buffer**. Use lesser volume of **EB buffer** if a higher concentration is required.



Notes

17

All centrifugations were performed using a benchtop centrifuge (Eppendorf - 5415 D)

18

The Nucleic acid quality is measured using a Nanodrop, and the A260/A230 ratio of 1.8-2.0 indicates good DNA purity. Qubit can be used to determine nucleic acid concentration.

19

Store DNA sample at -20°C until further use.

Reagent preparation

20

CTAB Extraction Solution

Note

🧪 11.69 g ⚗️ Sodium Chloride **Catalog #PubChem CID: 5234**

🧪 10 mL

⚗️ Tris HCl Buffer 1M Solution, Sterile pH 7.5 **Bio Basic Inc. Catalog #SD8124.SIZE.450ml**

🧪 5 mL ⚗️ EDTA (0.5 M), pH 8.0 **Life Technologies Catalog #AM9260G**

🧪 3.33 g ⚗️ CTAB

Water to 🧪 100 mL

20.1

CTAB Dilution Solution



Note

10 mL

 Tris HCl Buffer 1M Solution, Sterile pH 7.5 **Bio Basic Inc. Catalog #SD8124.SIZE.450ml** 4 mL EDTA (0.5 M), pH 8.0 **Life Technologies Catalog #AM9260G**

2 g CTAB

Water to 200 mL

20.2

0.4M NaCl in TE

Note

 16 mL **5M** NaCl

2 mL

 Tris HCl Buffer 1M Solution, Sterile pH 7.5 **Bio Basic Inc. Catalog #SD8124.SIZE.450ml** 400 μ L EDTA (0.5 M), pH 8.0 **Life Technologies Catalog #AM9260G**

Water to 200 mL

20.3

1.42 M NaCl in TE

Note

 24.8 mL **5M** NaCl

1 mL

 Tris HCl Buffer 1M Solution, Sterile pH 7.5 **Bio Basic Inc. Catalog #SD8124.SIZE.450ml** 200 μ L EDTA (0.5 M), pH 8.0 **Life Technologies Catalog #AM9260G**

Water to 100 mL