

Jul 30, 2023

Culturing and passaging of iPSC derived intestinal organoids

DOI

dx.doi.org/10.17504/protocols.io.4r3l22or3l1y/v1

rachel.bates Bates¹

¹UCL



rachel.bates Bates

ucl

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.4r3l22or3l1y/v1

Protocol Citation: rachel.bates Bates 2023. Culturing and passaging of iPSC derived intestinal organoids. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.4r3l22or3l1y/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 30, 2023



Last Modified: August 28, 2024

Protocol Integer ID: 85691

Keywords: ASAPCRN, stemdiff intestinal organoid kit, intestinal organoid, intestinal organoid kit, organoid, passaging of ipsc, using stemdiff

Disclaimer

Protocol Particulars Video

The video below is a supplement with extra context and tips, as part of the Aligning Science Across Parkinson's (ASAP) Protocol Particulars video interview series, featuring conversations with protocol authors.

<https://www.youtube.com/embed/IA1YeJurpug?si=nb3HPUyzMtBiAKq4>

Abstract

Culturing and passaging of iPSC derived intestinal organoids derived using STEMDIFF intestinal organoid kit. We usually use organoids after 5 passages once consistent growth has been established and until 15th passage.



Materials

Intestinal growth medium.

⊗ DMEM F12/HEPES **Gibco - Thermo Fischer Catalog #113300**

⊗ N-2 max supplement **R&D Systems Catalog #AR009**

⊗ N21-MAX Media Supplement (50X) **R&D Systems Catalog #AR008**

⊗ HEPES **Merck MilliporeSigma (Sigma-Aldrich) Catalog #H6147**

FC WNT

EGF

Noggin

SB202

Lacetylcysteine

Nicotinamide

L-Gastrin

A3801



Protocol materials

- ✕ N-2 max supplement **R&D Systems Catalog #AR009**
- ✕ N21-MAX Media Supplement (50X) **R&D Systems Catalog #AR008**
- ✕ HEPES **Merck MilliporeSigma (Sigma-Aldrich) Catalog #H6147**
- ✕ DMEM F12/HEPES **Gibco - Thermo Fisher Scientific Catalog #113300**
- ✕ STEMdiff intestinal organoid kit **STEMCELL Technologies Inc. Catalog #05140**
- ✕ DMEM F12/HEPES **Gibco - Thermo Fisher Scientific Catalog #113300**
- ✕ DMEM F12/HEPES **Gibco - Thermo Fisher Scientific Catalog #113300**
- ✕ Anti-Adherence Rinsing Solution **STEMCELL Technologies Inc. Catalog #07010**
- ✕ ACCUTASE™ 100 mL **STEMCELL Technologies Inc. Catalog #7920**
- ✕ DMEM F12/HEPES **Gibco - Thermo Fisher Scientific Catalog #113300**
- ✕ Cultrex® 3-D Culture Matrix™ Reduced Growth Factor Basement Membrane Extract, PathClear® **Merck MilliporeSigma (Sigma-Aldrich) Catalog #3445-001-01**
- ✕ Cultrex® 3-D Culture Matrix™ Reduced Growth Factor Basement Membrane Extract, PathClear® **Merck MilliporeSigma (Sigma-Aldrich) Catalog #3445-001-01**
- ✕ Nunc® Cell-Culture Treated Multidishes, 12 well **Thermo Fisher Catalog #150628**
- ✕ DMEM F12/HEPES **Gibco - Thermo Fisher Scientific Catalog #113300**
- ✕ Anti-Adherence Rinsing Solution **STEMCELL Technologies Inc. Catalog #07010**
- ✕ Cultrex® 3-D Culture Matrix™ Reduced Growth Factor Basement Membrane Extract, PathClear® **Merck MilliporeSigma (Sigma-Aldrich) Catalog #3445-001-01**
- ✕ DMEM F12/HEPES **Gibco - Thermo Fisher Scientific Catalog #113300**
- ✕ DMEM F12/HEPES **Gibco - Thermo Fisher Scientific Catalog #113300**
- ✕ DMEM F12/HEPES **Gibco - Thermo Fisher Scientific Catalog #113300**
- ✕ Cultrex® 3-D Culture Matrix™ Reduced Growth Factor Basement Membrane Extract, PathClear® **Merck MilliporeSigma (Sigma-Aldrich) Catalog #3445-001-01**
- ✕ Nunc® Cell-Culture Treated Multidishes, 12 well **Thermo Fisher Catalog #150628**



Establishing organoids.

- 1 Intestinal organoids were generated using
 STEMdiff intestinal organoid kit **STEMCELL Technologies Inc. Catalog #05140**
until protocol stage 6.2.2.4. For some cell lines yields of organoids can be very low using STEMcell methods. We usually harvest organoids at 9 days differentiation.
- 2 Prepare required number of 15mL falcon by washing with
 Anti-Adherence Rinsing Solution **STEMCELL Technologies Inc. Catalog #07010** .
Rinse with PBS.
- 3 Thaw an aliquot of
 Cultrex® 3-D Culture Matrix™ Reduced Growth Factor Basement Membrane Extract, PathClear® **Merck MilliporeSigma (Sigma-Aldrich) Catalog #3445-001-01**
on ice.  60 µL will be required for each well.
- 4 Monolayer cultures in 24 well plate displaying spheroid budding are washed 3 times with cold  DMEM F12/HEPES **Gibco - Thermo Fischer Catalog #113300** to remove spheroids and place in precoated falcon.
- 4.1 if yield is low add  1 mL 2m
 ACCUTASE™ 100 mL **STEMCELL Technologies Inc. Catalog #7920** to monolayer for  00:02:00 at  37 °C . Remove detached monolayer in accutase and add to falcon.
- 5 Make volume up  10 mL with
 DMEM F12/HEPES **Gibco - Thermo Fischer Catalog #113300** and allow organoids and fragments of tissue to sink
- 6 Carefully aspirate supernatant leaving  1 mL behind in falcon containing organoids. 5m
Add a new 1ml of  DMEM F12/HEPES **Gibco - Thermo Fischer Catalog #113300** to the spheroids. Centrifuge  300 x g, 4°C, 00:05:00 .



- 7 Carefully remove as much supernatant as possible. Place falcon on ice and add  60 μ L  Cultrex® 3-D Culture Matrix™ Reduced Growth Factor Basement Membrane Extract, PathClear® Merck MilliporeSigma (Sigma-Aldrich) Catalog #3445-001-01 for each well required. Recommend plating 100-200 domes per dome (2 domes per well).
- 8 Gently mix organoids into the Cultrex being careful not to introduce bubbles. Using a new pipette tip pick up  30 μ L of mix and transfer to one well of a  Nunc® Cell-Culture Treated Multidishes, 12 well Thermo Fisher Catalog #150628 30m
- . Repeat the process so that the well contains 2 domes. Allow domes to solidify at  37 °C in incubator for  00:30:00 . We find organoids distribute best if turn the plate upside down. We have found that nunc plates provide the best adherence for domes without detachment over the 10 day growth period.
- 9 Carefully add  1 mL of Intestinal organoid growth medium. Incubate at  37 °C . Feed cells every 2/3 days.

Passaging intestinal organoids.

- 10 Coat a 15ml Falcon with  Anti-Adherence Rinsing Solution STEMCELL Technologies Inc. Catalog #07010 before washing off with PBS.
- 11 Thaw an aliquot of  Cultrex® 3-D Culture Matrix™ Reduced Growth Factor Basement Membrane Extract, PathClear® Merck MilliporeSigma (Sigma-Aldrich) Catalog #3445-001-01 ready for passaging. Usually split a well from a 12 well plate 1:4 and allow  60 μ L for each new well.
- 12 After 10 days in culture domes are ready for passaging. Remove medium from well. Add  1 mL of cold  DMEM F12/HEPES Gibco - Thermo Fischer Catalog #113300 . Gently break up the pellet and until it comes loose from the plate. Add into prepared falcon tube. Combine up to 4 wells from 12 well plate.



- 13 Wash well with a further  1 mL
 DMEM F12/HEPES **Gibco - Thermo Fischer Catalog #113300** to collect any remaining organoids and add to falcon. Make up to  10 mL .
- 14 Keep falcon on ice for 5 minutes to allow organoids to sink. Carefully remove supernatant leaving  1 mL in falcon with organoids. Add a new  2 mL of  DMEM F12/HEPES **Gibco - Thermo Fischer Catalog #113300** . Using a 1ml pipettor pipette suspension up and down 5 times to gently break apart the domes. incubate on ice for  00:05:00 . 5m
- 15 Centrifuge  300 x g, 4°C, 00:05:00 . Carefully aspirate supernatant. Cultrex will form a cloudy layer on pellet at bottom of tube. try to remove as much of this as possible to clear single cells and smaller fragments. Add  2 mL  DMEM F12/HEPES **Gibco - Thermo Fischer Catalog #113300** to pellet. Using a 1ml pipette vigorously pipette organoids up and down 20-25 times. Incubate on ice for  00:05:00 10m
- 16 Centrifuge pellet  300 x g, 4°C, 00:05:00 . Carefully remove as much supernatant as possible. Place falcon on ice and add  60 µL  Cultrex® 3-D Culture Matrix™ Reduced Growth Factor Basement Membrane Extract, PathClear® **Merck MilliporeSigma (Sigma-Aldrich) Catalog #3445-001-01** per 12 well required. Carefully mix organoid fragments with cultrex ensuring not to introduce bubbles. Split into two  30 µL domes per well of  Nunc® Cell-Culture Treated Multidishes, 12 well **Thermo Fisher Catalog #150628** . Incubate at  37 °C for  00:30:00 . 35m
- 17 Add  1 mL of intestinal organoid growth medium per well.

Intestinal organoid growth medium

18

