



Mar 05, 2025

Culturing RBL-2H3 rat basophils

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.14egn67qpl5d/v1>

Justin Donnelly¹, Emily C.P. Weiss¹

¹Arcadia Science

Arcadia Science



Arcadia Science

Arcadia Science

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn67qpl5d/v1>

Protocol Citation: Justin Donnelly, Emily C.P. Weiss 2025. Culturing RBL-2H3 rat basophils. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.14egn67qpl5d/v1>

Manuscript citation:

Borges AL, Donnelly J, Morazan E, Rollins M. (2025). Compound 48/80 is toxic in HMC1.2 and RBL-2H3 cells.
<https://doi.org/10.57844/arcadia-3207-4695>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 10, 2024

Last Modified: March 14, 2025

Protocol Integer ID: 99551

Keywords: tissue culture, mast cell, basophil

Abstract

Guidelines for thawing, maintaining, subculturing, and freezing RBL-2H3 (rat basophil) line.

Materials

 RBL-2H3 rat basophil cells **ATCC Catalog #CRL-2256**

 T75 flasks for cell culture **Thermofisher**

 Mast Cell Degranulation Assay Kit **Merck MilliporeSigma (Sigma-Aldrich) Catalog #IMM001**

 Eagles Minimum Essential Medium **VWR International (Avantor) Catalog #MSPP-302003**

 Heat-inactivated fetal bovine serum **Thermofisher Catalog #A3840101**

 70% ethanol **VWR International (Avantor)**

 Dimethyl sulfoxide (DMSO) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D2438**

T25 flasks

 IMDM **Gibco - Thermo Fisher Scientific Catalog #12440053**

 PenStrep **Invitrogen - Thermo Fisher**

 1-Thioglycerol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M6145**

 Trypsin-EDTA (0.25%) **Thermo Fisher Scientific Catalog #25200056**

PBS

15 mL centrifuge tubes



Protocol materials

- ☒ 70% ethanol **VWR International (Avantor)**
- ☒ PenStrep **Invitrogen - Thermo Fisher**
- ☒ 1-Thioglycerol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M6145**
- ☒ RBL-2H3 rat basophil cells **ATCC Catalog #CRL-2256**
- ☒ T75 flasks for cell culture **Thermofisher**
- ☒ Trypsin-EDTA (0.25%) **Thermo Fisher Scientific Catalog #25200056**
- ☒ Dimethyl sulfoxide (DMSO) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D2438**
- ☒ Mast Cell Degranulation Assay Kit **Merck MilliporeSigma (Sigma-Aldrich) Catalog #IMM001**
- ☒ Eagles Minimum Essential Medium **VWR International (Avantor) Catalog #MSPP-302003**
- ☒ Heat-inactivated fetal bovine serum **Thermofisher Catalog #A3840101**
- ☒ IMDM **Gibco - Thermo Fisher Scientific Catalog #12440053**

Thawing the RBL-2H3 rat basophil cell line

8m

- 1 To make the complete culture medium for these cells, add Heat-inactivated fetal bovine serum **Thermofisher Catalog #A3840101** to a final concentration of 1M 15 % (v/v) in Eagles Minimum Essential Medium **VWR International (Avantor) Catalog #MSPP-302003**.

. Add PenStrep **Invitrogen - Thermo Fisher** to a final concentration of 1×. Sterile-filter the media under strict aseptic conditions.

It's important to avoid excessive alkalinity of the medium during recovery of the cells. Prior to adding cells, you should place the culture vessel containing the complete growth medium into the incubator for at least 00:15:00 to allow the medium to reach 7.0 – 7.6 .
- 2 The frozen cell culture should have been stored in liquid nitrogen and not at -70 °C , as storage at -70 °C will result in loss of viability. To use the cells, first thaw the vial by gentle agitation in a 37 °C water or bead bath. To reduce the possibility of contamination, keep the O-ring and cap out of the water (if using). Thawing should be rapid (approximately 00:02:00 for water bath, 00:05:00 for bead bath).
- 3 Remove the vial from the water bath as soon as the contents are thawed, and decontaminate the outside of the vial by dipping in or spraying with 70% ethanol **VWR International (Avantor)** . All of the operations from this point on should be carried out under strict aseptic conditions.
- 4 Working quickly to avoid loss of viability, transfer the vial contents to a centrifuge tube containing 9.0 mL complete culture medium and spin at approximately 125 x g, 00:05:00 .
- 5 Resuspend cell pellet with the recommended complete medium (see the specific batch information for the recommended dilution ratio) and dispense into a 75 cm² culture flask.

15m

7m

5m



- 6 Incubate the culture at  37 °C in a suitable incubator with a 5% CO₂ in air atmosphere.

Subcultivation

10m

- 7 Subculture the cells every 2–3 days. We recommend a subcultivation ratio of 1:4 to 1:8.

- 7.1 Remove and discard used culture medium.

- 7.2 Rinse the cell layer with  Trypsin-EDTA (0.25%) **Thermo Fisher Scientific Catalog #25200056** or PBS to remove traces of serum.

- 7.3 Add  2.0 mL to  3.0 mL of  Trypsin-EDTA (0.25%) **Thermo Fisher Scientific Catalog #25200056** solution to the flask. Incubate cells at 37 °C until they detach — for RBL-2H3, this is usually within  00:10:00 .

10m

- 7.4 Add  6.0 mL to  8.0 mL of complete growth medium ($\geq 2\times$ the trypsin volume) and dislodge remaining cells in the flask by gently pipetting.

- 7.5 Add appropriate aliquots of the cell suspension to new culture vessels and incubate the culture at  37 °C in a suitable incubator. We recommend 5% CO₂ in air atmosphere.

Cryopreservation

8m

- 8 Cryopreserve many aliquots of low-passage cells for future experiments.

- 8.1 Recover cells from flask as in 7.1–7.4.



- 8.2 Count cells. Prepare a mixture of  10 μL cells and  10 μL trypan blue, inverting cells gently before removing sample for counting. Mix well, then transfer  10 μL trypan-diluted cells to a Countess slide. Count using Countess to determine cell density and relative viability.
- 8.3 Freeze cells at ~1 million live cells in  1 mL cryo-media per vial. Ensure there are enough live cells for the desired number of aliquots. Centrifuge cells to be frozen at  125 x g, 00:05:00 . 5m
- 8.4 Meanwhile, prepare cryo-media by diluting  Dimethyl sulfoxide (DMSO) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D2438** in complete growth medium to a final concentration of  10 % (v/v) .
- 8.5 After centrifugation, aspirate supernatant. Resuspend cell pellet in cryo-media. Working quickly to avoid loss of viability, aliquot cells to pre-labeled screw-top cryotubes.
- 8.6 Place cells in a Mr. Frosty or similar container that controls freezing rate. Place container at  -80 °C  Overnight . 10m
- 8.7 Transfer frozen cells to vapor-phase liquid nitrogen for long-term storage.

Protocol references

HMC1.2 Manufacturer Guidelines:

https://www.emdmillipore.com/US/en/product/HMC-1.2-Human-Mast-Cell-Line,MM_NF-SCC062?ReferrerURL=https%3A%2F%2Fwww.google.com%2F#documentation

RBL-2H3 Manufacturer Guidelines:

<https://www.atcc.org/api/pdf/product-sheet?id=CRL-2256>