

Nov 01, 2019

# Culturing THP-1 Cells

DOI

dx.doi.org/10.17504/protocols.io.y6gfzbw

Teesha C Luehr<sup>1</sup>

<sup>1</sup>University of British Columbia

Leonard Foster's Lab



### Teesha C Luehr

University of British Columbia

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account





DOI: https://dx.doi.org/10.17504/protocols.io.y6gfzbw

Protocol Citation: Teesha C Luehr 2019. Culturing THP-1 Cells. protocols.io

https://dx.doi.org/10.17504/protocols.io.y6gfzbw

Manuscript citation:

https://www.atcc.org/Products/All/TIB-202.aspx



**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 13, 2019

Last Modified: November 01, 2019

**Protocol Integer ID: 21416** 

**Keywords:** human monocyte suspension cell line, human monocyte suspension cell line from peripheral blood, acute mnocytic leukemia, cell

### **Abstract**

THP-1 cells are a human monocyte suspension cell line from peripheral blood of a 1 year old infant who had acute mnocytic leukemia.

## **Troubleshooting**



### **Preparing Media**

- 1 The base medium for this cell line is RMPI-1640.
- 2 Reguired supplements:
  - [M] 1 % volume L-glutamine
  - [M] 10 % volume Fetal Bovine Serum

#### Note

Most catalog numbers of RMPI-1640 contain L-glutamine, however, some do not. Ensure that it is in the media before using it for culturing.

- 3 **Optional Supplements:** 
  - [M] 1 % volume PenStrep
  - [M] 0.05 millimolar (mM) 2-mercaptoethanol

#### Note

PenStrep is not required for THP-1 culturing, however, if you are having issues with bacterial contamination, it can be used at 1X.

#### Note

2-mercaptoethanol is stated as a required component for complete RPMI-1640 medium, however, in our laboratory it is not standard practice to add it.

## **Cell Storage**

4 Always store cells in liquid nitrogen. This is for both the original tube of cells from ATCC and any passages afterwards.

## **Preparation of Materials & Reagents**

5 Place the media bottle in the \$\\$37 \circ\$ water bath at least \( \frac{1}{10} \) 00:30:00 prior to using



6 Thaw cells at \$\mathbb{8}\$ 25 °C (room temperature) for \bigode{\infty} 00:10:00 or **\$** 37 °C in a water bath for **♦** 00:02:00

## Working from fozen cells

- 7 Sanitize all items going into the Biological Safety Cabinet with 70% ethanol
- 8 As soon as the cells are thawed, transfer the cells to a  $\perp$  15 mL conical tube and add

▲ 10 mL of complete media

Note

Cells are stored with 5% DMSO, which can lyse cells if they are left for too long.

- 9 Pellet cells for 00:03:00 at 500g
- 10 Discard supernatant
- 11 Resuspend cells by pipetting up and down 5X in 4 5 mL complete media
- 12 Transfer cells + media to a T-25 flask

### Incubation

13 Incubate cells at 37 °C and 5% CO2 and 80% humidity

## Feeding and Splitting



- 14 THP-1 cells replicate after ~26 hours. In practice, it takes 2 days for a true doubling.
- 15 Once cells have doubled OR when media has begun to change colour, it is time to add media, split cells into new flasks, or to spin down to remove all media

STEP CASE

#### Adding media 4 steps

- 16 If concerned about cell concentration, perform a cell count
- 17 Double the total media volume with new complete media
- 18 Carefully mix the new media in by rocking the flask back and forth
- 19 Place the flask back in the incubator