

Nov 04, 2024

Current clamp recordings

DOI

dx.doi.org/10.17504/protocols.io.n92ldr3n8g5b/v1

Chuyu Chen¹, Loukia Parisiadou¹, Yevgenia Kozorovitskiy²

¹Northwestern University, Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815;

²Department of Neurobiology, Northwestern University, Evanston, IL, USA



Chuyu Chen

Northwestern University, Aligning Science Across Parkinson's...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldr3n8g5b/v1>

Protocol Citation: Chuyu Chen, Loukia Parisiadou, Yevgenia Kozorovitskiy 2024. Current clamp recordings. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n92ldr3n8g5b/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: November 04, 2024

Last Modified: November 04, 2024

Protocol Integer ID: 111535

Keywords: ASAPCRN, Irrk2 effects on intrinsic ispns adaptation, intrinsic ispns adaptation, egfp bac reporter mice, Irrk2 effect, vivo striatal slice, current clamp recording, Irrk2, ispns in ex

Funders Acknowledgements:

Aligning Science Across Parkinson's [ASAP-020600] through the Michael J. Fox Foundation for Parkinson's Research (MJFF)
Grant ID: ASAP-020600

Abstract

To better understand LRRK2 effects on intrinsic iSPNs adaptations, *Drd2-eGFP* BAC reporter mice were injected with haloperidol, haloperidol + MLI-2, or vehicle controls for 14 days. Then, we performed whole-cell current-clamp recordings on identified iSPNs in ex vivo striatal slices.

Troubleshooting

Acute Slice Preparation

- 1 *Drd2-eGFP* mice were treated with either haloperidol or haloperidol plus LRRK2 inhibitor (MLi-2) for 14 days.
- 2 On day 15 or 16, mice were deeply anesthetized on isoflurane (3%) prior to perfusion with ice-cold aCSF containing in (mM): 127 NaCl, 2.5 KCl, 1.25 NaH₂PO₄, 25 NaHCO₃, 20 glucose, 2 CaCl₂, and 1 MgCl₂.
- 3 Brains were dissected and sectioned on a Leica vibratome (VT1000S) to 300 μ m thickness.
- 4 Slices were transferred to a holding chamber and incubated at 34°C for 20 minutes prior to recordings.

Current clamp recordings

- 5 Slices were placed in the recording chamber at RT and continually perfused with oxygenated aCSF with 1mM scopolamine to block cholinergic interneuron firing.
- 6 The dorsal striatum was located under a 10x air objective, then switched to a 60x water immersion objective was used to target neurons for recording (Olympus, Tokyo, Japan). Slices were visualized using DIC and a QIClick microscope camera (QImaging, Surrey, Canada). GFP-positive cells were selected using a CoolLED pe4000 system (CoolLED Ltd., Andover, UK).
- 7 Patch electrodes were pulled from borosilicate capillary glass to a resistance of 2.5–6 M Ω . Electrodes were filled with an intracellular solution comprised of (in mM) 135 K-gluconate, 4 KCl, 10 HEPES, 10 Na-phosphocreatine, 4 MgATP, 0.4 Na₂GTP, 1 EGTA, 20 μ M Alexa 594 (pH 7.28, 298–305 mOsm/L).
- 8 After break-in, cells were held for ~5–10 minutes prior to current clamp recordings.
- 9 Recordings were acquired using a Multiclamp 700B amplifier (Axon Instruments, Union City, CA), digitized at a rate of 20 kHz, and filtered at a rate of 3–4 kHz.
- 10 Five second-long sweeps were acquired using a version of the MATLAB-based (MathWorks, Natick, MA) acquisition suite, ScanImage.



- 11 A series of 500 ms-long current steps from -50 to 300 pA were pseudorandomized. Current steps occurred one second after the start of the sweep, with a 20 ms-long, -20 pA current step added at the end of the sweep to test input resistance throughout the recording.
- 12 Cells with input resistance that changed >10% over the recording were discarded. Sweeps occurred with an inter-sweep interval of 10-20 seconds. A minimum of five sweeps for each current step magnitude was obtained.

Current clamp recording analysis

- 13 Average firing rate at each current step was determined using custom built MatLab code.
- 14 All values are reported as Hz. Firing rate values were then exported to GraphPad Prism (GraphPad, LaJolla, CA) for data visualization and statistics.