

Jan 09, 2025

CUT and RUN - Cas9 crosslinked



Forked from [CUT and RUN](#)



[Cell Genomics](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l6dr5rvqe/v1>

anita.adami¹

¹Laboratory of Molecular Neurogenetics, Department of Experimental Medical Science, Wallenberg Neuroscience Center and Lund Stem Cell Center, BMC A11, Lund University, 221 84 Lund, Sweden.



Anita Adami

Lund University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l6dr5rvqe/v1>

External link: <https://doi.org/10.1016/j.xgen.2025.100979>

Protocol Citation: anita.adami 2025. CUT and RUN - Cas9 crosslinked . [protocols.io](#)

<https://dx.doi.org/10.17504/protocols.io.bp2l6dr5rvqe/v1>

Manuscript citation:

Adami A, Garza R, Gerdes P, Johansson PA, Dorazehi F, Koutounidou S, Castilla-Vallmanya L, Atacho DA, Sharma Y, Johansson JG, Tam O, Kirkeby A, Barker RA, Hammell MG, Douse CH, Jakobsson J (2025) LINE-1 retrotransposons mediate cis-acting transcriptional control in human pluripotent stem cells and regulate early brain development. Cell Genomics 5(10). doi: [10.1016/j.xgen.2025.100979](https://doi.org/10.1016/j.xgen.2025.100979)

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 09, 2025

Last Modified: January 09, 2025

Protocol Integer ID: 117996

Keywords: ASAPCRN, cas9, fresh cell, protocol, cell

Funders Acknowledgements:

Aligning Science Across Parkinson's through the Michael J. Fox Foundation for Parkinson's Research

Grant ID: ASAP-000520

Swedish Research Council

Grant ID: 2018-02694

Swedish Brain Foundation

Grant ID: FO2019-0098

Cancerfonden

Grant ID: 190326

Barncancerfonden

Grant ID: PR2017-0053

NIHR Cambridge Biomedical Research Centre

Grant ID: NIHR203312

Swedish Society for Medical Research

Grant ID: S19-0100

National Institutes of Health

Grant ID: HG002385

Swedish Research Council

Grant ID: 2021-03494

Swedish Research Council

Grant ID: 2020-01660

Abstract

This protocols describe how to perform crosslink CUT&RUN on fresh cells



Sample preparation

25m

- 1 Detach cultured cells with Accutase (37 °C incubation for 5 minutes, or until cells completely detach) and resuspend in washing buffer (10% knockout serum [KSR] in iPSBrew).

NOTE: For optimal outcome, it is essential to use fresh cells.

- 2 Pellet cell by centrifuging them @  400 x g for  00:05:00 .

5m

- 3 Count cells and harvest 500 000 cells per condition (and antibody).

EXAMPLE: If you are running a CRISPRi experiment (CRISPRi gRNA hiPSCs & non-targeting control hiPSCs) and you want to perform a Cas9 CUT&RUN on both conditions, you will harvest 500 000 CRISPRi and 500 000 control hiPSCs. If you are also including CUT&RUN controls (e.g. IgG control), you will need 500 000 per condition, per antibody- in this case: 500 000+500 000 CRISPRi hiPSCs & 500 000+500 000 control hiPSCs.

CUT&RUN (crosslinked)

2h 23m

- 4 Activate ConA-coat magnetic beads (Epicypheer) by washing twice in Bead Binding Buffer [20 mM HEPES pH 7.5, 10 mM KCl, 1 mM CaCl₂, 1 mM MnCl₂]. Place on ice until use.

- 5 Centrifuge cells at  600 x g for  00:03:00 and resuspend the pellet in 1 ml of 0.1% formaldehyde solution to allow light crosslinking. Incubate for  00:01:00 at  Room temperature .

4m



NOTE: Crosslinking is performed to capture more challenging targets such as Cas9, but it is not needed for the CUT&RUN control samples (e.g. IgG). To perform a regular CUT&RUN, refer to the the standard CUT and RUN protocol:

dx.doi.org/10.17504/protocols.io.j8nlkwb8dl5r/v1

CAUTION: Formaldehyde is classified as a carcinogen and can cause respiratory and skin irritation upon exposure. Handle under a chemical hood and with the relevant protection equipment.

- 6 Add 125 mM of glycine and centrifuge again at  600 x g for  00:03:00 .

3m



- 7 Remove supernatant, and resuspend cells in 800 μL of XL-Wash Buffer (20 mM HEPES pH 7.5; 150 mM NaCl; 0.5 mM spermidine; 1 Roche Complete Protease Inhibitor EDTA-free tablet; 1% Triton X-100; 0.05% SDS - all diluted in dH_2O).
- 8 Centrifuge cells at 1000 x g for 00:03:00 and resuspend again in XL-Wash Buffer. Repeat two more centrifugations at 1200 x g (3 minutes each). 3m
- 9 Bind nuclei to beads (10 μL) for 00:10:00 at RT with gentle rotation, and then split bead-bound nuclei into three equal volumes (here, corresponding to IgG control and Cas9 antibody treatment). 10m
- 10 Remove wash buffer and resuspend nuclei in 100 μL cold XL-Antibody Buffer (XL-Wash Buffer; 0.05% digitonin; 2 mM EDTA; 1:50 anti-Cas9 antibody).
- CAUTION:** Digitonin is toxic. Gloves should be worn when preparing and working with digitonin.
- 11 Incubate Overnight at 4 $^{\circ}\text{C}$ on gentle shake. 3m
- 12 On the following day, remove supernatant by placing the sample(s) on a magnetic stand, and resuspend/wash the beads twice in XL-Digitonin Buffer (XL-Wash Buffer + 0.05% digitonin).
- 13 After the final wash, add pA-MNase (final concentration 700 ng/ml) in XL-Digitonin Buffer and incubate with the cells at 4 $^{\circ}\text{C}$ for 01:00:00 . Wash nuclei twice, resuspend in 100 μL XL-Digitonin Buffer, and chill to 0 $^{\circ}\text{C}$ - 2 $^{\circ}\text{C}$ in a metal block sitting on wet ice. 1h
- 14 Stimulate genome cleavage by addition of 2 mM CaCl_2 at 0 $^{\circ}\text{C}$ for 00:30:00 . Quench the reaction by adding 100 μL 2x Stop Buffer (0.35 M NaCl; 20 mM EDTA; 4 mM EGTA; 0.05% digitonin; 50 ng/ml of RNase A; 50 ng/ml of glycogen; 10 ng/ml spike in DNA; all diluted in distilled H_2O) and vortex. 30m
- 15 Incubate at 37 $^{\circ}\text{C}$ for 00:30:00 to allow the release of the fragmented chromatin. 30m



- 16 Place bead-bound nuclei on the magnet stand and collect the supernatant in a new 1.5 ml Eppendorf tube.
- 17 Add 0.09% SDS and 0.22 mg/ml of Proteinase K to reverse the crosslinking.
- 18 Isolate the DNA using NucleuSping PCR Cleanup columns (Macherey-Nagel) and prepare libraries for sequencing.

Sequencing

- 19 Prepare Illumina sequencing libraries using the Hyperprep kit (KAPA) with unique dual-indexed adapters (KAPA).