

Mar 15, 2024

CUT and RUN

DOI

dx.doi.org/10.17504/protocols.io.j8nlkwb8dl5r/v1

anita.adami¹

¹Laboratory of Molecular Neurogenetics, Department of Experimental Medical Science, Wallenberg Neuroscience Center and Lund Stem Cell Center, BMC A11, Lund University, 221 84 Lund, Sweden.



Raquel Garza

Lund University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlkwb8dl5r/v1>

Protocol Citation: anita.adami 2024. CUT and RUN. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.j8nlkwb8dl5r/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 29, 2022

Last Modified: May 31, 2024

Protocol Integer ID: 69297

Keywords: ASAPCRN, human brain tissue, brain tissue, cut, tissue, brain

Funders Acknowledgements:

Aligning Science Across Parkinson's through the Michael J. Fox Foundation for Parkinson's Research

Grant ID: ASAP-000520

Swedish Research Council

Grant ID: 2018-02694

Swedish Brain Foundation

Grant ID: FO2019-0098

Cancerfonden

Grant ID: 190326

Barncancerfonden

Grant ID: PR2017-0053

NIHR Cambridge Biomedical Research Centre

Grant ID: NIHR203312

Swedish Society for Medical Research

Grant ID: S19-0100

National Institutes of Health

Grant ID: HG002385

Swedish Research Council

Grant ID: 2021-03494

Swedish Research Council

Grant ID: 2020-01660




Abstract

This protocols describe how to perform CUT&RUN on human brain tissue (frozen)









Troubleshooting



Sample extraction

- 1 Flash-freeze postmortem brain tissue
- 2 Sample  50 mg -  100 mg from human brain tissue and store at  -80 °C until use

CUT&RUN

- 3 Activate ConA-coat magnetic beads (Epicypther) by washing twice in bead binding buffer [20 mM HEPES pH 7.5, 10 mM KCl, 1 mM CaCl₂, 1 mM MnCl₂]. Place on ice until use.
- 4 **Isolate nuclei from frozen tissue** after incubating with Recombinant Alexa Fluor® 488 Anti-NeuN antibody [EPR12763] - Neuronal Marker (ab190195) at a concentration of 1:500 for 30 minutes on ice.
Run nuclei through the FACS at  4 °C with low flowrate using a 100 mm nozzle and isolate 300.000 nuclei Alexa Fluor – 488 positive nuclei.
- 5 Pellet the sorted nuclei at 1,300 x g for  00:15:00 and resuspend in  1 mL of ice-cold nuclear wash buffer (20 mM HEPES, 150 mM NaCl, 0.5 mM spermidine, 1x cOmplete protease inhibitors, 0.1% BSA) and  10 µL per antibody treatment of ConA-coated magnetic beads (Epicypther) added with gentle vortexing (Pipette tips for transferring nuclei were pre-coated with 1% BSA). 15m
- 6 Bind nuclei to beads for  00:10:00 at RT with gentle rotation, and then split bead-bound nuclei into three equal volumes (corresponding to IgG control, H3K4me3 and H3K9me3 treatments). 10m
- 7 Remove wash buffer and resuspend nuclei in  100 µL cold nuclear antibody buffer (20 mM HEPES pH 7.5, 0.15 M NaCl, 0.5 mM Spermidine, 1x Roche complete protease inhibitors, 0.02% w/v digitonin, 0.1% BSA, 2 mM EDTA) containing primary antibody at 1:50 dilution and incubate at  4 °C  Overnight with gentle shaking. 10m
- 8 Wash nuclei thoroughly with nuclear digitonin wash buffer (20 mM HEPES pH 7.5, 150 mM NaCl, 0.5 mM Spermidine, 1x Roche cOmplete protease inhibitors, 0.02% digitonin, 0.1%



BSA) on the magnetic stand.

- 9 After the final wash, add pA-MNase in nuclear digitonin wash buffer and incubate with the nuclei at 4 °C for 01:00:00 . Wash nuclei twice, resuspend in 100 µL digitonin buffer, and chill to 0 °C - 2 °C in a metal block sitting in wet ice. 1h
- 10 Stimulate genome cleavage by addition of 2 mM CaCl₂ at 0 °C for 30 min. Quench the reaction by adding 100 µL 2x stop buffer (0.35 M NaCl, 20 mM EDTA, 4 mM EGTA, 0.02% digitonin, 50 ng/µL glycogen, 50 ng/µL RNase A, 10 fg/µL yeast spike-in DNA) and vortex.
- 11 Incubate 00:30:00 at 37 °C to release genomic fragments. Place bead-bound nuclei on the magnet stand and purify fragments from the supernatant using a NucleoSpin clean-up kit (Macherey-Bagel). 30m

Sequencing

- 12 Prepare Illumina sequencing libraries using the Hyperprep kit (KAPA) with unique dual-indexed adapters (KAPA).