

Mar 28, 2025

CycloneSEQ library construction from DNA of isolated bacteria

In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vzk3k2vx1/v1>

Wang Mengmeng¹, Hewei Liang¹, Tao Zeng¹, Junyi Chen¹, Yuanqiang Zou¹, Liang Xiao¹

¹BGI

GigaScience Press

BGI



Mengmeng Wang

BGI



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzk3k2vx1/v1>

Protocol Citation: Wang Mengmeng, Hewei Liang, Tao Zeng, Junyi Chen, Yuanqiang Zou, Liang Xiao 2025. CycloneSEQ library construction from DNA of isolated bacteria. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.rm7vzk3k2vx1/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 21, 2024

Last Modified: March 28, 2025

Protocol Integer ID: 112526

Keywords: genomics, long-read sequencing, microbiology, CycloneSEQ, cycloneseq library construction from dna, isolated bacteria cycloneseq, cycloneseq library construction, cycloneseq library protocol, research this cycloneseq library protocol, cycloneseq platform, assembly of circular complete genome, circular complete genome, read sequencing, sequencing platform, genome, bacterial strain, mixed microbial community, bacteria, dna

Funders Acknowledgements:

The Shenzhen Municipal Government of China

Grant ID: XMHT20220104017, CXB201108250097A

Abstract

CycloneSEQ is a novel long-read sequencing platform developed by BGI-Research This CycloneSEQ library protocol is tailored for the assembly of circular complete genomes from mixed microbial communities. DNA is extracted from bacterial strains and after construction the resulting libraries can be sequenced on the CycloneSEQ platform.

Troubleshooting

CycloneSEQ library construction from DNA of isolated bacteria

50m

- 1 Approximately  2 mL of liquid medium with 10^8 or 10^9 CFU (colony forming units) of bacterial cells are extracted using a DNA extraction Kit for high-throughput applications.
- 2 CycloneSEQ library preparation and sequencing follows the manufacturer's guidelines. Each sample, containing 2 μg of input DNA (≥ 21 ng/ μL), is diluted with nuclease-free water to  192 μL , then mixed with  14 μL each of DNA repair buffers 1 and 2,  12 μL of DNA repair enzyme 1, and  8 μL μL of DNA repair enzyme 2.
- 3 The mixtures are then incubated in a thermocycler at  20 $^{\circ}\text{C}$ for  00:10:00,  65 $^{\circ}\text{C}$ for  00:10:00, and held at  4 $^{\circ}\text{C}$.
- 4 After incubation, the mixtures are purified with 1.0x DNA clean beads and eluted with  240 μL of nuclease-free water.
- 5 The end-repaired samples are then mixed with  10 μL of sequencing adaptors,  100 μL of 4x ligation buffer,  40 μL of DNA ligase, and  10 μL of nuclease-free water, and incubated at  25 $^{\circ}\text{C}$ for  00:30:00.
- 6 The ligated products are purified again with 1.0x DNA clean beads, resuspended with long fragment wash buffer, and recovered into  42 μL of elution buffer.
- 7 The resulting libraries are quantified using a Qubit fluorometer and sequenced on the CycloneSEQ WuTong02 platform according to the manufacturers sequencing protocols.

20m

30m