

Jun 12, 2024

Version 1

CYToF Staining V.1

 [Nature Communications](#)

DOI

dx.doi.org/10.17504/protocols.io.14egn69nml5d/v1

Meelad Amouzgar¹, Patricia Favaro¹, Daniel Ho¹, Trevor Bruce¹, Kausalia Vijayaragavan¹, Sean Bendall¹

¹Stanford University



Kausalia Vijayaragavan

Stanford University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn69nml5d/v1>

External link: <https://doi.org/10.1038/s41467-025-63883-4>

Protocol Citation: Meelad Amouzgar, Patricia Favaro, Daniel Ho, Trevor Bruce, Kausalia Vijayaragavan, Sean Bendall 2024. CYToF Staining. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.14egn69nml5d/v1>

Manuscript citation:

Amouzgar M, Favaro P, Ho D, Bruce T, Bendall SC (2025) A deep single cell mass cytometry approach to capture canonical and noncanonical cell cycle states. Nature Communications 16(). doi: [10.1038/s41467-025-63883-4](https://doi.org/10.1038/s41467-025-63883-4)



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 12, 2024

Last Modified: June 12, 2024

Protocol Integer ID: 101699

Keywords: flow cytometry, optical flow cytometry, mass cytometry, staining cytof, simultaneous analysis of diverse immune cell population, cell immunological assay, cytof, immunological assay, comprehensive characterization of immune population, diverse immune cell population, immune population, fluorescent label

Disclaimer

none

Abstract

CYToF has significantly advanced immunophenotyping, especially in exploratory settings where comprehensive characterization of immune populations is necessary but sample sizes are limited. It is a powerful tool for single-cell immunological assays, particularly for the detailed and simultaneous analysis of diverse immune cell populations.

By using metal isotopes instead of fluorescent labels, mass cytometry eliminates the spectral overlap issues seen in optical flow cytometry, allowing robust analysis of around 60 individual parameters simultaneously. The technique involves increased complexity in the design, execution, and interpretation of experiments.

Troubleshooting




Reagents

- 1
 - Paraformaldehyde (PFA) 16% (Electron Microscopy Sciences – catalog number 15710)
 - Cisplatin Cell-ID™ Cisplatin-198Pt—100 µL (#201198standard biotools)
 - 1x Low-Barium PBS
 - BSA
 - Sodium Azide
 - Cell Suspension Media (CSM)

Reagent Preparation

2 Cisplatin-198Pt Preparation

100mM stock Cisplatin-198Pt is stored at  -20 °C

Prepare 0.5uM working stock solution of Cisplatin-198Pt by carrying out two dilution

1. add 1ul of 100mM Cisplatin-198Pt stock + 199ul of Low-Barium PBS = 0.5mM (500uM) Cisplatin-198Pt
2. add 1ul of 0.5mM (500uM) Cisplatin-198Pt + 999ul of Low-Barium PBS = 0.005mM (0.5uM) Cisplatin-198Pt

For final live/dead cell staining with PBMCs is 0.025uM Cisplatin-198Pt/1million cell suspension.

Cisplatin concentration used does not reach saturation. So 0.5uM is plenty for 1-10 million cells. Cisplatin staining intensity increases with cell size (monocytes population is greater than lymphocyte population), neurons greater than astrocytes than microglial. Cisplatin binds covalently to any proteins in the dead/dying cells.

	A	B	C	D	E
	Stock (mM)	dilution 1	Working Conc (uM)	dilution 2	Final (uM)
	100	0.005	500	0.001	0.5

Cisplatin Concentrations

3 Cell Suspension Media (CSM)

1 L of 1x Low-Barium PBS
5g of BSA
200 mg of sodium azide
Mix ingredients in appropriate container



Filter using 2× 500 mL vacuum filter unit

CyToF Staining

4 **PFA Fixation:**

1. Add 100 μ L of 16% PFA per 1 mL of media and incubate cells for 10 minutes at room temperature.
2. Wash with 4mL of CSM 2x
3. Pellet the cells by centrifugation at 600xg for 5 min.
4. Store fixed cells in 1mL of CSM at -80C

5 **Barcoding protocol:**

1. Transfer cells to cluster tube (if not already).
2. Wash 1x with cold PBS, spin 5' 600G 4C, aspirate, vortex to re-suspend pellet.
3. Wash 1x with 1mL cold PBS-S, spin 5' 600G 4C, aspirate, vortex to re-suspend pellet.
4. Cut barcode plate cover to expose only the barcode wells to be used (keep cover on until ready to pipette)

The next 3 steps must be completed quickly to prevent weak/uneven staining of the barcoding metals

1. Using a multichannel pipette, remove 100 μ L **cold** PBS-S from reservoir, add to the corresponding barcode wells on the barcoding plate & pipette up and down several times.
2. Transfer the 100 μ L of now PBS-S + Barcode to the corresponding sample tubes, mix up & down several times (final sample + barcode volume ~180 μ L) (we skip the step where we add an additional 100 μ L PBS to make 280 μ L)
3. Use a multichannel pipette to repeat for all samples to be barcoded, then let incubate 15' at RT in dark.
4. After 15 minutes, add 700 μ L CSM to quench, Spin 5' 600G 4C, aspirate, then wash 1X with 1 mL CSM.
5. Aspirate & re-suspend cell pellets by vortex. For each cluster tube, add 100 μ L CSM, then adjust pipette volume to collect the entire volume (~160 μ L), and deposit into a FACS tube (use larger tube if total cell number in pooled sample is large). Repeat for each sample to pool all barcoded cells into the FACS tube for staining.

6 **Surface staining**

100 μ L total staining.

1. Bring volume to 55 μ L.
2. Add 5 μ L of FC block (**human**) and incubate for 10' at RT in the shaker/
3. Make surface antibody cocktail (~40 μ L total) and add to cells (100 μ L total).
4. Incubate for 30' at RT in the shaker.
5. Wash cells with 4 mL of CSM 600xg for 5 minutes.



7 Permeabilization with Methanol

1. Add 1mL of 4oC 100% methanol (keep on ice) (slowly)
2. Incubate for 10' on ice.
3. Add 3mL of CSM.
4. Pellet at 600xg for 5 min.
5. Wash more 2 times. Vortex the pellet and add the panel to a final volume of 100ul (1x test)

8 Intracellular staining

1. **Filter the intracellular antibody cocktails before adding to the cells.**
2. **10,000g 1 minute**
3. **prewet with csm (if volume is not too small. 20 ul is very little)**
4. Add the intracellular staining
5. Incubate for 30' at RT.
6. Add 3mL of CSM.
7. Pellet at 600xg for 5 min.

9 DNA intercalator

Iridium Intercalator: DNA1/DNA2 = Ir191/Ir192

For preparing Iridium DNA intercalator:

1. Prepare 5x the DNA intercalator solution:

DNA intercalator: 0.25μL of DNA intercalator (stored at 4°C) + 1mL of PBS +		
100uL of PFA 16%	+5 ul WGA	

Add 1ml of DNA intercalator solution to each tube and mix and incubate either 00:30:00 at Room temperature to run at cytof at same day or **overnight at 4°C**.

1. You can also store in fridge for ~ a week in the intercalation solution

10 Day of CyToF run

1. wash 1x CSM and 2x Water and run on CyTOF. 600xg/5min
1. Dilute EQ beads 1:10 in ddH₂O and add 1 mL of diluted beads per 1 million cells.
Usually, we add 1.5ml of water/beads per tube and run 500ul in one push at cytof. It is



enough to visualize the titration. If you thing you need more, run another push and concatenate the data afterwards.

OR

Long-term Storage

1. you Wash with 3ml CMS and spin down 600g/5min ,
2. resuspend in CSM+10% DMSO, and put right the way in -80oC.
3. You can keep it there for months if needed. Once ready, thaw and wash 1 time with CSM and 2 times with ddH2O and then acquire, 600g/5min

11