



May 08, 2025

Cytokine Array

 In 2 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.kxygx75yol8j/v1>

Lucas Blasco-Agell^{1,2}, Meritxell Pons-Espinal^{1,2}, Antonella Consiglio^{1,2}, Miquel Vila³

¹Department of Pathology and Experimental Therapeutics, Bellvitge University Hospital- IDIBELL, 08908 Hospitalet de Llobregat, Spain;

²Institute of Biomedicine of the University of Barcelona (IBUB), Carrer Baldori Reixac 15-21, Barcelona 08028, Spain.;

³VHIR-CIBERNED-ASAP

Vilalab Public



Miquel Vila

VHIR-CIBERNED-ASAP

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygx75yol8j/v1>

Protocol Citation: Lucas Blasco-Agell, Meritxell Pons-Espinal, Antonella Consiglio, Miquel Vila 2025. Cytokine Array.
protocols.io <https://dx.doi.org/10.17504/protocols.io.kxygx75yol8j/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 17, 2025

Last Modified: May 08, 2025

Protocol Integer ID: 130932

Keywords: immunology multiplex assay, cytokine array, based multiplex assay, multiplex assay, infy with bead, infy

Funders Acknowledgements:

ASAP

Grant ID: ASAP-020505

Abstract

The EMD MILLIPLEX® MAP Human Cytokine/Chemokine/Growth Factor Pane A – Immunology Multiplex Assay (Merck Millipore) is employed to simultaneously analyse TNF- α , IL-1 β , IL 6, and INF γ with Bead-Based Multiplex Assay using the Luminex® technology.

- 1 Culture supernatants from non-stimulated and 24-hour LPS-stimulated hMG are collected and stored at -80 °C for long storage.
- 2 Quality Controls (QC) and Human Cytokine Standard (STD) are prepared with deionized water by serial dilutions as manufacturer's instructions. Assay Buffer is used for Background or STD 0.
- 3 Wells are filled with 200 µL of Wash Buffer (WB), sealed and mixed at RT for 10 minutes on a Thermomixer Comfort plate shaker (EppendorfTM) at 500-800 rpm.
- 4 Plate was firmly tap upside down on absorbent paper and 25 µL of each STD and QC with 25 µL of RPMI medium in duplicates. 25 µL of sample with 25 µL of Assay Buffer are added in duplicates. Finally, 25 µL of the mixed antibody-bead preparation is added to every well and incubated at 2-8°C ON in a plate shaker.
- 5 The day after, well contents are gently removed by firmly attaching plate to a handheld magnet and washed with WB for 3 times.
- 6 Then, 25 µL of Detection Antibodies are added into each well and incubated at RT for 1 hour on a plate shaker.
- 7 After that, 25 µL of Streptavidin-Phycoerythrin are added into each well and incubated at RT for 30 minutes on a plate shaker. Well contents are gently removed with a handheld magnet and washed with WB for 3 times.
- 8 Antibody-beads are resuspended by adding 150 µL of Sheath Fluid to all wells and incubated at RT for 5 minutes on a plate shaker.
- 9 Plate is run on a multiplex system Luminex[®] 200TM (InvitrogenTM). Probe heights of every antibody are adjusted according to Luminex[®] recommended protocols employing the xPONENT[®] software (Luminex[®]).
- 10 Mean Fluorescent Intensity of every analyte is analyzed using a 5-parameter logistic or spline curve-fitting method for calculating analyte concentrations in samples.