

Feb 29, 2024

Cytotoxicity Assay Protocol

DOI

dx.doi.org/10.17504/protocols.io.14egn35mpl5d/v1

Roshni Jaffery^{1,2}, Ningbo Zheng^{1,2}, Jiakai Hou^{1,2}, Ashley Guerrero^{1,2}, Si Chen^{1,2}, Chunyu Xu^{1,2},
Nicholas A. Egan^{1,2}, Ritu Bohat^{1,2}, Weiyi Peng^{1,2}

¹Department of Biology and Biochemistry, University of Houston, Houston, TX, USA;

²Aligning Science Across Parkinson's Collaborative Research Network



Sarfraz Ahmed

Mass General Hospital

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn35mpl5d/v1>

Protocol Citation: Roshni Jaffery, Ningbo Zheng, Jiakai Hou, Ashley Guerrero, Si Chen, Chunyu Xu, Nicholas A. Egan, Ritu Bohat, Weiyi Peng 2024. Cytotoxicity Assay Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.14egn35mpl5d/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: February 29, 2024

Last Modified: February 29, 2024

Protocol Integer ID: 95903

Keywords: cytotoxicity assay, cancer cell lines, flow cytometry, cytotoxicity assay protocol, protocol for cytotoxicity, using flow cytometry, cytotoxicity, cancer cell line, protocol, assay, cell line

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-000312

Abstract

This is the protocol for cytotoxicity assay in cancer cell lines using flow cytometry.

Materials

1. Tumor cell line
2. Tumor-reactive T cells
3. Tumor cell medium: RPMI 1640 (Gibco 11875) + 10% Fetal Bovine Serum (FBS) complete mixture
4. T cell medium: 50ml complete MEM media (Gibco 11095080) + 500 IU/ml IL2 + add 50 μ L of 1000X beta-mercaptoethanol (55mM)
5. 1X Phosphate-buffered saline (PBS)
6. FACS Buffer: 2% FBS in PBS
7. APC-conjugated rabbit anti-active caspase 3 Ab
8. BD Cytotfix/CytopermTM Fixation/Permeabilization Kit
9. Fixation/permeabilization solution
10. BD Perm/WashTM Buffer

Troubleshooting



- 1 Begin T cell culturing a week prior. Grow and culture tumor cell lines. Ensure the cells are not over-confluent.
- 2 Split the tumor cells the day before the assay. Perform the cytotoxic assay as described in next steps.
- 3 Count tumor cells. Collect 3-5M cells for each cell line.
- 4 Centrifuge the cells and wash the pellet once with PBS. Ensure that all media is washed off as this can inhibit the staining.
- 5 Aspirate the supernatant and leave as little PBS as possible
- 6 Re-suspend cells in culture medium at the appropriate ratio for the experiment (1×10^6 / ml for 50,000 cells/50ul/well or about ~80% confluent) using 96 well round well plate with 50,000 cell tumor and in triplication.
- 7 Add 50 ul of T cell suspension diluted to achieve the desired E:T ratio. Below is an example of how to set up E:T ratios:

	A	B	C
	ET Ratio(s)	T cells/50ul /well	Tumor cells/50ul/well
	0:1	0 (only with T cell medium)	5×10^4 cells
	1:1	5×10^4 cells	5×10^4 cells
	3:1	1.5×10^5 cells	5×10^4 cells

- 8 Centrifuge the mixed cells in the tubes at 300 rpm for 5 minutes to gently pellet the cells. Incubate for 3 hours at 37°C .
- 9 Add 150ul PBS into the plate and then centrifuge at 2000rpm for 5min.



- 10 Remove the supernatant and add then 100ul/well Trypsin. The plate is incubated at 37C incubator for 3-5mins.
- 11 Add 100ul FACS buffer to stop the trypsinization and centrifuge at 2000rpm for 5 min.
- 12 Remove the supernatant carefully.
- 13 Fix and permeabilize cells using Cytofix/Cytoperm™ Kit (BD BioSciences) for 20 minutes at room temperature (100ul/well – **mix well**).
- 14 Wash with Perm/Wash™ Buffer (150ul/well) at 2000 rpm for 5 min.
- 15 Prepare caspase antibody dilution such that each well gets stained with 2.5ul of APC-conjugated rabbit anti-active caspase 3 Ab in 50 ul Perm/Wash™ Buffer. Mix gently and incubate for 30 min in the dark on ice or in the 4°C refrigerator.
- 16 Wash cells 1Xwith Perm/Wash™ Buffer (150ul).
- 17 Resuspend the cell pallets into 100ul FACS buffer (2% FBS in 1X PBS).
- 18 Analyze the samples using flow cytometry.