

Aug 08, 2020 Version 2

Cytotoxicity assay using LLC-MK2 cell line V.2

 [PLOS One](#)

DOI

dx.doi.org/10.17504/protocols.io.bjgbkjsn

Iuanaborba¹

¹Universidade Federal do Rio de Janeiro

L. P. Borba-Santos



Iuanaborba

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bjgbkjsn

External link: <https://doi.org/10.1371/journal.pone.0240658>

Document Citation: Iuanaborba 2020. Cytotoxicity assay using LLC-MK2 cell line. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bjgbkjsn>

Manuscript citation:

Borba-Santos LP, Vila T, Rozental S (2020) Identification of two potential inhibitors of *Sporothrix brasiliensis* and *Sporothrix schenckii* in the Pathogen Box collection. PLoS ONE 15(10): e0240658. doi: [10.1371/journal.pone.0240658](https://doi.org/10.1371/journal.pone.0240658)

License: This is an open access document distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Created: August 07, 2020

Last Modified: August 08, 2020

Document Integer ID: 40163

Abstract

Protocol that described the cytotoxicity assay using a mammalian epithelial cell line (LLC-MK2; ATCC CCL-7).

Cytotoxicity assay using LLC-MK₂ cell line

1. The mammalian epithelial cell line LLC-MK₂ (ATCC CCL-7) was cultivated in RPMI 1640 medium¹ supplemented with 2 mM L-glutamine and heat-inactivated 10% fetal bovine serum and buffered with sodium bicarbonate at pH 7.5;
2. Aliquots of 100 µl containing 2.5×10^5 cells/ml were added into flat-bottom 96-well microplates (TPP™) and incubated for 48 h at 37 °C, in a 5% CO₂ atmosphere, to obtain confluent monolayers of LLC-MK₂ cells;
3. The supernatant was gently discarded, samples were washed in sterile PBS, and 100 µl of compounds previously diluted² in RPMI 1640 medium were added on confluent monolayers;
4. Cells were treated for 48 h at 37 °C, in a 5% CO₂ atmosphere;
5. The supernatant was gently discarded and samples were washed in sterile PBS;
6. Cell viability was analyzed by the tetrazolium (XTT) reduction assay and 150 µl of XTT solution (1 mg/ml XTT¹ and 1mM menadiona¹ in PBS) were added in each well;
6. Microplates were incubated for 2 h at 37°C, in a 5% CO₂ atmosphere, in the dark;
7. Microplates were centrifugated at 4000 rpm for 5 min and the supernatant was added into a new microplate;
8. Spectrophotometric readings at 492 nm were performed using a microtiter plate reader³;
9. The absorbance value for each well was subtracted from the value for the negative controls⁴ and inhibition of cell growth (I) relative to positive controls⁵ was calculated according to the following equation: $I = 100 - (A \times 100/C)$, where A is the absorbance of treated wells, and C is the absorbance of untreated control wells;
10. The concentration of compounds that elicited 50% cytotoxicity (CC₅₀) was estimated by linear regression;
11. The diluent control containing 1% DMSO was included in experiments;
12. Experiments were performed in triplicate in two independent moments.

¹Sigma Chemical Co., USA.

² Stock solutions of compounds in dimethyl sulfoxide (DMSO) at 1 mM were diluted in RPMI 1640 medium supplemented with 2 mM L-glutamine to obtain concentrations of 0.1, 1, 2, 4, 5, 7, 8, 9, and 10 μ M.

³EMax Plus, Molecular Devices, USA.

⁴ Wells without cells containing only RPMI media that were incubated in the same conditions of wells with cells.

⁵ Wells with untreated cells.