



Aug 22, 2022

DAB immunostaining of thin, fixed mouse brain tissue sections using HNA or NCAM to characterize human iPSC-derived cell xenografts

DOI

dx.doi.org/10.17504/protocols.io.14egn7pb6v5d/v1

Benjamin Trist¹, Ashish Mathai¹, Asheeta Prasad¹

¹The University of Sydney



Benjamin Trist

The University of Sydney

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn7pb6v5d/v1>

Protocol Citation: Benjamin Trist, Ashish Mathai, Asheeta Prasad 2022. DAB immunostaining of thin, fixed mouse brain tissue sections using HNA or NCAM to characterize human iPSC-derived cell xenografts. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.14egn7pb6v5d/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 12, 2022

Last Modified: May 31, 2024

Protocol Integer ID: 68551

Keywords: NCAM, HNA, Human-to-mouse xenograft, Human iPSC, Immunohistochemistry, ASAPCRN, mouse brain tissue section series, fixed mouse brain tissue section series, mouse brain tissue section, brain of athymic mice, fixed mouse brain tissue section, growth of human ipsc, derived cell, immunohistochemistry, derived cell xenograft, living brain, human ipsc, athymic mice, cell xenografts this protocol, mice

Funders Acknowledgements:

Michael J Fox Foundation

Grant ID: ASAP-000497

Abstract

This protocol describes our use of chromogenic 3,3'-diaminobenzidine (DAB) immunohistochemistry to identify human iPSC-derived cells within thin, fixed mouse brain tissue section series'. We apply this workflow for post-mortem assessment of the survival and growth of human iPSC-derived cells which have been transplanted into the living brain of athymic mice.

Attachments



[it75bj7ap.docx](#)

100KB






Materials



Equipment:

- Horizontal rocker
- Vortex
- Microcentrifuge
- Glass petri dish
- Oven

Consumables:

-  20 mL scintillation vials
- Paint brushes
- Gelatin-Chrom Alum-coating microscope slides
 1. See related protocol - Coating superfrost microscope slides with gelatin-chromium potassium sulfate
- Microscope slide coverslips (no. 1.5,  25 mm x  75 mm)
- Glass pipettes
- Rubber teats
- Transfer pipettes

Key reagents:

- Optimal Cutting Temperature (OCT) compound
- Bovine Serum Albumin (BSA)
- Casein
- Sodium citrate
- Tween-20 and Triton X-100
- Ethanol
- Hydrogen peroxide (H₂O₂)
- DEPEX
- 3, 3'-diaminobenzidine (Sigma #D5905)
- Antibodies
 1.  NCAM **Abcam Catalog #ab75813**
 2. HNA (Novus #NOVNBP-313912)
 3. Biotinylated anti-rabbit secondary antibody (Vector Labs #BA-1000)
 4.  Avidin/Biotin HRP Complex **Vector Laboratories Catalog #PK-6100**

Solutions:

- 1x PBS, pH 7.4

A	B
Antigen retrieval (AR) buffer	
Sodium citrate	2.94 g (10 mM)
Tween-20	500 µL (0.05%)
Up to 1L with dH ₂ O, pH 6.0	

- Quenching solution
10mL (3.3%) 33% H₂O₂, 50mL (50%) ethanol up to 100mL with 1x PBS
- 1x PBST
500 µL (0.05%) Tween-20 in 1L 1x PBS

A	B
Blocking solution	
Casein	1 g (1% w/v)
Triton X-100	250 µL (0.25% v/v)
Glycine	1.5 g (1.5% w/v)
BSA	5 g (5% w/v)
Up to 100mL with 1xPBS	







Material input (animal, cell, tissue, fraction details):

Thin, fixed athymic mouse brain tissue sections prepared from whole mouse brains grafted with human iPSC-derived neural progenitor cells.





Troubleshooting



Day 1 (~4-6 hrs)

- 1 Pre-heat oven and AR buffer to 70 °C .
- 2 Label scintillation vials to match labels on section storage plates (mouse and section series IDs, name, date etc.).
- 3 Transfer sections into scintillation vials using a transfer pipette or fine paintbrush.
- 4 Remove anti-freeze solution and perform 3× 7 min washes in 1x PBS at Room temperature with gentle agitation.   
 - Slow shaking on an orbital rocker recommended for washes/incubations to ensure even contact with solutions.
 - Use a glass pipette and rubber teat to remove solution during wash changes.
 - Anti-freeze solution must be rinsed off prior to immunostaining.
- 4.1 Remove anti-freeze solution and wash in 1x PBS at Room temperature for 00:07:00 (1/3). 
- 4.2 Remove anti-freeze solution and wash in 1x PBS at Room temperature for 00:07:00 (2/3). 
- 4.3 Remove anti-freeze solution and wash in 1x PBS at Room temperature for 00:07:00 (3/3). 

Antigen retrieval (AR)

- 5 Incubate sections in AR buffer for 00:30:00 at 70 °C . 

- 6 Preheat AR buffer to 70 °C prior to addition to sections.
- 7 After antigen retrieval, allow sections to cool for 00:30:00 before proceeding with staining. 
- 8 Perform 3× 7 min washes in 1x PBST with agitation. 



8.1 Wash in 1x PBST with agitation for ⌚ 00:07:00 (1/3).

7m



8.2 Wash in 1x PBST with agitation for ⌚ 00:07:00 (2/3).

7m



8.3 Wash in 1x PBST with agitation for ⌚ 00:07:00 (3/3).

7m



Quenching step

9 Incubate sections in quenching solution for ⌚ 00:30:00 at 🌡 Room temperature with gentle agitation.

30m



10 Perform 3× 7 min washes in 1x PBST with gentle agitation.



10.1 Wash in 1x PBST with gentle agitation for ⌚ 00:07:00 (1/3).

7m



10.2 Wash in 1x PBST with gentle agitation for ⌚ 00:07:00 (2/3).

7m



10.3 Wash in 1x PBST with gentle agitation for ⌚ 00:07:00 (3/3).

7m



Blocking step

11 Incubate sections in blocking solution for ⌚ 01:00:00 at 🌡 Room temperature with gentle agitation.

1h



Primary antibody step



- 12 Incubate sections with NCAM (1:20,000) or HNA (1:15,000) primary antibodies diluted in blocking buffer Overnight at 4 °C with gentle agitation.



Day 2 (~8hrs)

- 13 Perform 3× 7 min washes in 1x PBST with gentle agitation.



- 13.1 Wash in 1x PBST with gentle agitation for 00:07:00 (1/3).

7m



- 13.2 Wash in 1x PBST with gentle agitation for 00:07:00 (2/3).

7m



- 13.3 Wash in 1x PBST with gentle agitation for 00:07:00 (3/3).

7m



Secondary antibody step

- 14 Incubate sections in anti-rabbit biotinylated secondary antibody 1:500 diluted in blocking buffer that has been diluted 2-fold for 02:00:00 at Room temperature .

2h



- 15 Perform 3× 7 min washes in 1x PBST with gentle agitation.



- 15.1 Wash in 1x PBST with gentle agitation for 00:07:00 (1/3).

7m



- 15.2 Wash in 1x PBST with gentle agitation for 00:07:00 (2/3).

7m





- 15.3 Wash in 1x PBST with gentle agitation for 00:07:00 (3/3).

7m






Tertiary complex step

16 Incubate sections in Avidin-Biotin Complex (ABC) kit solution (Vector Laboratories) for  02:00:00 at  Room temperature .

2h



17 Prepare tertiary complex  00:30:00 prior to use according to the manufacturer's instructions.

30m

-  100 μ L A +  100 μ L B +  9800 μ L 1x PBS (1:100 A + 1:100 B in 1x PBS).

18 Perform 3 \times 7 min washes in 1x PBST with gentle agitation.



18.1 Wash in 1x PBST with gentle agitation for  00:07:00 (1/3).

7m



18.2 Wash in 1x PBST with gentle agitation for  00:07:00 (2/3).

7m





18.3 Wash in 1x PBST with gentle agitation for  00:07:00 (3/3).

7m



Chromogen step

19 Perform DAB staining as follows;



20 Prepare DAB solution by dissolving  10 mg DAB tablet into  20 mL PBS ([M] 0.5 mg/mL).

21 Filter through Whatman paper #1 or 0.22 μ m syringe filter before use.


Note

NB: DAB is a suspected mutagen and should be handled with care.



22 Prepare DAB-H₂O₂ solution **immediately** prior to use by adding  10 μ L 30% H₂O₂ per  5 mL DAB solution and mix thoroughly.



23 Incubate sections with DAB-H₂O₂ solution for  00:08:00 .

8m



24 Perform 3× 7 min washes in 1x PBST with gentle agitation.



24.1 Wash in 1x PBST with gentle agitation for  00:07:00 (1/3).


7m



24.2 Wash in 1x PBST with gentle agitation for  00:07:00 (2/3).



7m



24.3 Wash in 1x PBST with gentle agitation for  00:07:00 (3/3).

7m



25 **Mount tissue sections** onto super-frost slides pre-coated with gelatin-chrome alum and allow to dry at  Room temperature  Overnight .




Day 3 (2 days later)

33m


26 Process slide-mounted tissue sections through the following solutions;

27 dH₂O  00:03:00 .


3m

28 50% ethanol  00:03:00 .

3m








29 70% ethanol  00:03:00 .

3m

30 95% ethanol  00:03:00 .

3m



- 31 100% ethanol  00:03:00 . 3m
- 32 100% ethanol  00:03:00 . 3m
- 33 Xylene  00:05:00 . 5m
- 34 Xylene  00:05:00 . 5m
- 35 Xylene  00:05:00 . 5m
- 36 **Coverslip slides with DEPEX mounting media** and allow to dry in the fume hood 
 Overnight before proceeding with microscopy.
- 37 Image sections using bright field microscopy for subsequent xenograft characterization.