



Feb 16, 2024

DAB Staining for GFP on Free-floating Fixed NHP Brain Tissue



Forked from [Standard DAB Staining for Free-floating Fixed NHP Brain Tissue](#)

DOI

dx.doi.org/10.17504/protocols.io.eq2lyjy4mlx9/v1

Aydin Alikaya¹, William Stauffer¹

¹Department of Neurobiology, University of Pittsburgh, Pittsburgh, PA, USA



Andreea Bostan

University of Pittsburgh

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2lyjy4mlx9/v1>

Protocol Citation: Aydin Alikaya, William Stauffer 2024. DAB Staining for GFP on Free-floating Fixed NHP Brain Tissue .
protocols.io <https://dx.doi.org/10.17504/protocols.io.eq2lyjy4mlx9/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: February 07, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 94871

Keywords: ASAPCRN, Immunostaining, DAB, NHP Brain Tissue, fixed nhp brain tissue, nhp brain tissue, green fluorescence protein, brain tissue, fixed brain tissue section, brain tissue section, biotin abc complex, tissue, avidin, gfp


Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-020519

Abstract

This protocol details the procedure for immunohistochemical 3,3'-Diaminobenzidine (DAB) staining of free-floating fixed brain tissue sections using the avidin/biotin ABC complex to detect Green Fluorescence Protein (GFP).

This protocol has been tested with free-floating non-human primate (NHP) brain tissue that has been fixed with 4% paraformaldehyde, cryoprotected with sucrose gradients, and cryo-sectioned  50 μm .

Guidelines

When using 6 well tissue culture plates [Falcon, 353046] to react individual sections, you will need **2+ mL** solutions for *each* well plate.

When using circular staining nets [e.g., Brain Research Laboratories #4115] to react multiple series of sections, you will need **50 mL** solutions for *each* staining net.



Materials

Tissue:

NHP brain tissue sections (50 μm).

Materials/Equipment:

- Tissue culture plates or circular staining nets
- Orbital shaker
- Fume hood
- Nitrile Gloves
- Glass slides (charged or subbed)

Reagents:

- Phosphate-buffered saline (PBS)
- Hydrogen Peroxide: H_2O_2 (30%)
- Distilled water: dH_2O
- Primary Antibody: Thermo Fisher, Molecular Probes Cat# A11122, RRID: AB_221569
- Vectastain ABC-HRP Kit, Peroxidase (Rabbit IgG) (PK-4001, Vector Laboratories)
- Peroxidase (HRP) with Nickel (3,3'-diaminobenzidine) (SK-4100) (Vector Laboratories)

Troubleshooting

Safety warnings

- ! Use appropriate care when using hydrogen peroxide (reactive, can cause skin/eye damage) and DAB (suspected carcinogen). Collect DAB solution for chemical waste disposal.



Part I (Day 1)

3h

1 Bring tissue to Room temperature in Phosphate Buffered Saline (**PBS**, pH 7.2-7.4) on an orbital shaker for 30 minutes. 00:30:00

30m

2 Prepare **Peroxide Solution (0.3 % H₂O₂)** in **dH₂O**.

5m

For 10 mL **0.3% H₂O₂** use:

- 100 µL **30% H₂O₂**
- 9900 µL **dH₂O**

3 Prepare **Normal Goat Serum Blocking Solution** in **PBS**.

5m

To 10 mL **PBS** add:

- 150 µL **Normal Goat Serum** (= 3 drops of serum from Vectastain ABC-HRP Kit, Peroxidase Rabbit IgG PK-4001)

4 Prepare **Primary Antibody Solution (rabbit anti-GFP)** at **1:10000** dilution in **PBS**:

5m

- 1 µL **rabbit anti-GFP** (Thermo Fisher, Molecular Probes Cat# A11122, RRID:AB_221569)
- 9999 µL **PBS**.

5 **Rinse** in PBS on a shaker at Room temperature : **3 × 5 minutes**. 00:05:00

15m

6 Quench endogenous peroxide in **Peroxide Solution (0.3 H₂O₂)** on a shaker at

30m

Room temperature : **30 minutes**. 00:30:00

7 **Rinse** in PBS on a shaker at Room temperature : **3 × 5 minutes**. 00:05:00

15m

8 Incubate in **Normal Goat Serum Blocking Solution** on a shaker at RT: **1 hour**.

1h

01:00:00

DO NOT RINSE after blocking serum.




- 9 Incubate in **Primary Antibody Solution** on a shaker at 4 °C Overnight . 12h
- Part II (Day 2)** 4h
- 10 Bring tissue (in the **Primary Antibody Solution**) to Room temperature on a shaker (30 minutes). 00:30:00 30m
- 11 Prepare **ABC Solution** in PBS (at least 30 minutes before use). 00:30:00 5m
To 10 mL **PBS** add:
▪ **2 drops A** from Vectastain ABC-HRP Kit, Peroxidase Rabbit IgG PK-4001.
▪ **2 drops B** from Vectastain ABC-HRP Kit, Peroxidase Rabbit IgG PK-4001.
- 12 Prepare **Secondary Antibody Solution (1:200)** in PBS. 5m
To 10 mL **PBS** add:
▪ 150 µL (= 3 drops from Vectastain ABC-HRP Kit, Peroxidase Rabbit IgG PK-4001) of **Normal Goat Serum**.
▪ 50 µL (= 1 drop from Vectastain ABC-HRP Kit, Peroxidase Rabbit IgG PK-4001) of **biotinylated goat anti-rabbit IgG secondary antibody**.
- 13 **Rinse** in **PBS** on a shaker at Room temperature : **3 × 5 minutes**. 00:05:00 15m
- 14 Incubate in **Secondary Antibody Solution** on a shaker at Room temperature : **30 minutes**. 00:30:00 30m
- 15 **Rinse** in **PBS** on a shaker at Room temperature : **3 × 5 minutes**. 00:05:00 15m
- 16 Incubate in **ABC Solution** on a shaker at Room temperature : **60 minutes**. 01:00:00 . 1h
- 17 **Rinse** in **PBS** on a shaker at Room temperature : **3 × 5 minutes**. 00:05:00 15m



18 Prepare **DAB Peroxide Substrate Solution** in **dH₂O**.

5m

To use the Vector Labs DAB Peroxidase Substrate Kit (SK-4100):

In  5 mL **dH₂O**:

- **2 drops Reagent 1**
- **4 drops Reagent 2**
- **2 drops Reagent 3**
- [optional] 2 drops of Reagent 4 (Nickel) if a black reaction product is desired



Note: Mix well before use. Use immediately.

19 Incubate in **DAB Peroxide Substrate Solution** on a shaker at  Room temperature :

6m

 00:03:00 -  00:06:00 .

Note: Watch the tissue closely to avoid high background staining.

20 **Rinse** in **PBS** on a shaker at  Room temperature : **3 × 5 minutes.**  00:05:00

15m

21 Mount tissue on charged glass slides in **1:8 PBS** in **dH₂O** and let air dry.

22 Rinse slides with **dH₂O** and let air dry (preferably in a hood).

23 Coverslip clean and dry slides with Cytoseal 60 (Thermo Fisher #830-16).

Protocol references

https://vectorlabs.com/productattachments/protocol/VL_SK-4100_UserGuide_LBL02267.pdf

<https://vectorlabs.com/products/vectastain-abc-hrp-kit-rabbit-igg>