



Aug 27, 2024

DAB Staining for Tyrosine Hydroxylase (TH) on Free-floating Fixed NHP Brain Tissue

 Forked from [Standard DAB Staining for Free-floating Fixed NHP Brain Tissue](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l2q1dql1y/v1>

Andreea Bostan¹

¹Department of Neurobiology, University of Pittsburgh



Andreea Bostan

University of Pittsburgh

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l2q1dql1y/v1>

Protocol Citation: Andreea Bostan 2024. DAB Staining for Tyrosine Hydroxylase (TH) on Free-floating Fixed NHP Brain Tissue . protocols.io <https://dx.doi.org/10.17504/protocols.io.4r3l2q1dql1y/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: August 27, 2024

Last Modified: August 27, 2024

Protocol Integer ID: 106519

Keywords: ASAPCRN, Immunostaining, DAB, NHP Brain Tissue, TH, fixed nhp brain tissue, nhp brain tissue, tyrosine hydroxylase, tyrosine hydroxylase, brain tissue, fixed brain tissue section, brain tissue section, biotin abc complex, tissue section

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-020519

Abstract

This protocol details the procedure for immunohistochemical 3,3'-Diaminobenzidine (DAB) staining for Tyrosine Hydroxylase (TH) of free-floating fixed brain tissue sections using the avidin/biotin ABC complex and the Rabbit Polyclonal Anti-TH NB300-109 (Novus).

This protocol has been tested with free-floating non-human primate (NHP) brain tissue that has been fixed (10% formalin or 4% paraformaldehyde), cryoprotected (sucrose or glycerol gradients), and cryo-sectioned

→ 50 μ m .

Guidelines

This protocol has been optimized for using 6 well tissue culture plates [Falcon, 353046] to react individual sections. You will need **2+ mL** solutions for **each** well plate.

It is recommended to use new culture plates for antibody incubation and discard plates used for DAB.

Plates used for other steps can be washed and re-used.



Materials

Tissue:

Brain tissue sections (20 - 50 μm).

Materials/Equipment:

- Tissue culture plates or circular staining nets
- Orbital shaker
- Fume hood
- Nitrile Gloves
- Glass slides (charged or subbed)

Reagents:

- Phosphate-buffered saline (PBS)
- Phosphate-buffered saline (PBS) with 0.2% Triton-X
- Akoya Blocking Reagent Powder (SKU FP1012) to make TNB Blocking Buffer/
- Hydrogen Peroxide: H_2O_2 (3% or 30%)
- Distilled water: dH_2O
- Primary Antibody: Novus Biologicals Tyrosine Hydroxylase (TH) Polyclonal Rabbit.
- Secondary Antibody (to match the host of the primary antibody): biotinylated Goat-Anti Rabbit.
- Normal Serum Blocking Solution (e.g., Normal Goat Serum, S-1000-20)
- Vectastain Elite ABC Peroxidase Kit (Standard) (PK-6100) (Vector Laboratories)
- Vectastain ABC-HRP Kit, Peroxidase (Rabbit IgG) (PK-4001, Vector Laboratories)
- DAB Substrate Kit

Examples:

Peroxidase (HRP) with Nickel (3,3'-diaminobenzidine) (SK-4100) (Vector Laboratories)

ImmPACT DAB (SK-4105)

Safety warnings

- ⚠ Use appropriate care when using hydrogen peroxide (reactive, can cause skin/eye damage) and DAB (suspected carcinogen). Collect DAB solution for chemical waste disposal.



Part I (Day 1)

3h

- 1 Bring tissue to Room temperature in **phosphate-buffered saline (PBS)** on an orbital shaker for **30 minutes**.
 00:30:00 .
- 2 Prepare **Peroxide Solution (3 % H₂O₂)** in **dH₂O**.
E.g., for 10 mL 3% H₂O₂ use:
 - 1000 µL 30% H₂O₂
 - 9000 µL dH₂O
- 3 Prepare **Blocking Serum Solution: Normal Goat Serum (NGS) in TNB Blocking Buffer**.
E.g., in 10 mL buffer (TNB Blocking Buffer) add:
 - 150 µL NGS (or 3 drops of normal serum if using an ABC kit, e.g. Vectastain ABC-HRP Kit Rabbit IgG PK-4001)
- 4 Prepare **Primary Antibody Solution: anti-TH NB300-109 at 1:1000 TNB Blocking Buffer**.
E.g., for 10 mL **Primary Antibody Solution** use:
 - 9990 µL TNB Blocking Buffer
 - 10 µL anti-TH NB300-109
- 5 **Rinse in PBS with 0.2% Triton X (PBS-Tx)** on a shaker at Room temperature :
3 × 3 minutes. 00:09:00
- 6 Quench endogenous peroxide in **Peroxide Solution (3 % H₂O₂)** on a shaker at
 Room temperature :
2 × 10 - 15 minutes. 00:20:00 - 01:00:00
- 7 **Rinse in PBS with 0.2% Triton X (PBS-Tx)** on a shaker at Room temperature :
3 × 3 minutes. 00:09:00
- 8 Incubate in **Blocking Serum Solution** on a shaker at RT: **1 hour**.
 01:00:00 .

DO NOT RINSE after blocking serum.



- 9 Incubate in **Primary Antibody Solution** on a shaker at 4 °C : **overnight x 3 (60 - 72 hours)**.

72:00:00 .

3d

Part II (Day 2)

4h

- 10 Bring tissue (in the **Primary Antibody Solution**) to Room temperature on a shaker (**45 minutes**). 00:30:00 - 01:00:00

30m

- 11 Prepare **ABC Solution** in **PBS with 0.2% Triton X (PBS-Tx)** (at least 30 minutes before use). 00:30:00 .

5m

- 12 Prepare **Secondary Antibody Solution (1:200)** in **TNB Blocking Buffer**.

5m

In 10 mL TNB Blocking Buffer add:

- 150 µL Normal Goat Serum (NGS) (= 3 drops of normal serum if using an ABC kit, e.g. Vectastain ABC-HRP Kit Rabbit IgG PK-4001).
- 50 µL Secondary biotinylated goat-anti rabbit (= 1 drop of normal goat serum if using an ABC kit, e.g. Vectastain ABC-HRP Kit Rabbit IgG PK-4001).

- 13 **Rinse** in **PBS with 0.2% Triton X (PBS-Tx)** on a shaker at Room temperature : **3 × 3 minutes**. 00:09:00 .

9m

- 14 Incubate in **Secondary Antibody Solution** on a shaker at Room temperature : **30 minutes**. 00:30:00 .

30m

- 15 **Rinse** in **PBS with 0.2% Triton X (PBS-Tx)** on a shaker at Room temperature : **3 × 3 minutes**. 00:09:00 .

9m

- 16 Incubate in **ABC Solution** on a shaker at Room temperature : **60 minutes**. 01:00:00 .

1h



- 17 **Rinse in PBS with 0.2% Triton X (PBS-Tx)** on a shaker at Room temperature : 00:09:00 . 9m
- 18 Prepare **Peroxide Substrate Solution** in dH₂O. 5m
To use the Vector Labs DAB Peroxidase Substrate Kit (SK-4100):
In 5 mL dH₂O:
- 2 drops Reagent 1
 - 4 drops Reagent 2
 - 2 drops Reagent 3
 - [optional] 2 drops of Reagent 4 (Nickel) if a black reaction product is desired
- Note: Mix well before use. Use immediately.**
- 19 Incubate in **Peroxide Substrate Solution** on a shaker at Room temperature : 00:03:00 - 00:06:00 . 6m
- Note:** Watch the tissue closely to avoid high background staining.
- 20 **Rinse** in buffer (e.g. PBS) on a shaker at Room temperature : 00:09:00 . 9m
- 21 Mount tissue on glass slides (subbed or charged) in 1:8 PBS in dH₂O and let air dry.
- 22 Rinse slides with dH₂O and let air dry (preferably in a hood).
- 23 Coverslip clean and dry slides with Cytoseal 60 (Thermo Fisher #830-16).

Protocol references

https://vectorlabs.com/productattachments/protocol/VL_SK-4100_UserGuide_LBL02267.pdf

<https://vectorlabs.com/products/vectastain-elite-abc-hrp-kit-standard>

https://my.akoyabio.com/ccrz_ProductDetails?sku=FP1012&cclcl=en_US