



Nov 17, 2025

DAB staining in PFA-fixed grafted mouse brain slices



Forked from [Immunofluorescence in PFA-fixed grafted mouse brain slices](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vz9585gx1/v1>

Roberto Garcia Swinburn¹, Ernest Arenas¹

¹Karolinska Institute Stockholm

SOX6 mDA differentiation



Roberto Garcia Swinburn

Karolinska Institute Stockholm

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vz9585gx1/v1>

Protocol Citation: Roberto Garcia Swinburn, Ernest Arenas 2025. DAB staining in PFA-fixed grafted mouse brain slices .
protocols.io <https://dx.doi.org/10.17504/protocols.io.rm7vz9585gx1/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 11, 2025

Last Modified: November 17, 2025

Protocol Integer ID: 227002

Keywords: immunohistochemistry, immunofluorescence, animal tissue, dab labeling of protein, grafted mouse brain slice, fixed grafted mouse brain slice, dab labeling, μ m mouse brain slice, conjugated antibody, protein expression, fixed brain slice, protein, dab, brain slice

Abstract

This protocol was used to tag protein expression using HRP-conjugated antibodies on 20 μ m mouse brain slices. DAB labeling of proteins and nuclei was succesful on PFA-fixed brain slices.

Materials

PBS1X

PBTx0.3% [Triton-X 0.3% in PBS]

PBTx0.1% [Triton-X 0.1% in PBS]

Blocking Serum Solution (BSS) [Donkey Serum 5%, BSA 1 mg/ml, PBTx0.1%]

Primary antibodies

Secondary antibodies

ReadyProbes Autofluorescence quenching agent

DAPI

Mounting medium

Protocol materials

 3,3'-Diaminobenzidine tetrahydrochloride **Merck MilliporeSigma (Sigma-Aldrich)** Catalog #D5905



Day one: Primary antibody

19h 35m

- 1 Place up to eight $\pm$ 20-30 μm slices in wells filled with PBS1X.
For free-floating, 12 well plates were used and slow constant shaking to ensure coating without compromising sample structural integrity.
- 2 Wash tissue with PBS1X.
PBS1X can be commercial or self-made, but make sure that it is $\text{pH } 7.4$
- 3 Wash tissue with PBS1X
- 4 Wash with PBTx0.3%
- 5 Wash with PBTx0.1%
- 6 Wash with PBTx0.1%
- 7 Incubate samples in BSS [Goat serum 5%, BSA 1 mg/ml, PBTx0.1%], preferably in plates of smaller wells (e.g 24 well) to optimize volume used.
If 2ndary antibody is hosted in donkey, use donkey serum
- 8 Incubate samples in primary antibody solution (in BSS) at 4°C
We used mouse-anti STEM121 for our protocol.

5m



5m



15m



5m



5m



1h



18h



Day two: Secondary antibody

2h 50m

- 9 Recover slices to 12 well plates. Wash with PBTx0.1%
- 10 Wash with PBTx0.1%

5m



5m



11 Incubate in HRP-conjugated secondary solution (in BSS) Room temperature
We used Goat anti-mouse

2h



12 Wash with PBTx0.1%

5m

13 Wash with PBTx0.1%

5m

DAB incubation (in hood)

18m

14 Prepare DAB solution (fresh, protect against light) as manufacturer's instruction



14.1 Dissolve 1 tablet DAB in 15 ml of Tris-buffered saline, pH 7.6



3,3'-Diaminobenzidine tetrahydrochloride **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D5905**

DAB is carcinogenic, so until solution is discarded, all steps should be performed in a protective hood.

14.2 Add 12 µl of fresh 30% hydrogen peroxide

15 Incubate the samples in DAB solution for 8 minutes

8m

16 Wash with PBS1X

5m

17 Wash with PBS1X

5m

Mounting

18 Prepare mounting media (PBS1X: Glycerol 1:1, 0.2% sodium azide). This media lasts up to a year, although it should be discarded if any growths/deposits are observed.



- 19 Mount stained samples on glass slides using the mounting media.