

Mar 29, 2024

DAB Staining of Fixed Mouse Brain Tissue Sections

DOI

dx.doi.org/10.17504/protocols.io.n92ldm127l5b/v1

Katerina Rademacher¹, Ken Nakamura¹

¹Gladstone Institute of Neurological Disease



Haru Yamamoto

UCSF

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldm127l5b/v1>

Protocol Citation: Katerina Rademacher, Ken Nakamura 2024. DAB Staining of Fixed Mouse Brain Tissue Sections. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n92ldm127l5b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 07, 2024



Last Modified: May 31, 2024

Protocol Integer ID: 93035

Keywords: ASAPCRN, dab staining of fixed mouse brain tissue section, mouse brain tissue section, fixed mouse brain tissue section, dab staining, dab

Funders Acknowledgements:

ASAP

Grant ID: 020529

Abstract

This protocol describes steps for immunohistochemical 1 3,3'-Diaminobenzidine (DAB) staining of free floating fixed mouse brain tissue sections.



Materials

- 24-well netwell insert (Fisher NC9979302)
- dPBS (UCSF Stem Cell Core)
- Tris buffered saline (TBS)
- Methanol
- 30% H₂O₂ (Millipore Sigma H1009)
- Endogenous peroxidase quenching buffer:
 - 10% MeOH, 3% H₂O₂, in 1X TBS
- TBS++++ blocking buffer:
 - 10mL 10X TBS
 - 10mL fetal bovine serum (UCSF Stem Cell Core)
 - 3g BSA (Millipore Sigma A3803)
 - 1g glycine
 - 30mg sodium azide (Millipore Sigma S2002)
 - Bring to 96mL with ddH₂O
 - Store @-20°C in 12mL aliquots
 - Add 500 µL of 10% Triton x-100 (Millipore Sigma 93443) right before use
- Primary antibody (e.g. rabbit anti-TH, Abcam AB152)
- ABC complex - VECTASTAIN[®] Elite ABC-HRP Kit, Peroxidase (Rabbit IgG)
 - Manufacturer: Vector Laboratories
 - Cat#: PK-6101
 - Includes biotinylated goat anti-rabbit IgG secondary antibody (Vector Laboratories BA-1000-1.5)
- 50X 1,3,3'-diaminobenzidine (DAB) dissolved in 0.1M Tris
 - Millipore Sigma D12384
- 0.1 M Tris, pH 8.0
- Ethanol
- Xylene (ThermoFisher Scientific X51)
- Permout mounting medium (ThermoFisher Scientific SP15)
- Superfrost Plus slides (Fisher 12 550 15)

Troubleshooting

Day 1

- 1 Blocking and primary antibody incubation.
 - 1.1 Pull all sections into 1x dPBS in 24-well netwells.
 - 1.2 Wash 3 times with 1x dPBS, 10 mins each.
 - 1.3 Wash 3 times with 1x TBS, 10 mins each.
 - 1.4 Quench endogenous peroxidase activity with quenching buffer, 5 mins at RT. Just covering the sections is sufficient, and drain excess buffer by dabbing on kimwipes before and after.
 - 1.5 Wash 3 times with 1x TBS, 10 mins each.
 - 1.6 Block in TBS++++ for 75 mins at RT on rocker, 2 mL/well.
 - 1.7 Incubate in primary antibody diluted in TBS++++ on rocker at RT O/N, 500 μ L/well in a new 24-well plate without netwell insert.

Day 2

- 2 Secondary antibody incubation and DAB development
 - 2.1 Wash 3 times with 1x TBS, 5 mins each in netwell insert.
 - 2.2 Incubate in secondary antibody at a 1:300 dilution in TBS++++ for 2 hrs, 500 μ L/well in a 24-well plate without netwell insert.
 - 2.3 15 mins prior to end of incubation, prepare ABC complex:
 - Both buffer A and B are 1:300 dilution in TBS++++



- Incubate in RT beaker of water for 30 mins, in drawer

2.4 Wash 3 times with 1x TBS, 5 mins each in netwell insert.

2.5 Incubate 2 hours in ABC complex in TBS++++ at RT on rocker without netwell insert.

2.6 Wash 3 times with 0.1 M Tris buffer, pH 8.0, 5 mins each in netwell insert.

2.7 Prepare 1x DAB solution:

- Thaw aliquots of 50x DAB, 1 tube is around 200 μ L
- Add 50x DAB into 0.1M Tris buffer
- Vortex for 5 sec to mix well, sit in dark for 10 mins
- Right before development, add 30% H₂O₂ (1:10000 dilution), vortex to mix well
- Distribute 2 mL to each well on a 12-well plate

2.8 Prepare an additional 12-well plate with Tris buffer to start washes immediately after development.

2.9 Develop DAB reaction, minimizing time difference between wells. Drain excess buffer by dabbing on kimwipes before and after.

2.10 Wash 2 times with Tris, 5 mins each.

2.11 Wash 2 times with TBS, 5 mins each.

2.12 Wash 2 times with PBS, 10 mins each.

2.13 Sections can be stored in PBS at 4C.

2.14 Mount onto Superfrost Plus slides labeled with pencil and allow to dry completely.

Day 3



- 3 Dehydration and clearing.
- 3.1 Wash 2 × 2min in ddH₂O.
- 3.2 Wash 2 × 2min in 70% EtOH.
- 3.3 Wash 2 × 2min in 95% EtOH.
- 3.4 Wash 2 × 2min in 100% EtOH.
- 3.5 Wash 2 × 2min in Xylene.
- 3.6 Coverslip using Permount mounting medium.
- 3.7 Let dry in fume hood overnight then store in drawer in slide box until imaging.