

May 17, 2019

DABCO Mounting Slide Protocol for *Drosophila melanogaster* embryos

DOI

<https://dx.doi.org/10.17504/protocols.io.2yjgfun>

Ashley Albright¹

¹University of California, Berkeley



Ashley Albright

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.2yjgfun>

Protocol Citation: Ashley Albright 2019. DABCO Mounting Slide Protocol for *Drosophila melanogaster* embryos. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.2yjgfun>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 16, 2019

Last Modified: May 17, 2019

Protocol Integer ID: 23275

Keywords: mounting drosophila melanogaster embryo, drosophila melanogaster embryos protocol, dabco mounting slide protocol, drosophila melanogaster embryo, pbt onto slide, confocal microscopy, slide, conventional fluorescence, fluorescence

Abstract

Protocol for mounting *Drosophila melanogaster* embryos stored in PBS or PBT onto slides for conventional fluorescence and/or confocal microscopy.

Adapted from:

Kosman, D., Mizutani, C.M., Lemons, D., Cox, W.G., McGinnis, E., and Bier, E. (2004)

<https://science.sciencemag.org/content/305/5685/846.long>

Materials

MATERIALS

⊗ 1X PBS (Phosphate-buffered saline)

⊗ Microscope slides

⊗ Glycerol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #G5516**

⊗ 14-Diazabicyclo[2.2.2]octane **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D27802**

⊗ clear nail polish

Safety warnings

! Read SDS for DABCO prior to starting this protocol!



Make Mountant

- 1 In light-shielded  50 mL conical (I cover mine in foil) add  1.25 g of DABCO crystals,  15 mL of 1X PBS, and  35 mL of glycerol.

Safety information

Read SDS for DABCO before starting.

- 2 Mix on rocking platform until the solution is homogeneous. Store at  -20 °C .

Mount Slides

- 3 Take off as much liquid as you can, resuspend embryos in an arbitrary amount of DABCO mountant (depends on how dilute you want the embryos to be, but generally a couple hundred ul). Allow the embryos to settle  01:00:00 to  03:00:00 at  Room temperature or overnight at  4 °C before mounting. Embryos can be stored for months/years in microcentrifuge tubes in light-shielded boxes at  -20 °C .
- 4 Transfer approximately  35 µL of embryos with using pipette tip cut with scissors onto a clean slide between two 22 × 22 mm coverslips approximately 0.5 cm apart in the middle of the slide. Cover with another 22 × 22 mm cover slip. Seal the edges with a quick-drying, transparent nail polish. Slides can be stored for a few weeks at  4 °C in the dark (or a few months at  -20 °C) .




Image source:

https://www.uwo.ca/sci/research/biotron/pdf/microscopy/LSM_SampleGuide.pdf

Note

Keep in mind if you're DAPI staining, that tends to diffuse away. It is better to mount a fresh aliquot of embryos if samples are reanalyzed at a later date.

- 5 Analyze embryos by conventional fluorescence or confocal microscopy.