



Jan 05, 2026

Daphnia magna: DNA extraction for downstream whole genome NGS

DOI

<https://dx.doi.org/10.17504/protocols.io.8epv556j6v1b/v1>



Albano Pinto¹, Joana Luísa Pereira¹

¹CESAM - Centre for Environmental and Marine Studies, Department of Biology, University of Aveiro, Portugal



Joana Luísa Pereira

University of Aveiro

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv556j6v1b/v1>

Protocol Citation: Albano Pinto, Joana Luísa Pereira 2026. *Daphnia magna*: DNA extraction for downstream whole genome NGS. protocols.io <https://dx.doi.org/10.17504/protocols.io.8epv556j6v1b/v1>

Manuscript citation:

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 23, 2025

Last Modified: January 05, 2026

Protocol Integer ID: 235705

Keywords: kit-based DNA extraction, Daphnia, microcrustaceans, high quality DNA, dna from the freshwater microcrustacean daphnia magna, freshwater microcrustacean daphnia magna, dna extraction for downstream whole genome ng, dna extraction, extraction of dna, masterpure complete dna, level dna methylation assessment, downstream whole genome ng, genome sequencing, rna kit, sufficient dna yield, procedures for the extraction, sequencing application, dna, biosearch technology, extraction, rna, sequencing

Funders Acknowledgements:

European Union

Grant ID: 101078991

FCT

Grant ID: 2022.10817.BD

Abstract

This protocol provides optimised procedures for the extraction of DNA from the freshwater microcrustacean *Daphnia magna*, based on the use of a commercial kit (MasterPure Complete DNA and RNA Kit, Epicentre, Biosearch Technologies), added indications for preliminary quality control. By running this protocol, one should obtain sufficient DNA yield towards downstream next-generation whole-genome sequencing applications, e.g. base-level DNA methylation assessment (Methyl-seq); as well, samples should withstand good quality levels for the purpose, namely low-to-minimal contaminants and high molecular weight (low fragmentation).

Materials

Apart from the MasterPure Complete DNA and RNA Kit (Epicentre, Biosearch Technologies), the following materials will be needed to carry out the protocol:

 Isopropanol, 99.5%, molecular biology grade **Thermo Fisher Scientific Catalog #327272500**

 Ethanol 70%

 SYBR Safe DNA Gel Stain **Invitrogen - Thermo Fisher**

 10X TBE Buffer **Thermo Fisher Scientific**

 6X DNA Loading Dye **Thermo Fisher Scientific**

 GeneRuler 1 kb Plus DNA Ladder **Thermo Fisher Catalog #SM1331**

Sterile 1.5-mL Eppendorf tubes

Sterile 15 mL Falcon tubes

Sterile pipette tips

Sterile pestles matching 1.5-mL Eppendorf tubes

Ice

Thermoblock set to 65°C

Refrigerated centrifuge (4°C)

Microcentrifuge for short spins

Vortex

Qubit fluorometer (Qubit 3.0, Invitrogen)

Nanodrop spectrophotometer (NanoDrop 1000 spectrophotometer, Nano-Drop Technologies, Wilmington, DE, USA)

Electrophoresis system

 Daphnids

Protocol materials

 Isopropanol, 99.5%, molecular biology grade **Thermo Fisher Scientific Catalog #327272500**

 Ethanol 70%

 SYBR Safe DNA Gel Stain **Invitrogen - Thermo Fisher**

 6X DNA Loading Dye **Thermo Fisher Scientific**

 GeneRuler 1 kb Plus DNA Ladder **Thermo Fisher Catalog #SM1331**

 10X TBE Buffer **Thermo Fisher Scientific**



Before start

Note the materials, equipment and conditions required for the protocol, especially those items that are required apart from the kit. Ensure also the proper environment to carry out the protocol, including disinfected surfaces, established flows in handling and transfer of samples to prevent cross-contamination, availability of general labware, including disposal vessels, sterile vessels, tube holders as needed, labelling material, gloves, etc.



Tyssue lysis and digestion

5h 41m 20s

- 1 Prepare the lysis solution: dilute 1 μl of Proteinase K (50 $\mu\text{g}/\mu\text{l}$) into 300 μl of Tissue and Cell Lysis Solution (both supplied by the kit) per sample.
- 2 Add 150 μl of Tissue and Cell Lysis solution with proteinase K to each tube with the  Daphnids and grind the daphnids using a pestle.
- 3 Add more 150 μl of Tissue and Cell Lysis solution with proteinase K to each  Daphnids tube and mix by inverting the tubes several times.
- 4 Incubate the  Daphnids tubes at 65°C for 15 min in the thermoblock, and vortex the tubes every 5 minutes to mix.
- 5 Cool the  Daphnids down to 37°C
- 6 Add 1 μl of 5 $\mu\text{g}/\mu\text{l}$ RNase A (supplied by the kit) to each sample, then mix by inverting the tubes several times. Spin-down the samples.
- 7 Incubate the  Daphnids tubes at 50°C for 3 hours in the thermoblock.
- 8 Incubate the  Daphnids tubes on ice for 5 minutes.
- 9 Add 150 μl of MPC protein precipitation solution (supplied by the kit) to each  Daphnids tube. Mix by vortexing for 10 seconds.
- 10 Centrifuge at 4 °C for 10 minutes at 10000xg.
- 11 Discard the pellet (do not disturb the pellet) and transfer the supernatant to new sterile 1.5-mL eppendorf tubes.

5m



20s



30s



20m



20m



2m



3h



5m



30s



15m





12 Add 500 μ l of



Isopropanol, 99.5%, molecular biology grade **Thermo Fisher Scientific Catalog #327272500**

. Mix by gently inverting 40 times.

3m



13 Centrifuge at 4 °C for 10 minutes at 13,000 rpm.

15m



14 Discard the supernatant (isopropanol) **using a micropipette**, carefully to not disturb the pellet. Centrifuge the tubes again for 10 seconds at room temperature and remove the remaining supernatant with a smaller bore pipette tip. Check if the DNA pellet (white) is visible.



15 Add 500 μ l of refrigerated  Ethanol 70% , **gently**, to avoid lifting the pellet. Spin-down the tubes at 4 °C if the pellet is loose. Discard the supernatant (ethanol) **using a micropipette**. Repeat this step 2 times.

5m



16 Dry the pelleted DNA at RT for 30 minutes in the tubes with the lids open.

30m

17 Resuspend the DNA with 35-37 μ l of autoclaved UP water. Leave the DNA at 4 °C for 30 minutes to aid in DNA suspending.

35m



18 Aliquot the samples ($\approx$ 15 μ l) for quality control checks.

5m

19 Store DNA at -80°C or keep it at 4°C if preservation for a short period is intended before further steps.

gDNA quantity and quality check

1h 25m 20s

20 Run DNA quantification using fluorometry (recommended)

20.1 Prepare a 1-20 μ L aliquot of sample (commonly 2 μ L) and add the proprietary kit reagent for DNA quantification (dsDNA BR/HS) with Qubit to complete a final reading volume of 200 μ L. For samples starting from higher biomass dilution might be needed before this step (check the expected DNA concentration vs. reading range).

5m





20.2 Read the solution in the Qubit for DNA quantification. The outcome will be the DNA concentration of the sample in ng/μL. If sample dilution applied, mind it to record the DNA concentration of the sample. For the reference samples as indicated in this protocol ( Daphnids), the expected DNA concentration is generally within 5-10 ng/μL.

20s



21 Run DNA quality screening using spectrophotometry

21.1 Set the Nanodrop spectrophotometer software for nucleic acids assessment

2m

21.2 Pipette a drop (2 μL recommended volume) of the sample into the holder in the reader. No dilution is needed if using similar biomass as indicated for the reference sample in this protocol ( Daphnids)

30s



21.3 Read the sample. The outcome will be:
(i) Total nucleic acids concentration, typically within the range of 150-300 ng/μL if using samples similar to the reference sample in this protocol ( Daphnids). As there is no consistent quantification of gDNA concentration only, this concentration is not particularly informative.
(ii) Absorbance ratio 260/280 (~1.8 should be the target; appreciably lower ratios may indicate the presence of protein, phenol or other contaminants that absorb strongly at or near 280 nm). With this protocol, similar samples to that indicated as a reference ( Daphnids) and using the indicated kit, ratios around 1.9 are commonly reached.
(iii) Absorbance ratio 260/230 (2.0-2.2 should be the target; appreciably lower ratios indicate the presence of contaminants which absorb at 230 nm, such as EDTA, carbohydrates, phenol and some extraction reagents, e.g. TRIzol). With this protocol, similar samples to that indicated as a reference ( Daphnids) and using the indicated kit, ratios around 2.0 can be reached, but more commonly the ratios are within 1.7-1.9; this poses little to no issues in most downstream NGS applications.

30s



21.4 Clean and repeat the read at least once for consistency check. Some variability in readings is relatively common.

2m



22 Run DNA quality screening via electrophoresis

22.1 Prepare a 0.8% agarose gel stained with Invitrogen

 SYBR Safe DNA Gel Stain **Invitrogen - Thermo Fisher** , following the manufacturer

20m





instructions.

- 22.2 Place the gel in an electrophoresis tank filled with an appropriate electrophoresis buffer (e.g.  10X TBE Buffer **Thermo Fisher Scientific** , diluted to 1x) and load the samples with an appropriate loading dye (e.g.,  6X DNA Loading Dye **Thermo Fisher Scientific** , diluted to 1x). In the first slot add the ladder (e.g.,  GeneRuler 1 kb Plus DNA Ladder **Thermo Fisher Catalog #SM1331**).

10m



- 22.3 Run an electrophoresis at 100V for 30 min.

40m

- 22.4 Take your gel to a gel documentation platform and check the bands pattern. Ideally, integer DNA fragments should be on the top of the gel, above the 10,000 bp ladder band.

5m



Acknowledgements

This protocol was optimised and elaborated under the scope of the project EPIBOOST, funded by the European Union (Grant 101078991; doi/org/10.3030/101078991). Views and opinions expressed are however those of the authors only and do not necessarily reflect those of the EU or the European Research Executive Agency; neither the EU or the granting authority can be held for them. Albano Pinto is grateful to FCT for the awarded individual PhD grant with the reference 2022.10817.BD and DOI identifier: <https://doi.org/10.54499/2022.10817.BD>.