

Mar 22, 2023

DAPI-Based Polyphosphate Estimation with Extraction Sufficiency Validation: A Method for Quantifying Polyphosphate from Microalgae Samples

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov1ox87lr2/v1>

Ying-Yu Hu¹, Zoe V. Finkel¹


¹Dalhousie University

Marine Microbial Macroecology Lab
Tech. support email: ruby.hu@dal.ca



Ying-Yu Hu

Dalhousie University



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov1ox87lr2/v1>

Protocol Citation: Ying-Yu Hu, Zoe V. Finkel 2023. DAPI-Based Polyphosphate Estimation with Extraction Sufficiency Validation: A Method for Quantifying Polyphosphate from Microalgae Samples. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.ewov1ox87lr2/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 17, 2023

Last Modified: March 22, 2023

Protocol Integer ID: 77245

Keywords: DAPI, polyphosphate, microtiter plate, microplate, microalgae, fluorescence, fluorometric estimation for polyphosphate, based polyphosphate estimation with extraction sufficiency validation, method for quantifying polyphosphate, quantifying polyphosphate, fluorescence of dapi, based polyphosphate estimation, microalgae samples the utilization, samples in quartz, microalgae sample, polyphosphate, analysis from microalgae, significant background fluorescence, accurate polyp measurement, black microtiter plate, measuring polyp, stained sample, fluorescence, polyp measurement, dapi photobleaching, minimizing dapi photobleaching, significant hurdle in polyp measurement, high background fluorescence, based fluorometric estimation, quartz, dapi, polyp from the sample undetectable, polyp from field sample, insufficient polyp extraction, microalgae, matrix effects in microplate, microplate, extraction, spectrofluorometer, collecting sample

Funders Acknowledgements:

CBIOMES

Grant ID: 549937

SCOPE-Gradients

Grant ID: 723789



Abstract

The utilization of DAPI-based fluorometric estimation for polyphosphate (polyP) analysis from microalgae has become increasingly prevalent in field samples since its publication by Martin P. et al. This technique involves evaluating the fluorescence of DAPI-stained samples in quartz cuvettes using a spectrofluorometer. To reduce the consumption of reagent, time, and labor while minimizing DAPI photobleaching, we have adapted this method to a 96-well black microtiter plate with a black film-covered lid. Additionally, the calculation method has been modified to account for matrix effects in microplates.

Testing the number of treatment rounds necessary to extract all polyP is crucial. However, even when collecting samples from the same field location or cultivation condition, there can be high variability in treatment rounds among replicates, leading to significant background fluorescence and rendering the polyP from the sample undetectable. This challenge is especially prominent when measuring polyP from field samples. Limited sample availability and insufficient polyP extraction, combined with high background fluorescence, make the laborious measurement unpredictable and hinder accurate polyP measurement. This obstacle is a significant hurdle in polyP measurement.

In our assay, we overcome the challenge by validating the sufficiency of extraction for each sample and then measuring the polyP values.

To conduct the assay, roughly 400 uL RNase, 400 uL DNase, and 700 uL proteinase are required for four samples.

Citation

Martin, Patrick & Van Mooy, Benjamin (2026)

. Fluorometric Quantification of Polyphosphate in Environmental Plankton Samples: Extraction Protocols, Matrix Effects, and Nucleic Acid Interference.

Applied and Environmental Microbiology.

<http://doi.org/10.1128/AEM.02592-12>

LINK

Guidelines

1. Extracted polyphosphate must be measured on the same day. Polyphosphate loss has been observed if the extraction is processed days after.
2. The polyphosphate standard aliquot can only be thawed and used once. Do not refrozen and thawed multiple times.



Materials

Chemicals

⊗ Tris Buffer 1M pH 7.0 **Fisher Scientific Catalog #BP1756-500**

⊗ Sodium phosphate glass type 45 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S4379-500MG**

⊗ Proteinase-K **Fisher Scientific Catalog #BP1700-500**

⊗ RNase A: 500 U/mL; RNase T1: 20000 U/mL **Fisher Scientific Catalog #AM2288**

⊗ TURBO DNase 2 U/μL **Fisher Scientific Catalog #AM2239**

⊗ DAPI: 4'6-Diamidino-2-phenylindole dihydrochloride **Fisher Scientific Catalog #D1306**

Sample collection

- 1 Filter microalgae in liquid media onto GFF or PC filters, using gentle vacuum pressure (5 inches Hg).

Equipment

Filter forceps

NAME

blunt end, stainless steel

TYPE

Millipore

BRAND

XX6200006P

SKU

- 2 Rinse sample with filtered artificial seawater (no nutrients)
- 3 Place sample filters in cryogenic vials
- 4 Filter blank media (without cells) through GFF or PC filter as blank.
- 5 Flash freeze filters and stored at  -20 °C
- 6 Freeze dry before measurement.



Equipment

FreeZone® 2.5 L Benchtop Freeze Dryers

NAME

Labconco®

BRAND

700202000

SKU

Preparation of reagents

7 Tris buffer [M] 20 Mass Percent  7.0

Note

Budget:
About 250 mL per four samples

7.1 In a 1 L volumetric flask, top  20 mL [M] 1 Mass Percent  7.0 Tris buffer to 1 L with MilliQ

7.2 Store at  Room temperature

8 PolyP primary standard stock

8.1 Weigh one glass pellet of polyP (45) and write down the weight.



Equipment

Microbalance

NAME

Cubis series

TYPE

Sartorius

BRAND

MSE6.6S-000-DM

SKU

8.2 Transfer the pellet into a 100 mL graduated cylinder.

8.3 Dilute to 100 mL with Tris IMJ 20 Mass Percent pH 7.0

8.4 Aliquot primary stock into 10~50 uL per microtube with Stepper and store at 🌡️ -20 °C

9 PolyP secondary standard stock

If the pellet is far more than 10 mg, dilute primary to secondary to bring down the concentration before preparing working standard

10 Proteinase K IMJ 20 Mass Percent

10.1 Add 🧪 25 mL MilliQ directly into the original package of Proteinase K, vortex to mix

10.2 Aliquot 700 uL to microtubes and keep frozen at 🌡️ -20 °C

11 DAPI primary stock IMJ 14.3 Mass Percent

Add 🧪 2 mL MilliQ directly into the original package and keep frozen at 🌡️ -20 °C

Preliminary extraction efficiency test

12 Prepare boiling bath.

Equipment

VWR® Advanced Hot Plates	NAME
VWR	BRAND
97042-658	SKU

Equipment

Hollow Polypropylene (PP) Ball Bath Covers, 20 mm	NAME
Cole-Parmer	BRAND
UZ-06821-04	SKU

Equipment

Tube rack	NAME
Simport MultiRack™	BRAND
CA48648-606	SKU

13 Prepare  37 °C incubator/shaker.



14 Transfer sample into glass centrifuge tube

Equipment

Disposable Glass Screw-Cap Centrifuge Tubes^{NAME}

10 mL

TYPE

Corning®

BRAND

99502-10

SKU

15 Label centrifuge tube for different samples, place one Pasteur pipet into the tube for transferring extract from the same sample

16 Label 15 mL Falcon tube from 1 to 15 for each one sample.

17 Add  4 mL Tris buffer [M] 20 Mass Percent  7.0 , vortex and then sonicate.

15s

Equipment

Specific Pipette Tips 5mL

NAME

Thermo Scientific™ Finntip™

BRAND

21-377-304

SKU

18 Keep in boiling bath.

5m

**Note**

Make sure the tube rack is in the middle of the boiling bath and covered with PP balls. Tris solution in the tube should be boiling during the 5 minutes' incubation.

19 Sonicate

15s

20 Vortex and then transfer extract to 15 mL Falcon tube, according to the extract number.

Note

Do not push filter to the bottom. Use Pasteur pipet, gently lift the filter upwards, and then transfer as much extract as possible. Gently press the extract out of the filter.

Equipment

Disposable Soda-Lime Glass Pasteur Pipets	NAME
5 3/4"	TYPE
Fisherbrand	BRAND
13-678-6A	SKU

21 Repeat Step 17 to Step 20 until complete 15 times' extraction in total.

22 Centrifuge the extract

5m

 3200 rpm, Room temperature, 00:05:00

23 Use forward pipetting, load black microtitre plate with  200 µL supernatant from the extract (one well for one extract, no need to load replicates).

Tris buffer  20 Mass Percent  7.0 is used as blank.



	1	2	3	4	5	6	7	8	9	10	11	12
A	1-1	1-9	2-1	2-9	3-1	3-9	4-1	4-9				
B	1-2	1-10	2-2	2-10	3-2	3-10	4-2	4-10				
C	1-3	1-11	2-3	2-11	3-3	3-11	4-3	4-11				
D	1-4	1-12	2-4	2-12	3-4	3-12	4-4	4-12				
E	1-5	1-13	2-5	2-13	3-5	3-13	4-5	4-13				
F	1-6	1-14	2-6	2-14	3-6	3-14	4-6	4-14				
G	1-7	1-15	2-7	2-15	3-7	3-15	4-7	4-15				
H	1-8	Tris	2-8		3-8		4-8					

Equipment

96-Well Black Microplates

NAME

Polystyrene

TYPE

Greiner Bio-One

BRAND

655076

SKU

- 24 Prepare DAPI working solution [M] 100 Mass Percent
Dilute  12.6 μL of [M] 14.3 Mass Percent DAPI stock with  1800 μL MilliQ in a foil wrapped microtube and vortex.
- 25 In the dimmed room with only red light bulb on add  24 μL [M] 100 Mass Percent DAPI to each sample in the plate.
- 26 Adhere black film on the top of a microplate lid and cover the plate with this lid.

**Equipment****Black Vinyl Films for Fluorescence and Photoprotection** NAME

VWR

BRAND

89087-692

SKU27 Shake at room temperature for  00:07:00

7m

28 Read fluorescence: excitation at 410 nm and emission at 550 nm

Equipment**Varioskan LUX Multimode Microplate Reader** NAME

Thermo Fisher

BRAND

VL0L00D0

SKU

29 Plot fluorescence intensity versus number of extraction.
The number of extract (N) is the stationary point where the fluorescence of stained extract stops decreasing or the derivative of the fluorescence after that point is close to zero.





If $RFU(15) - RFU(Tris) > 1$, proceed to extract five additional times. And then measure the stained extract following the previous steps.

- 30 Combine Extraction 1 to Extraction N into a falcon tube.

Note

Try to transfer all solution including debris from each tube.
If the total volume is over 50 mL, use a beaker instead.

	Sample code	N	V(Tris) per extract (mL)

- 31

$$C_{extract} = \frac{\sum_{i=1}^N C_i V_i}{V_T} = \frac{\sum_{i=1}^N C_i \times V_{Tris}}{N \times V_{Tris}} = \frac{\sum_{i=1}^N C_i}{N}$$

$$C'_{extract} = \frac{\sum_{i=1}^N C_i (V_i - V)}{V_T - N \times V} = \frac{\sum_{i=1}^N C_i \times (V_{Tris} - V)}{N \times V_{Tris} - N \times V} = \sum_{i=1}^N C_i \times \frac{V_{Tris} - V}{N \times V_{Tris} - N \times V} = \frac{\sum_{i=1}^N C_i}{N}$$

$\Rightarrow C_{extract} = C'_{extract}$ Assuming the volume of extract from each vial is precisely removed for preliminary test.

Enzyme treated extract

- 32 Well mix 1~N extract, transfer 12 mL into 15 mL falcon tube, centrifuge

3200 rpm, Room temperature, 00:05:00

5m

- 33 Transfer 1.8 mL supernatant to a 2 mL tube (Set S).

Note

Sample is triplicated into S1a, S1b and S1c; S2a, S2b, S2c...etc.

- 34 Centrifuge extract "N+1" 3200 rpm, Room temperature, 00:05:00

5m

Note

Blank is duplicated into B1a and B1b; B2a and B2b... etc.

35 Transfer  1.5 mL supernatant into a 2 mL tube (Set B).

36 In Set S, add  18 µL RNase and  18 µL DNase

Note

RNase tends to leave residue in the tip. However one package has only 1 mL RNase, it will be a waste to use reverse pipetting. After dispensing RNase into the vial, use the same tip to draw the solution and gently dispense it back into the solution for about three time, so that there is no residue remaining in the tip. Replace a new tip for the next vial.

Note

Require ~400 uL RNase and ~400 uL DNase.

37 In Set B, add  15 µL RNase and  15 µL DNase

38 Incubate at  37 °C , shake continuously

10m

Equipment

SHAKING INCUBATOR

NAME

71L

TYPE

Corning® LSE™

BRAND

6753

SKU

**Note**

Start the timer when temperature reaches 37 °C

39 Thaw proteinase (~700uL)

40 In Set S, add 36 μ L Proteinase

41 In Set B, add 30 μ L Proteinase

42 Incubate at 37 °C , shake continuously.

30m

Note

Start the timer when temperature reaches 37 °C

Enzyme treated standard amended extract

43 Prepare PolyP working standard [PO₃]⁻ 7.6 Mass Percent

Based on the actual concentration of PolyP (45) primary or secondary standard stock, dilute a certain volume of stock with Tris buffer 20 Mass Percent 7.0

For a final concentration 7.6 Mass Percent

Note

Total volume = 160 X N (ul), where N = sample number

Note

FW(45Na₂O.55P₂O₅)=10600
Mol of PO₃ per mol of PolyP (45) = 110

44 Transfer  840 μL of enzyme treated extract (1~N) into 2 mL tubes (Set A).

Note

Forward pipetting, aspire and dispense for three times to mix.

45 Add  160 μL  7.6 Mass Percent polyP working standard to  840 μL of enzyme treated extract, vortex.

46 Prepare DAPI working solution  100 Mass Percent

Dilute  12.6 μL of  14.3 Mass Percent DAPI stock with  1800 μL MilliQ in a foil wrapped microtube and vortex.

Load microtiter plate

7m

47 Load  200 μL blanks (B: N+1), samples (S: 1~N) and amended samples (A: Amended 1~N) to the microplate. Organize samples as shown in the following scheme:

	1	2	3	4	5	6	7	8	9	10	11	12
A	B1a	B1a	B2a	B2a	B3a	B3a	B4a	B4a	B1a (UN)	B2a (UN)	B3a (UN)	B4a (UN)
B	B1b	B1b	B2b	B2b	B3b	B3b	B4b	B4b	B1b (UN)	B2b (UN)	B3b (UN)	B4b (UN)
C	S1a	S1a	S2a	S2a	S3a	S3a	S4a	S4a	S1a (UN)	S2a (UN)	S3a (UN)	S4a (UN)
D	S1b	S1b	S2b	S2b	S3b	S3b	S4b	S4b	S1b (UN)	S2b (UN)	S3b (UN)	S4b (UN)
E	S1c	S1c	S2c	S2c	S3c	S3c	S4c	S4c	S1c (UN)	S2c (UN)	S3c (UN)	S4c (UN)
F	A1a	A1a	A2a	A2a	A3a	A3a	A4a	A4a	A1a (UN)	A2a (UN)	A3a (UN)	A4a (UN)
G	A1b	A1b	A2b	A2b	A3b	A3b	A4b	A4b	A1b (UN)	A2b (UN)	A3b (UN)	A4b (UN)
H	A1c	A1c	A2c	A2c	A3c	A3c	A4c	A4c	A1c (UN)	A2c (UN)	A3c (UN)	A4c (UN)



Note

Reverse pipetting

- 48 In a dimmed room with only red bulb on, add  24 μL DAPI working solution  100 Mass Percent to each sample in the microplate **except for those labelled with (UN).**

Note

Forward pipetting

- 49 Adhere black film on the top of a microplate lid and cover the plate with this lid.

- 50 Shake at room temperature for  00:07:00

7m

- 51 Shake duration: 1 min
Shaking type: continuous
Shaking speed and force: 600 rpm/High
Fluorescence: excitation at 410 nm and emission at 550 nm
Measurement time: 300 ms
Excitation bandwidth: 5 nm

Calculation

- 52



In the 1~N extraction	In the N+1 extraction	In the amended 1~N extraction
$PolyP_{1-N}(uM) = \frac{\sum_{i=1}^N C_i \times V}{V \times N} = \frac{1}{N} \sum_{i=1}^N C_i$ $Blank_{1-N}(uM) = \frac{\sum_{i=1}^N C_{bi} \times V}{V \times N} = \overline{C_{bi}}$	$Blank_{N+1}(uM) = C_{b(N+1)}$	$PolyP_{amend}(uM) = \frac{840}{1000} \times \frac{1}{N} \sum_{i=1}^N C_i$ $Blank_{amend}(uM) = \frac{840}{1000} \times \overline{C_{bi}}$
Tube set S: 1.8 mL extract 18 uL RNase 18 uL DNase 36 uL Proteinase	Tube set B: 1.5 mL extract 15 uL RNase 15 uL DNase 30 uL Proteinase	Tube set A: 840 uL from Tube set S 160 uL ~ 2 nmol standard
SignalA(DAPI) SignalA(unstained)	SignalB(DAPI) SignalB(unstained)	SignalA(DAPI+std) SignalA(+std, unstained)

$$Signal(extract + enzyme) = [SignalA(DAPI) - SignalB(DAPI)] - [SignalA(unstained) - SignalB(unstained)]$$

$$Signal(extract + enzyme + std) = \left[SignalA(DAPI + std) - \frac{840}{1000} \times SignalB(DAPI) \right] - \left[SignalA(+std, unstained) - \frac{840}{1000} \times SignalB(unstained) \right]$$

$$Signal(std) = Signal(extract + enzyme + std) - \frac{840}{1000} \times Signal(extract + enzyme)$$

53

$$Conc(std)_{uM} = \frac{160}{1000} \times C_{PO3_2nd}$$

$$\frac{Signal(std)}{Conc(std)} = \frac{Signal(extract + enzyme)}{Conc(extract + enzyme)}$$

$$\Rightarrow Conc(extract + enzyme)_{uM} = \frac{Signal(extract + enzyme)}{Signal(std)} \times Conc(std)$$

$$Conc(extract)_{uM} = Conc(extract + enzyme) \times \frac{1800 + 18 + 18 + 36}{1800}$$

$$polyP_{umol/filter} = Conc(extract) \times 0.001 \times V_{Tris/extraction_mL} \times N_{extraction}$$

$$NaPO3_{ug/filter} = polyP_{umol/filter} \times 101.96$$

Citations

Martin, Patrick & Van Mooy, Benjamin. Fluorometric Quantification of Polyphosphate in Environmental Plankton Samples: Extraction Protocols, Matrix Effects, and Nucleic Acid Interference

<http://doi.org/10.1128/AEM.02592-12>