



Mar 11, 2019

Version 3

ddRADSeq in a Field Setting V.3



In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.y2efybe

Mrinalini Watsa^{1,2,3}, Gideon Erkenwick^{1,3}, Aaron Pomerantz^{4,5}, Stefan Prost^{6,7}

¹University of Missouri - Saint Louis; ²Washington University in Saint Louis; ³Field Projects International;

⁴University of California, Berkeley, CA, USA, Department of Integrative Biology;

⁵Marine Biological Laboratory, Woods Hole, MA - USA;

⁶LOEWE-Center for Translational Biodiversity Genomics, Senckenberg Museum, 60325 Frankfurt, Germany;

⁷South African National Biodiversity Institute, National Zoological Garden, Pretoria 0184, South Africa



Mrinalini Watsa

San Diego Zoo Wildlife Alliance

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.y2efybe>

Protocol Citation: Mrinalini Watsa, Gideon Erkenwick, Aaron Pomerantz, Stefan Prost 2019. ddRADSeq in a Field Setting. [protocols.io https://dx.doi.org/10.17504/protocols.io.y2efybe](https://dx.doi.org/10.17504/protocols.io.y2efybe)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 11, 2019

Last Modified: March 11, 2019

Protocol Integer ID: 21286

Keywords: ddrad, ddradseq, minion, nanopore, restriction digest, wildlife, fieldwork, sequencing of genomic dna, sequencing, minion sequencer, illumina sequencer, sequencing run, genomic dna, miseq platform, dna, field laboratory in the southeastern peruvian amazon, single minion, warbler sample, miseq, field laboratory, methods such as blue pippin

Abstract

This protocol begins with blood stored in Longmire's solution and has as its goal, the sequencing of genomic DNA from several individuals onto a single MinION sequencing run. It is based on a protocol outlined by Thrasher et al. (2018) for the same process on a large number of warbler samples, but it differs from this original protocol in a few key ways:

1. Since we ultimately were targeting sequencing by a MinION sequencer, the ultimate processing steps reflect library preparation for runs on this device, instead of Illumina sequencers.
2. For the same reason, we also size selected larger fragments than can be sequenced on a MiSeq platform.
3. This protocol is also meant to only analysed a few individuals (20 or less) since the MinION read coverage is lower than that of a MiSeq.
4. The entire protocol was carried out in a field laboratory in the southeastern Peruvian Amazon, and therefore many of the luxuries a institutional laboratory can afford were not present to us
5. We used gel-based size selection in the absence of methods such as Blue Pippin.
6. We stored blood on FTA cards and in Longmire's solution in our sampling program. We discovered that yields from FTA cards, regardless of including large numbers of hole punches per extraction, simply were too low for the sequencing goals of this project. Therefore, we include the protocol of a WGA amplification step for these low samples. Longmire's solution produced higher yields, but some times, multiple extractions might need to be pooled together with an SPRI bead cleanup to be fully effective.

Materials

STEP MATERIALS

 Buffer AE [Qiagen Catalog #19077](#)

 Buffer AE [Qiagen Catalog #19077](#)



Protocol materials

⊗ Buffer AE Qiagen Catalog #19077

⊗ Buffer AE Qiagen Catalog #19077



⊗ Buffer AE Qiagen Catalog #19077

⊗ Buffer AE Qiagen Catalog #19077

Troubleshooting



P1 Adaptor Ligation

- 1 Resuspend P1 oligos to 100uM in 0.5X AE Buffer (Qiagen), and follow instructions provided by supplier to do so if needed.
- 2 In a PCR strip tube, mix equal volumes of each oligo together to get a final concentration of 50uM of adapter. E.g.  20 μ L of P1F_X and  20 μ L of P1R_X 11.

- 3 In a

Equipment

MiniOne PCR System

NAME

PCR thermocycler

TYPE

MiniOne


BRAND



M4000


SKU

<https://theminione.com/product/minione-pcr-system/>^{LINK}



heat the mixture at  80 °C for  00:01:00



- 4 Remove tube from the thermocycler and place on a rack at  24 °C or room temperature .

- 5 Allow to cool to  24 °C or room temperature for at least  00:30:00 .



- 6 Further dilute the adapter to 5uM with 0.5X  Buffer AE **Qiagen Catalog #19077**



in 1.5ml tube. E.g Add  40 μL of annealed adapter and  360 μL 0.5X AE.

- 7 Store adapters at  4 °C or in freezer ( -20 °C minimum) for long term storage.

P1 1:20 Dilution

- 8 This can be done in a strip of PCR tubes to allow you to use the multichannel pipette to dispense the P1 adapters. (You will need 20 tubes if using all the P1 adapters)
- 9 Label tubes.
- 10 Add  19 μL of nuclease free water to each tube.
- 11 Add  1 μL of P1 adapter to corresponding one tube above.
- 12 Repeat for all P1 adapters being used.
- 13 Seal tubes, vortex briefly to mix, centrifuge.











FTA Card DNA Extractions

- 14 2. Rinse each sample in ddH₂O by vortexing three times for five seconds each
- 15 3. Using a steril tip, transfer each disk to a pre-labeled PCR tube.
- 16 4. Add 45 to 50 μL of ddH₂O (should cover the disk)
- 17 5. Spin down tubes and make sure disk is submerged




- 18 6. Heat to 95°C for 30 minutes in a PCR machine (
- 19 7. Remove tube and pulse vortex or gently tap 60 times.
- 20 8. Centrifuge for 30 seconds
- 21 9. Transfer all liquid into a pre-labelled tube – this will be Elution A
- 22 10. Add 45 to 50 uL (Confirm that you did use 30 ul for each elution) to the disk that remains in the PCR tube and repeat steps 4 – 9 to get Elution B
- 23 11. Store at 4°C if using it soon, or store at -20°C otherwise

P2 Adaptor Ligation

- 24 5In a PCR tube mix equal volumes of each oligo together to get a final concentration of 50uM of adaptor. E.g.  50 µL of P2F_Mspl and  50 µL of P2R_Mspl.
- 25 In a thermocycler heat the mixture at  80 °C for  00:01:00
- 26 Remove tube from the thermocycler and place on a rack at  24 °C or room temperature .
- 27 Allow to cool to  24 °C or room temperature for at least  00:30:00 .
- 28 Further dilute the adaptor to 25uM with 0.5 X  Buffer AE **Qiagen Catalog #19077** in 1.5ml tube. E.g Add  100 µL of annealed adaptor and  100 µL 0.5X AE.



29 Store adapter at  4 °C or in freezer for long term storage.