

Sep 10, 2025

DEFND cryopreserved colon tissue

DOI

<https://dx.doi.org/10.17504/protocols.io.kxygx47wdl8j/v1>

Kathy Plaisance¹, Antonina Mikorska¹, Thierry Voet¹

¹KULeuven Lab of Reproductive Genomics



Kathy Plaisance

KULeuven Lab of Reproductive Genomics

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygx47wdl8j/v1>

Protocol Citation: Kathy Plaisance, Antonina Mikorska, Thierry Voet 2025. DEFND cryopreserved colon tissue. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kxygx47wdl8j/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We tested this protocol, data analysis is still ongoing

Created: August 19, 2025



Last Modified: September 10, 2025

Protocol Integer ID: 226580

Keywords: Cryostat, Homogenization, Permeabilization, colon tissue this protocol, colon tissue, description of defnd, defnd, tissue

Funders Acknowledgements:

Aligning Science Across Parkinson's: ASAP

Grant ID: ASAP-025179

Abstract

This protocol details the description of DEFND cryopreserved colon tissue.

Materials

	A	B	C	D
	Colon lysis buffer	Stock concentration	Final concentration	V (μL)
	TrisHCl pH 7,4	1000 mM	10 %	10.1
	NaCl	5000 mM	10	2.0
	MgCl ₂	1000 mM	3	3.0
	DTT	100 mM	1	10.1
	RNase inhibitor RNasin (Promega)	40 U/μL	1	25.3
	Nuclease-free water			959.5
	Total volume			1010.0

	A	B	C	D
	Wash buffer Promega	Stock concentration	Final concentration	V (μL)
	TrisHCl pH 7,4	1000 mM	10 mM	22.2
	NaCl	5000 mM	10 mM	4.4
	MgCl ₂	1000 mM	3 mM	6.7
	BSA	10 %	1 %	222.0
	Tween-20	10 %	0.1 %	22.2
	DTT	100 mM	1 mM	22.2
	RNase inhibitor RNasin (Promega)	40 U/μL	1 U/μL	55.5
	Nuclease-free water			1864.8
	Total volume			2220.0

	A	B	C	D
	DEFND buffer	Stock concentration	Final concentration	V (μL)



	A	B	C	D
	TrisHCl pH 7,4	1000 mM	10 mM	2.2
	NaCl	5000 mM	10 mM	0.4
	MgCl ₂	1000 mM	3 mM	0.7
	Igepal (NP40 substitute)	10 %	0.1 %	2.2
	cOmplete protease inhibitor	25 X	1 X	8.8
	Pefabloc	20 mg/mL	1 mg/mL	11.0
	RNase inhibitor Protector (Sigma)	40 U/μl	0.5 U/μl	2.8
	Lithium diiodosalicylate (LIS)	100 mM	12.5 mM	27.5
	Nuclease-free water			164.5
	Total volume			220.0

	A	B	C	D
	Nuclei isolation buffer (NIB)	Stock concentration	Final concentration	V (μL)
	TrisHCl pH 7,4	1000 mM	10 mM	10.1
	NaCl	5000 mM	10 mM	2.0
	MgCl ₂	1000 mM	3 mM	3.0
	Igepal (NP40 substitute)	10 %	0.1 %	10.1
	cOmplete protease inhibitor	25 X	1 X	40.4
	RNasin (Promega)	20 U/μl	0.2 U/μl	10.1
	Nuclease-free water			934.3
	Total volume			1010.0

	A	B	C	D
	1X diluted nuclei buffer	Stock concentration	Final concentration	Volume (μL)
	20X nuclei buffer	20 X	1 X	5.0
	DTT	100 mM	1 mM	1.0



	A	B	C	D
	RNase inhibitor Protector (Sigma)	40 U/ μ L	1 U/ μ L	2.5
	Nuclease-free water			91.5
	Total volume			100.0

	A	B	C	D
	MNase reaction buffer	Stock concentration	Final concentration	Volume (μL)
	20X nuclei buffer	20 X	1 X	5.0
	CaCl ₂	1000 mM	5 mM	0.5
	DTT	100 mM	1 mM	1.0
	RNase inhibitor Protector (Sigma)	40 U/ μ L	1 U/ μ L	2.5
	Nuclease-free water			91.0
	Total volume			100.0

	A	B	C	D
	MNase 0,2U/μL	Units in powder	Stock concentration	Volume (μL)
	Lyophilised MNase in vial	50 U	0.2 U/ μ L	
	Nuclease-free water			250.0
	Total volume			250.0

	A	B	C	D
	Stop solution	Stock concentration	Final concentration	V (μL)
	EGTA	500 mM	100 mM	0.6
	Nuclease-free water			2.4
	Total volume			3.0

Remark: CaCl₂ concentration may not be optimal, 5 μ L of 0,1M CaCl₂ added to 500 μ L nuclei and 5 μ L enzyme is recommended. pH may also not be optimal, as it should be between 9 and 10.

Before start

Book necessary instruments and pre-chill centrifuge to 4°C.



Sectioning

- 1 Cool cryostat to -18 °C (blade and chamber).
- 2 Take sample out of -80 °C and place on dry ice.
- 3 Make 6-10 sections of 50 µm and place in DNA LoBind tube.

Homogenization

17m

4

Note

Perform all following steps on wet ice.

1m



Add 500 µL colon lysis buffer and incubate for 00:01:00 .

- 5 Pipet up and down 20 times.



- 6 Add an additional 500 µL colon lysis buffer.



- 7 Incubate On ice for 00:01:00 and pipet up and down 20 times.

1m



- 8 Triturate the sample 20 times using a pellet pestle.

- 9 Spin for 00:05:00 at 500 x g and 4 °C , remove the supernatant and resuspend in 1 mL wash buffer Promega.

5m





10 Transfer the sample with wide-bore pipette onto $\rightarrow$ 40 μm strainer (pluriStrainer mini) placed on Erlenmeyer 2 mL DNA LoBind tube.



11 Wash the Erlenmeyer 1.5 mL tube with additional Erlenmeyer 500 μL wash buffer Promega, transfer onto the filter.



12 Spin the Erlenmeyer 2 mL tube for clock 00:05:00 at rotor 500 x g and temp 4 $^{\circ}\text{C}$.

5m



13 Remove the supernatant after the centrifugation.

14 Resuspend the pellet in Erlenmeyer 700 μL of wash buffer Promega.



15 Filtrate through flowmi and move to a $\rightarrow$ 20 μm filter (pluriStrainer mini) placed on Erlenmeyer 2 mL DNA LoBind tube.

16 Centrifuge clock 00:05:00 at rotor 500 x g and temp 4 $^{\circ}\text{C}$.

5m



17 Remove the supernatant after the centrifugation.

Permeabilization

10m

18 Resuspend in Erlenmeyer 200 μL 12.5 mM DEFND buffer (start incubation time when resuspending).



19 Incubate temp On ice for clock 00:05:00.

5m



20 Immediately following this incubation, add Erlenmeyer 1 mL of NIB and centrifuge the nuclei at temp 4 $^{\circ}\text{C}$ for clock 00:05:00 at rotor 500 x g.

5m



21 Resuspend the pellet in either Erlenmeyer 25 μL 1X diluted nuclei buffer or Erlenmeyer 50 μL MNase reaction buffer (CaCl₂), depending on choice in step 25 (full Multiome experiment or





nucleosome depletion MNase QC).

22 Dilute further if clumping is present and/or the concentration is too high.

23 Count nuclei using Propidium Iodide and Acridine Orange (AO/PI).

24 Nuclei count

Total cell concentration (cells and nuclei):

Live cell concentration (cells):

Dead cell concentration (nuclei):

Viability:

25 Either continue with the full 10X Genomics Chromium Next GEM Single Cell Multiome ATAC + Gene Expression (CG000338 Rev F) protocol without alterations or perform a nucleosome depletion quality control using MNase (Merck), as instructed below.

Note

Adapted from Olsen TR, Talla P, Sagatelian RK, Furnari J, Bruce JN, Canoll P, Zha S, Sims PA. Scalable co-sequencing of RNA and DNA from individual nuclei. Nat Methods. 2025 Mar;22(3):477-487. doi: 10.1038/s41592-024-02579-x. Epub 2025 Feb 12. PMID: 39939719; PMCID: PMC12272635.

Digestion

2m

26 Transfer  50 μL (or other volume) of nuclei resuspend in reaction buffer to different eppendorf tubes, so split conditions from nucleosome depletion into two. 

27 Preheat ThermoShaker to  28 $^{\circ}\text{C}$.

28 Add desired volume of micrococcal nuclease solution (1 unit per  5 μL) to obtain 7 U of enzyme per 1×10^6 cells. 

29 Mix gently and constantly for  00:02:00 in a  28 $^{\circ}\text{C}$ ThermoShaker. 

2m



30 At the end of 2 minutes, add  1 μL (or other volume equal to 1/50 from sample volume) of  0.1 Molarity (M) EGTA to stop. 

31 Negative control reaction buffer and reaction buffer + stop solution without MNase.

DNA extraction

32 Add  100 μL RLT+ buffer to lyse nuclei and vortex. 

33 Extract DNA using DNA Clean & Concentrator® -5 column following manufacturer's protocol.

Protocol references

Olsen TR, Talla P, Sagatelian RK, Furnari J, Bruce JN, Canoll P, Zha S, Sims PA. Scalable co-sequencing of RNA and DNA from individual nuclei. Nat Methods. 2025 Mar;22(3):477-487. doi: 10.1038/s41592-024-02579-x. Epub 2025 Feb 12. PMID: 39939719; PMCID: PMC12272635.