

Jul 10, 2018

# Dengue virus type 2 (DENV-2) capsid-Thai TaqMan assay (no longer in use; see Guidelines)

DOI

[dx.doi.org/10.17504/protocols.io.q4ydyxw](https://dx.doi.org/10.17504/protocols.io.q4ydyxw)

Ian M Mackay<sup>1</sup>, Judy A Northill<sup>1</sup>, Alyssa Pyke<sup>1</sup>


<sup>1</sup>Public Health Virology, Forensic and Scientific Services

Public Health Virology, F...



Ian M Mackay

Public Health Virology, Forensic and Scientific Services



OPEN ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.q4ydyxw](https://dx.doi.org/10.17504/protocols.io.q4ydyxw)

**Protocol Citation:** Ian M Mackay, Judy A Northill, Alyssa Pyke 2018. Dengue virus type 2 (DENV-2) capsid-Thai TaqMan assay (no longer in use; see Guidelines). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.q4ydyxw>

## Manuscript citation:

Development and Evaluation of Serotype- and Group-Specific Fluorogenic Reverse Transcriptase PCR (TaqMan) Assays for Dengue Virus

<http://jcm.asm.org/content/39/11/4119.abstract>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Other

Archived former assay.

**Created:** June 18, 2018

**Last Modified:** July 10, 2018

**Protocol Integer ID:** 13176

## Abstract

This protocol was designed and developed at this laboratory and incorporates a previously published oligoprobe (see below).

The protocol specifically aims to amplify DENV02 viruses and not other dengue viruses. The assay targets the capsid region and is designed as a qualitative test for investigating suspected human cases of DENV-2 infections. This assay has been superseded by the [Dengue virus type 2 \(DENV-2\) MGB TaqMan \(DENV2-2016MGB\) assay](#).

## Guidelines

This was a past assay that we no longer in use.

For the best DENV-2 TaqMan assay, please refer to our recommended protocol:

Dengue virus type 2 (DENV-2) MGB TaqMan (DENV2-2016MGB) assay

<https://www.protocols.io/view/dengue-virus-type-2-denv-2-mgb-taqman-denv2-2016mg-n7kdhw>

## Materials

### MATERIALS

 SuperScript™ III Platinum™ One-Step qRT-PCR Kit **Life Technologies Catalog #11732088**

### STEP MATERIALS

 SuperScript™ III Platinum™ One-Step qRT-PCR Kit **Life Technologies Catalog #11732088**

 SuperScript™ III Platinum™ One-Step qRT-PCR Kit **Life Technologies Catalog #11732088**

## Protocol materials

 SuperScript™ III Platinum™ One-Step qRT-PCR Kit **Life Technologies Catalog #11732088**

 SuperScript™ III Platinum™ One-Step qRT-PCR Kit **Life Technologies Catalog #11732088**

 SuperScript™ III Platinum™ One-Step qRT-PCR Kit **Life Technologies Catalog #11732088**

 SuperScript™ III Platinum™ One-Step qRT-PCR Kit **Life Technologies Catalog #11732088**

## Before start

If using a different brand or model of real-time thermocycler, check the concentration of ROX is adequate. Method assumes the user is familiar with the thermocycler and software used to run the protocol and with PCR in general.

## Oligonucleotide sequences

1	<b>Name</b>	<b>Sequence 5'-3'</b>
	D2TaqC(b)-f <sup>1</sup>	TTCATGGCCCT <b>K</b> GTGGC
	D2Cor05r <sup>1</sup>	CCCCATCT <b>Y</b> TT <b>Y</b> ARTATCCCTG
	D2TaqCor-FAM <sup>2</sup>	<i>FAM</i> - TCCTTCGTTTCCTAACAATCC– <i>TAMRA</i>

1

Designed by Alyssa Pyke, Public Health Virology, 2005;

2

Previously published, Callahan  
*et al*

<http://jcm.asm.org/content/39/11/4119.abstract>

## Reagents

2



SuperScript™ III Platinum™ One-Step qRT-PCR Kit **Life**  
Technologies Catalog #11732088

## Reaction set-up

3 The assay has been used on both a Rotor-Gene 6000 and a Rotor-Gene Q real-time thermocycler

Prepare sufficient mix for the number of reactions.

Include a suitable 'dead volume' as necessary if using a robotic dispenser.

Reagent	Volume (μl) x1	Final reaction concentration
Nuclease-free water	4.44	N/A
D2TaqC(b)-f 150pmol/μl	0.04	300nM
D2Cor05r 150pmol/μl	0.04	300nM
D2TaqCor-FAM 75pmol/μl	0.04	150nM
2X Reaction Mix <sup>1</sup>	10	1X
SuperScript® III/Platinum® Taq Mix <sup>1</sup>	0.4	1X
ROX Reference Dye (25μM)	0.04	0.05μM
Template	5	N/A
<b>TOTAL</b>	20	

1

Superscript

TM

III Platinum

TM

One-step qRT-PCR kit

- Dispense 15µL to each reaction well.
- Add 5µL of template (extracted RNA, controls or NTC [nuclease-free water] ).
- Total reaction volume is 20µL

## Amplification

4


50°C	5min	1X
95°C	2min	1X
95°C	3sec	40X
60°C	30sec <sup>1</sup>	

1

Florescence acquisition step

## Result Analysis

- 5 The definition used for a satisfactory positive result from a real-time fluorogenic PCR should include each of the following:
1. A **sigmoidal curve** – the trace travels horizontally, curves upward, continues in an exponential rise and followed by a curve towards a horizontal plateau phase
  2. A **suitable level of fluorescence** intensity as measured in comparison to a positive control (y-axis)
  3. A **defined threshold (C<sub>T</sub>) value** which the fluorescent curve has clearly exceeded (Fig.1 arrow), which sits early in the log-linear phase and is <40 cycles
  4. A flat or non-sigmoidal curve or a curve that crosses the threshold with a C<sub>T</sub> >40 cycles is considered a negative result.
  5. NTCs should not produce a curve



**Figure 1.** Examples of satisfactory sigmoidal amplification curve shape when considering an assay's fluorescent signal output. The crossing point or threshold cycle ( $C_T$ ) is indicated (yellow arrow); it is the value at which fluorescence levels surpass a predefined (usually set during validation, or arbitrary) threshold level as shown in this normalized linear scale depiction. LP-log-linear phase of signal generated during the exponential part of the PCR amplification; TP-a slowing of the amplification and accompanying fluorescence signal marks the transition phase; PP-the plateau phase is reached when there is little or no increase in fluorescent signal despite continued cycling.