

Aug 04, 2025

DENV NS2B-NS3 fluorescence single point screening assay



Forked from [SARS-CoV-2 nsp3 macrodomain Time-Resolved FRET peptide displacement assay](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvwb6e9vmk/v1>

Haim Barr^{1,2}, Noa Lahav^{1,2}, Jiyun Zhu^{3,2}

¹The Weizmann Institute of Science; ²ASAP Discovery Consortium; ³Stanford University

Haim Barr: General acknowledgement: The Wohl Drug Discovery institute, The Nancy and Stephen Grand Israel National Center for Personalized Medicine.;

ASAP Discovery



Jiyun Zhu

Stanford University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwb6e9vmk/v1>

Protocol Citation: Haim Barr, Noa Lahav, Jiyun Zhu 2025. DENV NS2B-NS3 fluorescence single point screening assay. protocols.io <https://dx.doi.org/10.17504/protocols.io.e6nvwb6e9vmk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 01, 2024

Last Modified: August 04, 2025

Protocol Integer ID: 111386

Keywords: Fluorescence assay, Assay, Inhibitor, Fragment, DENV, Screening, NS2B-NS3 protease, single point, protease inhibitors against denv ns2b, fluorescent assays for zikv ns2b, ns3 fluorescence, ns3 cleavage of substrate bz, ns3 fluorescence single point, screening assay purpose, denv ns2b, ns3 cleavage, screening assay, fluorescent assay, protease inhibitor, peptide substrate, zikv ns2b, assay purpose, enzyme, containing serine, enzyme activity, candidates for medicinal optimization general description, ns3, fluorescence, fluorescence signal, medicinal optimization general description, excitation of amc, enzymatic reaction, lower enzyme activity, result of lower enzyme activity

Funders Acknowledgements:

National Institutes of Health/National Institute Of Allergy and Infectious Diseases (NIH/NIAID)

Grant ID: U19AI171399

Disclaimer

The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.

Abstract

Purpose: performing single point test of selected compounds from screening libraries containing serine-targeting fragments and serine protease inhibitors against DENV NS2B-NS3 to find warheads and candidates for medicinal optimization

General description: This protocol details the fluorescent assays for ZIKV NS2B-NS3 cleavage of substrate Bz-Nle-KRR-AMC peptide. This method measures the fluorescence from the released AMC product as a result of the enzymatic reaction. When hydrolyzed, AMC is liberated from the peptide substrate. Excitation of AMC at 350 nm emits a resonant energy of 450 nm. When the enzyme is inhibited, the fluorescence signal will decrease as a result of lower enzyme activity. The screening method is validated by calculating the Z prime number of each plate, and the percentage of inhibition is calculated and then evaluated for inhibitory efficacy.

Outcome: hits were selected from screening assays with more than 50% inhibition



Materials

Assay Buffer Reagents (Concentration listed are from Stock Solutions)

1. HEPES(Fisher Scientific). Dissolving HEPES powder in MilliQ water and adjust pH to 8.5 for a final concentration of 0.5 M. Filter with 0.2 um filter
2. Sodium chloride (sigma aldrich). Dissolving crystal into MilliQ water for a final concentration of 5 M. Filter with 0.2 um filter
3. glycerol (Sigma Aldrich)
4. Igepal (Sigma Aldrich) dissolving one part of Igepal in 200 part of MilliQ water for a final concentration of 0.5 %
5. TCEP (GoldBio). Dissolving in MilliQ water into final concentration of 1 M. Store in -20 °C

Additional Reagents:

[M] 242140 nanomolar (nM) **DENV NS2B/NS3 Enzyme**

- The Enzyme stock was originally [M] 242140 nanomolar (nM) and was diluted to [M] 10000 nanomolar (nM) before every experiment in **freshly made Assay Buffer**. The final assay concentration is [M] 12.5 nanomolar (nM)

[M] 20000000 nanomolar (nM) **Substrate Bz-Nle-KRR-AMC**

- Substrate stock was dissolved in DMSO to the stock concentration. Before each experiment, the substrate stock was diluted to [M] 10000 nanomolar (nM) in freshly made Assay Buffer. The final assay concentration is [M] 500 nanomolar (nM)

Safety warnings

⚠ Please be sure to wear proper Personal Protective Equipment (PPE) while performing this experiment.

Before start

Thaw TCEP solution on ice to make sure it is fresh



DENV NS2B-NS3 expression and purification

- 1 The protein used in this assay was expressed and purified and received from batch QQ01DVNS2B -c001 -p007 from Centre for Medicines Discovery

Prepare Reagents

- 2 **PREPARE** all of the reagents/buffers required for this experiment.

Assay Buffer

	A	B	C	D	E
	Reagent	Stock	Final	Units	Note
	HEPES pH 8.5	500	10	mM	
	NaCl	5000	50	mM	
	glycerol	100	5	% v/v	
	Igepal	0.5	0.05	%	
	TCEP	1000	0.5	mM	add freshly

Reagents (dilute reagents in assay buffer for required volume)

	A	B	C	D	E
	Reagent	Stock	Prep (2x)	Final in assay plate	Units
	DENV-NS2B/NS3	242140	10000	5000	nM
	Substrate (Bz-Nle-KRR-AMC)	20000000	10000	5000	nM

Prepare 384-well Plate

16m



- 3 **DILUTE** Dilute protein and substrate using the assay buffer
 - Protein dilution: 150 μ L protein stock solution is added into 3480 μ L Assay Buffer
 - Substrate dilution: 1 μ L x 20 mM substrate is added into 2 mL Assay buffer
- 4 **MIX** Add 10 μ L enzyme stock solution into 384-well plates containing inhibitor stocks (column 2-22)
add 20 μ L reaction buffer to A1-H1, I24-P24 (blank control)
add 10 μ L reaction buffer to I1-P1, A24-H24 (substrate control)
MIX 180 μ L diluted enzyme solution with 3.6 μ L x 10 mM DCI (3,4-Dichloroisocoumarin) and aliquot 10 μ L into I2-P2, A23-H23 (positive control inhibitor)
MIX 180 μ L diluted enzyme solution with 3.6 μ L DMSO and aliquot 10 μ L into A2-H2, I23-P23 (no inhibitor control)

Incubate at room temperature for 1 hour
- 5 **REACT** Add 10 μ L substrate solution into enzyme solution in the plate and mix well.

Read Plate Fluorescence


- 6 **READ** and **RECORD** the plate Relative fluorescence units (RFU) using Cytation 3 Multi-Mode Reader (BioTek, Winooski, VT)

Expected result

AMC product will give RFU signal at ex 350/em 450

Experimental Design

- 7 Plate Layout



- Buffer
- Substrate
- Enzyme + Substrate (negative control)
- Enzyme + Substrate + Inhibitor (positive control)
- Enzyme + Substrate + Library compounds

DATA PROCESSING

Calculate the Z prime number of the screening plates to validate the screening:

$$Z' = 1 - 3 * (\text{Std.Dev of positive} + \text{Std.Dev of negative}) / |(\text{Average of positive} - \text{Average of negative})|$$

Calculate the percentage of inhibition:


$$\text{Inhibition \%} = ([\text{RFU no inhibitor}] - [\text{RFU with compounds}]) / (\text{RFU no inhibitor}) * 100\%$$

PLOT DATA

In Prism software, plot percentage of inhibition with corresponding number of compounds

Result

8 Exemplar results are shown



Acknowledgements

The buffer recipe is provided by Noa Lahav from Weizmann Institute of Science