

Apr 29, 2022

DESI imaging mass spectrometry on liver tissue

DOI

dx.doi.org/10.17504/protocols.io.ewov1nze7gr2/v1

Presha Rajbhandari¹, Brent R. Stockwell¹, Brent R Stockwell¹

¹Columbia University, New York, NY

Human BioMolecular Atlas Program (HuBMAP) Method Development Community
Tech. support email: Jeff.spraggins@vanderbilt.edu



Presha Rajbhandari

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.ewov1nze7gr2/v1

Protocol Citation: Presha Rajbhandari, Brent R. Stockwell, Brent R Stockwell 2022. DESI imaging mass spectrometry on liver tissue. [protocols.io https://dx.doi.org/10.17504/protocols.io.ewov1nze7gr2/v1](https://dx.doi.org/10.17504/protocols.io.ewov1nze7gr2/v1)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 28, 2022

Last Modified: April 29, 2022

Protocol Integer ID: 61621

Keywords: DESI, imaging mass spectrometry

Abstract

Acquire DESI imaging mass spectrometry data on liver tissue sections at 40 μ m spatial resolution - to assess distribution of lipids and metabolites within the tissue.

Materials

Major Mix IMS/Tof Calibration Kit (Waters, SKU:186008113)

Leucine Enkephalin (Waters, SKU:186006013)

- 1 Perform mass calibration of the instrument using Sodium formate and ion mobility using CCS Major Mix using electrospray ionization source.

Equipment	
Synapt G2-Si	NAME
Mass Spectrometer	TYPE
Waters	BRAND
176003191	SKU

- 2 Prepare spray solvent: 98:2 methanol:water with 0.1% formic acid and 20pg/μl leucine-enkephalin as lockmass
- 3 Fill the syringe and start the flow at 1.5μl/min
- 4 Set the DESI sprayer angle of 75°, nebulizing gas (N2) pressure of 0.3 PSLM
- 5 Dry the sample slide with tissue section in a vacuum desiccator for 10 minutes
- 6 Set the sprayer distance from the slide and the inlet tube and optimize signal intensity on a spare tissue section
- 7 Mark the corners of the slide with a colored marker pen and scan the slide
- 8 Place the slides on the DESI slide holder
- 9 Open HDI Imaging software (Waters Corp.) and import the slide image



- 10 Select the instrument type, slide holder in use, and mark the slide corners to define image coordinates
- 11 Define the m/z range, polarity and analyzer mode
- 12 Define the region of interest to image, pixel size and scan rate and export the experiment file
- 13 Load the samples on Masslynx software and start acquisition using appropriate data acquisition parameters on MS tune page