

Dec 10, 2024

Designing Overexpression Plasmids for Cellular Studies

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.36wgqd3y3vk5/v1>

Mukesh Kumar¹, Timothy A. Ryan^{1,2}

¹Department of Biochemistry, Weill Cornell Medicine, New York, NY 10065;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, Maryland 20815, USA



Mukesh Kumar

Weill Cornell Medical College

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36wgqd3y3vk5/v1>

Protocol Citation: Mukesh Kumar, Timothy A. Ryan 2024. Designing Overexpression Plasmids for Cellular Studies. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.36wgqd3y3vk5/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: November 26, 2024

Last Modified: December 10, 2024

Protocol Integer ID: 114390

Keywords: ASAPCRN, designing overexpression plasmids for cellular study, designing overexpression plasmid, function in primary cultured neuron, primary cultured neuron, cellular study, protein localization, studies of protein localization, using gibson assembly, neuron, gibson assembly

Funders Acknowledgements:

NIH

Grant ID: NS036942

NIH

Grant ID: NS11739

ASAP

Grant ID: ASAP-024404

Abstract

This protocol provides a detailed stepwise method for generating fluorescently tagged constructs using Gibson Assembly, enabling studies of protein localization and function in primary cultured neurons.



Materials

Plasmids and Genes:

-  CMV mRuby3-Synapsin1a **addgene Catalog #187896**
-  DDHD2 **Genscript Catalog #OHu23306**
-  CPT2 **Genscript Catalog #OHu17939D**

Vectors and Promoters:

-  plenti hSyn PGK1-HALO **addgene Catalog #220910**
- **Restriction Enzymes:** AgeI, Sall, XhoI, PmeI (Acquired from NEB)
-  QIAquick Gel Extraction Kit (250) **Qiagen Catalog #28706**
- **Gibson Assembly:**
 -  NEBuilder HiFi DNA Assembly Master Mix - 50 rxns **New England Biolabs Catalog #E2621L**
- **PCR Reagents:** High-fidelity DNA polymerase

NEB Stable Competent E. coli:

-  NEB Stable Competent E.coli (High Efficiency) - 20×0.05 ml **New England Biolabs Catalog #C3040H**
- PCR machine
- Gel electrophoresis apparatus
- UV gel imager
- Heat block or thermal cycler



PCR Amplification of Inserts

- 1 Design primers for PCR amplification of target genes, ensuring the inclusion of overlapping sequences compatible with Gibson Assembly and restriction sites as needed:
 - 1.1 eGFP-DDHD2: From DDHD2_OHu23306C_pcDNA3.1(+)-N-eGFP, add XhoI and PmeI.
 - 1.2 CPT2: From CPT2_OHu17939D_pcDNA3.1+/C-(k)DYK, add XhoI and PmeI.
- 2 Amplify target genes using a high-fidelity DNA polymerase to minimize errors.
- 3 Verify the PCR product sizes via gel electrophoresis and purify the bands using gel extraction and cleanup kit.

Restriction Digestion of Vectors

- 4 Digest the pLV vector with appropriate restriction enzymes:
 - 4.1 For mRuby-Synapsin: AgeI and SalI.
 - 4.2 For eGFP-DDHD2 and CPT2: XhoI and PmeI.
- 5 Resolve the digested vector on an agarose gel and purify the linearized vector using Qiagen gel extraction kit.

Gibson Assembly Reaction

1h

- 6 Prepare the Gibson Assembly reaction using NEBuilder® HiFi DNA Assembly Master Mix:



6.1 Mix linearized pLV vector and purified PCR-amplified inserts at a 1:3 molar ratio.



6.2 Add  10 μL of NEBuilder® HiFi DNA Assembly Master Mix.



6.3 Adjust the total reaction volume to  20 μL with nuclease-free water.

7 Incubate the reaction at  50 °C for  01:00:00 .

1h



Transformation of Competent E. coli

1h

8 Transform  5 μL -  10 μL of the Gibson Assembly reaction into competent E. coli.

9 Plate the transformed cells on LB agar plates with appropriate antibiotics.

10 Incubate the plates  Overnight at  37 °C .

1h



Verification of Clones

11 Screen colonies by colony PCR or restriction digestion to confirm successful assembly of the constructs.

12 Sequence positive clones to verify the integrity of the constructs.

Plasmid Preparation

1h

13 Grow positive colonies in LB broth with antibiotics  Overnight .

1h

14 Extract plasmid DNA using a plasmid purification kit.



Note

- Primer design is critical for successful Gibson Assembly; ensure adequate overlapping sequences (>15 bases) for seamless assembly.
- Use high-fidelity enzymes for PCR to prevent mutations in amplified inserts.
- Handle restriction enzymes and Gibson Assembly reagents on ice to maintain activity.
- Verify plasmid sequences to ensure correct insert orientation and integrity.