



Nov 24, 2023

# Detection of accessible cholesterol in primary cilia using purified His-ALOD4-mNeon in 3T3 Fibroblasts

 Proceedings of the National Academy of Sciences of the United States of America

DOI

[dx.doi.org/10.17504/protocols.io.rm7vzx2qngx1/v1](https://dx.doi.org/10.17504/protocols.io.rm7vzx2qngx1/v1)

Sreeja V Nair<sup>1</sup>, Ayan Adhikari<sup>1</sup>, Suzanne R Pfeffer<sup>1</sup>

<sup>1</sup>Department of Biochemistry, Stanford University School of Medicine



Suzanne R Pfeffer

Stanford University School of Medicine

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzx2qngx1/v1>

External link: <https://doi.org/10.1073/pnas.2504774122>

**Protocol Citation:** Sreeja V Nair, Ayan Adhikari, Suzanne R Pfeffer 2023. Detection of accessible cholesterol in primary cilia using purified His-ALOD4-mNeon in 3T3 Fibroblasts. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.rm7vzx2qngx1/v1>

**Manuscript citation:**

Nair SV, Jaimon E, Adhikari A, Nikoloff J, Pfeffer SR Lysosomal glucocerebrosidase is needed for ciliary Hedgehog signaling: A convergent pathway contributing to Parkinson's disease. Proceedings of the National Academy of Sciences of the United States of America 122(31). doi: [10.1073/pnas.2504774122](https://doi.org/10.1073/pnas.2504774122)

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** November 23, 2023

**Last Modified:** May 31, 2024

**Protocol Integer ID:** 91374

**Keywords:** ASAPCRN, detection of accessible cholesterol, cilia in rpe cell, accessible pool of cholesterol, cholesterol, accessible cholesterol, different pools of cholesterol, amount of accessible cholesterol, ciliary marker, ciliary hedgehog, primary cilia, fibroblast, cilia, rpe cell, cell, plasma membrane

**Funders Acknowledgements:**

**Aligning Science Across Parkinson's disease**

Grant ID: ASAP-000463

## Abstract

There exist at least three different pools of cholesterol in the plasma membrane: the essential pool, the sphingomyelin-sequestered pool, and the accessible pool (Radhakrishnan et al, 2020). Ciliary Hedgehog signaling is regulated by the accessible pool of cholesterol (Kinnebrew et al, 2019), and His-mNeon-FLAG-ALOD4, a toxin-based probe, can be used to visualize this accessible pool. Here, we present a method for staining and measuring the amount of accessible cholesterol on cilia in 3T3 fibroblasts stably expressing Somatostatin Receptor (SSTR3)-mApple, a ciliary marker. The protocol described here is based upon previously established methods (Kinnebrew et al., 2019; Johnson and Radhakrishnan, 2021) and has also been used successfully to label cilia in RPE cells.

## Guidelines

Purified ALOD4 is best used fresh and will lose activity if frozen; store at 4°C



## Materials

Purified His-tagged mNeon-FLAG-ALOD4 protein (prepared fresh as described here)

BL21 DE3 Rosetta plysS bacterial cells (Sigma # 70956-3)

1mL HiTrap Talon column (Cytiva #28953766)

PD-10 desalting column (Cytiva #17085101)

InstantBlue Coomassie Protein Stain (Abcam #ab119211)

His-tagged mNeon-FLAG-ALOD4 plasmid (Johnson and Radhakrishnan, 2021)

Biorad, 4–20% Mini-PROTEAN<sup>®</sup> TGX<sup>™</sup> Precast Protein Gels, 12-well, 20 µl #4561095

NIH 3T3 cells (ThermoFischer #R76107)

Rat tail collagen (Gibco by Life Technologies # A10483-01)




















Glass coverslips (Fisher Scientific #12-545-81)

Paraformaldehyde 16% Solution (Electron Microscopy Sciences # 15710) diluted fresh in 1X PBS

## Troubleshooting



## His-mNeon ALOD4 purification

- 1 Make a  50 mL starter culture of BL21 DE3 Rosetta plysS cells expressing pRSET+ His-mNeon-FLAG-ALOD4, starting from a fresh colony picked from freshly transformed plates containing carbenicillin (100µg/ml) and chloramphenicol (34µg/ml). Grow  Overnight with shaking in an Erlenmeyer flask at  37 °C .
- 2 Day 2. Using a 6 liter flask, prepare 2 X 2 liters LB by mixing  50 g LB powder with  2 L water. Autoclave for  00:20:00 using a liquid cycle. Allow the media to cool to  Room temperature 20m
- 3 Add carbenicillin(  100 µg /ml )and chloramphenicol (  34 µg /ml ) to each flask and shake at  37 °C until an OD 0.5-0.6 is reached. This usually takes ~2.5 hr.
- 4 Transfer flask to an  18 °C shaker and equilibrate temperature for 10-20 min.
- 5 Add IPTG to achieve a concentration of  0.3 millimolar (mM)
- 6 Induce expression overnight with shaking at  180 rpm at  18 °C
- 7 Day 3. Transfer bacterial culture to 4,  1 L centrifuge bottles.
- 8 Pellet bacteria by spinning at 4000 RPM for 20 min 4°C; discard supernatant.
- 9 Add 1 Roche protease inhibitor tablet to  50 mL of lysis buffer, rotate at  4 °C to dissolve.
- 10 Set up a ~  200 mL beaker with a stir bar at  4 °C in the cold room.




- 11 Resuspend pelleted bacteria On ice or in cold room using ~ 10 mL of lysis buffer for each bottle.
- 12 Vortex the bottles to help resuspend pellet; transfer resuspended pellets to the beaker with a stir bar.
- 13 Take 40 mL of this suspension to the Emulsiflex cell breakage device and have ready, 4 × 65ml polycarbonate ultracentrifuge tubes on ice for the homogenate from the next step. Pass 10 mL lysis buffer through the Emulsiflex once at 0-30K PSI to equilibrate the machine. Pass suspension through the Emulsiflex once at ~50K PSI to lyse cells. Collect directly into polycarbonate tubes On ice .
- 14 Balance tubes and then spin at 40,000 RPM for 00:45:00 in a 45Ti rotor at 4 °C to pellet cell debris. Collect supernatant and filter through a 500 mL , 0.45µm Millipore filter. Collect the filtered supernatant in a 100 mL glass beaker On ice .
- 15 Take a 1 mL HiTrap Talon column (Cytiva #28953766) and wash it with 10 column volumes (CV) of distilled, degassed water. Equilibrate column with 10CV lysis buffer (50mM HEPES pH 8, 500 mM NaCl, 5mM MgCl<sub>2</sub>, 0.5mM TCEP, 10% glycerol, EDTA-free protease tablet)
- 16 Run the lysate through the column. Collect the flow through.
- 17 Wash the column with 30CV wash buffer (50mM HEPES pH 8, 500 mM NaCl, 5mM MgCl<sub>2</sub>, 0.5mM TCEP, 10% glycerol, 20mM imidazole). Collect the wash fractions. Elute with an imidazole gradient of 100mM-500mM. Collect the elution fractions and analyze them by SDS PAGE.
- 18 Pool the fractions having the desired protein and buffer exchange the pool using a PD-10 desalting column (Cytiva #17085101) into the storage buffer (50mM HEPES pH 8, 150mM NaCl, 5mM MgCl<sub>2</sub>, 0.5mM TCEP, 10% glycerol).

45m

## Note

ALOD4 will lose significant activity if frozen. Use fresh or store at 4°C and use within two weeks of purification

19



Purification of His mNeon FLAG ALOD4 using a Hi Trap Talon column (Affinity Chromatography)

Shown here is an example of SDS-PAGE analysis of the Hi Trap Talon purification. Samples were analyzed on 4–20% Precast Protein Gels. The yield is ~7mg from 4L culture. Depending on the purity of the protein after affinity chromatography, we either purify it further by size exclusion chromatography or more commonly, simply carry out buffer exchange using a PD-10 desalting column.

## Cilia Staining

1d 1h 10m

- 20 Seed  $0.3 \times 10^6$  cells onto collagen-coated coverslips. Once attached, starve the cells in media without serum or in low serum for 24:00:00 to induce ciliogenesis.
- 21 Set up an incubation chamber for ALOD4-mNeon staining: place a metal tray (covered with a layer of parafilm) on ice for 20-30 minutes before paraformaldehyde (PFA) fixation.

1d


- 22 Aspirate media and fix cells in 4% PFA in 1X PBS for 00:10:00 at Room temperature .
- 23 Wash the cells three times in 1X PBS; transfer the coverslips to the incubation chamber and incubate them for 01:00:00 in 4 $\mu$ M ALOD4-mNeon diluted in 1% BSA 1X PBS, sterilized by passage through a 0.2  $\mu$ m filter.
- 24 Wash the cells twice in 1X PBS; fix the cells again in 2% PFA diluted in 1X PBS for 5 minutes On ice . Wash the cells again twice in 1X PBS.
- 25 Optional: stain with DAPI (1:1000) for 4 minutes on ice. Wash the cells twice in 1X PBS and mount them onto clean glass slides using 4  $\mu$ l Mowiol. Air-dry the coverslips and proceed to imaging.

10m

1h

## Example of ciliary accessible cholesterol staining

26



Shown here is an example 3T3 cell expressing SSTR3-mApple, starved for 24h to trigger ciliogenesis and stained with purified His-tagged mNeon-FLAG-ALOD4 protein according to this protocol.



## Protocol references

Johnson KA and Radhakrishnan A (2021) The use of anthrolysin O and ostreolysin A to study cholesterol in cell membranes. *Methods in Enzymology* 649:543–566.

Kinnebrew, M., Iverson, E. J., Patel, B. B., Pusapati, G. V., Kong, J. H., Johnson, K. A., ... & Rohatgi, R. (2019). Cholesterol accessibility at the ciliary membrane controls hedgehog signaling. *Elife*, 8, e50051.

<https://doi.org/10.7554/eLife.50051>

Radhakrishnan, A., Rohatgi, R., & Siebold, C. (2020). Cholesterol access in cellular membranes controls Hedgehog signaling. *Nature chemical biology*, 16(12), 1303–1313. <https://doi.org/10.1038/s41589-020-00678-2>