

Jun 23, 2025

Detection of alpha-synuclein aggregates in floating sections using antibody 5G4

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l6yo15vqe/v1>

Bryan Killinger¹, atousa bahrami¹

¹Rush University



Bryan Killinger

Rush University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l6yo15vqe/v1>

Protocol Citation: Bryan Killinger, atousa bahrami 2025. Detection of alpha-synuclein aggregates in floating sections using antibody 5G4. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bp2l6yo15vqe/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 16, 2025

Last Modified: June 23, 2025

Protocol Integer ID: 220724

Keywords: Rehydration, Antigen, Peroxidase quenching , synuclein aggregate, synuclein protein, micron human brain section, disease aggregate, antibody with an epitope, using antibody, disease aggregates from oligomer, including thicker peripheral tissue biopsy, antibody, biomarker, thicker peripheral tissue biopsy, thicker tissue format, use of this antibody, human brain section

Funders Acknowledgements:

Co-Pathologies Drive Neuroinflammation and Progression in PD

Grant ID: ASAP-021030

PROTEOMIC MAPPING OF PATHOLOGICAL INTERACTIONS ACROSS THE SYNUCLEINOPATHY BRAIN

Grant ID: 5R01NS128467

Abstract

5G4 is an antibody with an epitope in the region of aa44-57 of human alpha-synuclein protein, with a high specificity for aggregates. 5G4 reacts with a wide variety of disease aggregates from oligomers to fibrillar type, and shows very little reactivity for physiologic monomeric asyn, making it an excellent generalized "marker" for misfolded pathogenic asyn species. Use of this antibody has been limited to paraffin-embedded sections. Here, we developed a working protocol for fixed floating 40-micron human brain sections. Using the thicker tissue format may help apply 5G4 to a wide variety of specimens, including thicker peripheral tissue biopsies, which are the preferred format for use as a biomarker.

Materials

5G4: Antibody can be purchased from millipore-sigma (Millipore Cat# MABN389, RRID:AB_2716647)

 Anti- α -Synuclein (SNCA) Antibody **Merck MilliporeSigma (Sigma-Aldrich) Catalog #MABN389**

Before start

Mount tissue on slide:

Mount tissue on superfrost plus adhesion slides and dry overnight in oven at 37 (make sure it's very dry).



Day 1: Dehydrate tissues

41m

- 1 Di water for 00:02:00 . 2m
- 2 50% ethanol for 00:03:00 . 3m
- 3 70% ethanol for 00:03:00 . 3m
- 4 95% ethanol for 00:03:00 . 3m
- 5 100% ethanol 2×5 min each in a chemical hood. (100% ethanol is denature) (absolute ethanol is the anhydrous).
- 5.1 100% ethanol 00:05:00 each in a chemical hood. (100% ethanol is denature) (absolute ethanol is the anhydrous. (1/2) 5m
- 5.2 100% ethanol 00:05:00 each in a chemical hood. (100% ethanol is denature) (absolute ethanol is the anhydrous. (2/2) 5m
- 6 Xylenes 2×10 min in the chemical hood.
- 6.1 Xylenes 00:10:00 in the chemical hood. (1/2) 10m
- 6.2 Xylenes 00:10:00 in the chemical hood. (2/2) 10m

Day 1: Rehydration

41m

- 7 Xylenes 2×10 min in the chemical hood.



- 7.1 Xylenes 00:10:00 in the chemical hood. (1/2) 10m
- 7.2 Xylenes 00:10:00 in the chemical hood. (2/2) 10m
- 8 100% ethanol 2×5 min each in a chemical hood. (100% ethanol is denature) (absolute ethanol is the anhydrous).
- 8.1 100% ethanol 00:05:00 each in a chemical hood. (100% ethanol is denature) (absolute ethanol is the anhydrous. (1/2) 5m
- 8.2 100% ethanol 00:05:00 each in a chemical hood. (100% ethanol is denature) (absolute ethanol is the anhydrous. (2/2) 5m
- 9 95% ethanol for 00:03:00 . 3m
- 10 70% ethanol for 00:03:00 . 3m
- 11 50% ethanol for 00:03:00 . 3m
- 12 Di water for 00:02:00 . 2m

Day 1: Antigen Retrieval (~45min)

1h 15m

- 13 Wash sections with sodium citrate buffer for 00:05:00 . 5m
- 14 Microwaving: 00:10:00 in citrate buffer (pH 6). 10m
- 1000 L H₂O
 - 2.94 g sodium citrate (dihydrate)
 - PH to 6.0





-  0.5 mL Tween -20

15 Let citrate acid to cool ~  00:30:00 .

30m

16  00:30:00 in 88% formic acid.

30m

17 Quick wash with DM.



DM preparation:

	A	B
	Tris-HCl	50mM
	NaCl	150mM
	Triton X-100	0.05%

Day 1: Blocking and peroxidase quenching (1h 40 min)

1h 40m

18 Make blocking buffer: DM + 3% BSA + 3% serum+ Triton.



1. For 100 mL:

	A	B
	DM	100 mL
	BSA	3 g
	Serum	3 mL
	Triton	0.5 mL

2. Serum must be from host species of secondary antibody

3. Take  50 mL out for blocking buffer and leave the other 50 for primary antibody

4. Add peroxidase quenching to the blocking buffer solution

- For 50ml:



	A	B
	30% H ₂ O ₂	0.5 mL
	Sodium azide (stock 10%)	0.5 mL

19 Dip slides in blocking buffer for  01:30:00 at  Room temperature .

1h 30m



20 Dip in TBS 2×5 min each.

Note

This removes excess DM which interferes with the hydrophobic barrier.

20.1 Dip in TBS  00:05:00 each. (1/2)

5m

20.2 Dip in TBS  00:05:00 each. (2/2)

5m

Day 1: Primary antibody

21 5G4 should be diluted 1:5000.

22 Cover the tissue section in primary antibody solution.  Overnight at  4 °C .



Day 2: Secondary antibody (1h 50 min)

1h 50m

23 Dip in DM 3×5 min.



23.1 Dip in DM  00:05:00 . (1/3)

5m

23.2 Dip in DM  00:05:00 . (2/3)

5m

23.3 Dip in DM  00:05:00 . (3/3)

5m

24 Dip in TBS 2×5 min.

24.1 Dip in TBS  00:05:00 . (1/2)

5m

24.2 Dip in TBS  00:05:00 . (2/2)

5m

25 Make secondary blocking buffer:



- For 50 mL working solution:

	A	B
	DM	50 mL
	BSA	1 g
	Serum	2 mL

26 Dilute secondary antibody (Horse Anti-Mouse IgG Antibody (H+L), Biotinylated (BA-2000-1.5) from VectorLabs) in secondary blocking buffer (for biotinylated secondary 1:200).

27 Cover sections in the secondary antibody for  01:00:00 at  Room temperature .

1h

28 Dip in DM 3×5 min.



28.1 Dip in DM 00:05:00 . (1/3)

5m

28.2 Dip in DM 00:05:00 . (2/3)

5m

28.3 Dip in DM 00:05:00 . (3/3)

5m

29 Dip in TBS 2×5 min.

29.1 Dip in TBS 00:05:00 . (1/2)

5m

29.2 Dip in TBS 00:05:00 . (2/2)

5m

Day 2: ABC (35 min)

1h 20m

30 Prepare ABC reagent. Must be prepared 00:30:00 prior to use. Stable for 24:00:00 after prep.

1. For 25mL:

- 2.5 mL blocking buffer + 1 drop A + 1 drop B
- After 30 minutes, fill to 25 mL

31 Cover sections with prepared ABC reagent for 01:15:00 at Room temperature .

1h 15m



32 Dip in DM 00:05:00 .

5m

Day 2: DAB (~1 hr)

1h 7m

33 Make imidazole buffer (IB).





- For  1 L : 6.8 sodium acetate + 0.68 imidazole.
- Once IB is homogeneous, separate  100 mL into a separate beaker with stirbar.
- For  900 mL of IB, pH to 7.2-7.4.
- For  100 mL of non-pH'd IB, add  50 mg DAB followed by  2.5 g nickel. Dissolve completely.

34 Dip slides in pH'd IB (i.e. pH 7.2-7.4) 3×10 min each.

34.1 Dip slides in pH'd IB (i.e. pH 7.2-7.4)  00:10:00 each. (1/3)

10m

34.2 Dip slides in pH'd IB (i.e. pH 7.2-7.4)  00:10:00 each. (2/3)

10m

34.3 Dip slides in pH'd IB (i.e. pH 7.2-7.4)  00:10:00 each. (3/3)

10m

35 Immediately before use add hydrogen peroxide to DAB solution.



- For 100mL DAB solution add  16 μ L of 30% H_2O_2 .

36 Dip sections in DAB solution for development.

17m

- Strong signals  00:00:30 -  00:02:00
- Medium signals,  00:02:00 -  00:05:00
- Weaker signals  00:05:00 -  00:10:00 .
- Development up to 10 min is typically free of background. Over 10 min and background may develop.

37 Monitor color development. (Nickel DAB = jet black, DAB = brown).

- To check color development slides can be placed in a coplin jar with pH'd IB and then checked under microscope.
- If color intensity is too weak, place sections back into DAB solution for additional time.
- Can repeat this procedure until desired color intensity is achieved.



38 Once done developing wash sections in PBS 2×10min.

38.1 Once done developing wash sections in PBS  00:10:00 . (1/2)

10m

38.2 Once done developing wash sections in PBS  00:10:00 . (2/2)

10m

39 Deactivate DAB solution with bleach and dispose in DAB waste container.