



Aug 09, 2024

# Determination of free and protein-bound DA and NE, 5HT and Ach and their metabolites and oxidation products by UPLC-MS/MS method

 In 1 collection

DOI

[dx.doi.org/10.17504/protocols.io.8epv5rmj6g1b/v1](https://dx.doi.org/10.17504/protocols.io.8epv5rmj6g1b/v1)

Miquel Vila<sup>1</sup>, Marta Gonzalez-Sepulveda<sup>1</sup>

<sup>1</sup>VHIR-CIBERNED-ASAP

Vilalab Public



Miquel Vila

VHIR-CIBERNED-ASAP

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5rmj6g1b/v1>

**Protocol Citation:** Miquel Vila, Marta Gonzalez-Sepulveda 2024. Determination of free and protein-bound DA and NE, 5HT and Ach and their metabolites and oxidation products by UPLC-MS/MS method. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.8epv5rmj6g1b/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 08, 2024

**Last Modified:** August 09, 2024

**Protocol Integer ID:** 104953

**Keywords:** oxidation products by uplc, oxidation product, metabolite, ms method protocol for the determination

## Abstract

Protocol for the determination of free and protein-bound DA and NE, 5HT and Ach and their metabolites and oxidation products by UPLC-MS/MS method

## Materials

### External standards

**DA, DOPA, NE, DOPAC, DOMA, DOPE, VMA** → up to 3000 nM in 25 mM FA in water

**3MT, AC, VMA, 5-SCDA, 5-SCD** → up to 1000 nM in 25 mM FA in water




### Internal standard (IS)

**DA-4d** → 500 nM in 25 mM FA in water

## Troubleshooting

## Preparation of the aminochrome (AC) external standard

1m

- 1 Mix  500  $\mu$ L of [M] 1 millimolar (mM) dopamine (DA) with  500  $\mu$ L of [M] 2 millimolar (mM) KIO<sub>4</sub> dissolved in [M] 100 micromolar ( $\mu$ M) aqueous ammonium acetate buffer  5.8 at RT with vigorous shaking for 1 min.

### Note

Following oxidation, aminochromes are placed on ice to prevent further decomposition. Significant degradation of all aminochromes occurs at both RT and 4 °C within 24 h and -20°C at 24-48h (Ochs 2005; Lemos-Amado 2001).

## Preparation of calibration curves

1h

- 2 Prepare a stock solution of the IS in [M] 25 millimolar (mM) FA in water and store it at -80 °C.

Prepare fresh solutions of each metabolite in [M] 25 millimolar (mM) FA in water and use them to make four mixtures: MIX1 [dopamine (DA), noradrenaline (NA), 3-methoxytyramine (3MT), 3,4-dihydroxyphenylalanine (L-DOPA) and aminochrome (AC)], MIX2 [3,4-Dihydroxyphenylacetic acid (DOPAC), 3,4-dihydroxymandelic acid (DOMA) and vanillylmandelic acid (VMA)], MIX3 [5SCD and 5SCDA] and MIX4 [serotonin(5-HT), tryptophan (Trp) and 5-hydroxyindole-3-acetic acid (5-HIAA)].

- 3 Serially dilute mixtures with [M] 25 millimolar (mM) FA in water to obtain the concentration series used in calibration curves.

### Note

Typically, final calibration levels cover a range of 1.72–3000 nM for DA, NE, and MIX2 and 0.39–1000 nM for L-DOPA, 3MT, AC and MIX3.



- 3.1 Serially dilute acetylcholine (ACh) standard with [M] 0.1 % volume FA in acetonitrile (ACN) to obtain a calibration curve covering a range of 0.2-1000 nM
- 4 Homogenize control samples (i.e brain, intestines, heart, blood serum, cells...) in the appropriate volume of [M] 250 millimolar (mM) FA
- 5 Distribute the sample into 90  $\mu\text{L}$  aliquots prior to the addition of 30  $\mu\text{L}$  of the appropriate working mixture (MIX1, MIX2, MIX3 or MIX4), 96  $\mu\text{L}$  of [M] 25 millimolar (mM) FA and 24  $\mu\text{L}$  of [M] 8 micromolar ( $\mu\text{M}$ ) IS (DA-4d + 5HT-4d) .
- 5.1 For ACh calibration curve, distribute the sample into 90  $\mu\text{L}$  aliquots prior to the addition of 30  $\mu\text{L}$  of the working mixture and 120  $\mu\text{L}$  of [M] 200 nanomolar (nM) IS (ACh-4d) in 0.1% FA in ACN.
- 6 20000 rcf, 4°C, 00:10:00 10m
- 7 Transfer supernatant to an Ostro protein precipitation and phospholipid removal plate (Waters, USA) to filter it.
- Save the pellet for protein-bound determinations (see below)
- 8 Finally, inject 7  $\mu\text{L}$  into the UPLC-MS/MS system.

## Sample preparation








2h

- 9 Add 300  $\mu\text{L}$  of [M] 250 millimolar (mM) FA to each brain, intestine, heart or cell pellet sample prior homogenization. Dilute blood serum samples 1:10

### Note


Due to the poor stability of aminochrome, usually a maximum of 50 samples can be analyzed at a time



- 10 Take a  20  $\mu\text{L}$  20  $\mu\text{L}$  sample for protein determination (diluted 1/5 in  25 millimolar (mM) FA )
- 11 Take  55  $\mu\text{L}$  of sample for ACh determination: dilute 1:4 with  0.1 % volume FA in ACN containing a final concentration of 100 nM of ACh-4d as IS.
- 12 Take  216  $\mu\text{L}$  for catecholaminergic and serotonergic determination and add  24  $\mu\text{L}$  of  8 micromolar ( $\mu\text{M}$ ) IS (DA-4d + 5HT-4d) .

#### Note


Important!!!: ensure the concentration of IS is exactly the same in both calibration curves and samples

- 13  20000 rcf, 4°C, 00:10:00

10m

#### Note

The supernatant is used to determine free neurotransmitters and metabolites (that is, those present in the deproteinated supernatant) while the pellet is used for protein-bound determinations (that is, those present in the acid-insoluble pellet and released by HCl hydrolysis)

- 14 Transfer supernatant to an Ostro protein precipitation and phospholipid removal plate (Waters, USA) to be filtered.
- 15 Inject  7  $\mu\text{L}$  of filtered supernatant samples into the UPLC-MS/MS system


## Reductive HCl hydrolysis of resulting pellets


18h

- 16




**Safety information**

Work in fume hood during all the procedure



After removal of the supernatant, wash the pellet (from both calibration curves and samples) with  1 mL of chloroform: methanol (1: 1, v/v) by vortex mixing



17  20000 rcf, 4°C, 00:10:00



10m

18 Transfer the resulting pellets to a sealed-capped tube with  6 Molarity (M) HCl containing  5 % volume thioglycolic acid and  1 Mass Percent phenol



**Note**

- **Calibration curves** --> add  280 µL of the mixture and  40 µL of the corresponding calibration curve working mixture


- **Samples** --> add  288 µL of the mixture and  32 µL of IS

19 Purge tubes with a stream of nitrogen, seal them and heat them at  110 °C for  16:00:00

16h

20 Let tubes cool at  4 °C for at least  00:30:00

30m

21  20000 rcf, 4°C, 00:10:00















10m

22 Treat the supernatant with with acid-washed alumina to extract catecholic compounds

**Alumina extraction of catecholic compounds**

1h



- 23 Transfer a  100  $\mu\text{L}$  aliquot of each hydrolysate into a new Eppendorf tube containing  50 mg of acid-washed alumina and  200  $\mu\text{L}$  of [M] 1 Mass Percent  $\text{Na}_2\text{S}_2\text{O}_5$  - [M] 1 Mass Percent  $\text{EDTA}\cdot 2\text{Na}$
- 24 Add  500  $\mu\text{L}$  of [M] 2.7 Molarity (M) Tris. HCl - [M] 2 Mass Percent  $\text{EDTA}\cdot 2\text{Na}$   9 to the mixture
- 25  1100 rpm, 22°C, 00:05:00 on a microtube mixer 5m
- 26  20000 rcf, 00:10:00 10m
- 27 Remove the aqueous layer by aspiration and wash alumina with  1 mL of Milli-Q water
- 28  20000 rcf, 00:10:00 10m
- 29 Remove the aqueous layer by aspiration and wash alumina with  1 mL of Milli-Q water
- 30  20000 rcf, 00:10:00 10m
- 31 Remove the aqueous layer by aspiration and wash alumina with  1 mL of Milli-Q water
- 32  20000 rcf, 00:10:00 10m
- 33 Remove the aqueous layer by aspiration
- 34 Elute catechols from alumina with  100  $\mu\text{L}$  of [M] 0.4 Molarity (M)  $\text{HClO}_4$  by shaking for 2 min 2m



- 35 Collect all liquid into the injection plate without taking alumina

**Note**

Alumina is discarded after extraction

- 36 Finally, inject  7  $\mu\text{L}$  into the UPLC-MS/MS system.

## UPLC-MS/MS analysis for catecholaminergic, serotonergic and cholinergic determination

- 37 A Waters Acquity™ UPLC system is coupled with a Xevo TQ-S triple quadrupole mass spectrometer with electrospray ionization interface (Waters). Instrument control, data acquisition, and analysis is performed using MassLynx V4.1 (Waters).
- 37.1 An Acquity HSS T3 (1.8  $\mu\text{m}$ , 2.1 mm  $\times$  100 mm) column coupled to an Acquity HSS T3 VanGuard (100 Å, 1.8  $\mu\text{m}$ , 2.1 mm  $\times$  5 mm) pre-column is used to detect MIX1-3 analytes. Column temperature is set at 45 °C and samples are maintained at 6 °C in the thermostatic autosampler.
- 37.2 An Acquity UPLC BEH C18 (1.7 $\mu\text{m}$ , 2.1 $\times$ 100mm) column coupled to a Acquity BEH C18 1.7 $\mu\text{m}$  VanGuard pre-column is used to detect MIX4. Column temperature is set at 55 °C and samples are maintained at 6°C in the thermostatic autosampler.
- 37.3 A Cortecs UPLC HILIC (1.6  $\mu\text{m}$ ; 2.1 $\times$ 75 mm) column coupled to a Cortecs UPLC HILIC VanGuard pre-column (Waters) is used for ACh determination. Column temperature is set at 50 °C and samples are maintained at 6 °C in the thermostatic autosampler.
- 38 The mobile phase for MIX1-4 consisted of solvent A (methanol 100%) and solvent B (25 mM FA in MQ water) at a flow of 0.4 mL/min with the following gradient profiles (depending on the MIX):

MIX1 and MIX2:

0.5% B maintained for 0.5 min, 5% B at 0.9 min and maintained for 2.1min, 50% B at 2.8 min and maintained for 1.2 min, 0.5% B at 4.1 min followed by 0.2 min of equilibration. Total run time 4.3 min.

0.5



### MIX3:

0.5% B maintained for 0.5 min, 8% B at 2.6 min, 50% B at 2.9 min and maintained for 0.6 min, 0.5% B at 3.7 min maintained 0.2 min for equilibration. Total run time 3.7 min

### MIX4:

1%B maintained for 0.5 min, 25 % B at 3 min, 50 % B at 3.1 min and maintained for 0.5 min, 1 % B at 3.6 min maintained 0.4 min for equilibration.

- 38.1 The mobile phase for ACh determination consisted of solvent A (0.1% FA in ACN) and solvent B (10 mM ammonium acetate in MilliQ water) at a flow of 0.5 mL/min with isocratic 70% A- 30% B conditions during 2.2 min.
- 39 The mass spectrometer detector operates under the following parameters: source temperature 150 °C, desolvation temperature 450 °C, cone gas flow 50 L/hr, desolvation gas flow 1100 L/hr and collision gas flow 0.15 mL/min. Argon is used as the collision gas. The capillary voltage was set at: 0.5 kV (MIX1, MIX3 and MIX3-PB), 2 kV (MIX2) or 3kV (MIX4, ACh). The electrospray ionization (ESI) source was operated in both positive and negative modes, depending on the analyte.
- 40 Multiple Reaction Monitoring (MRM) acquisition settings for the targeted metabolites are summarised in the following Table

A	B	C	D	E	F	G
<b>Table 1. MRM acquisition settings a</b>						
Analyte	MRM transition (m/z)	MIX	RT (min)	CV (V)	CE (eV)	CpV (kV)
<b>ACh</b>	145,98 > 86,80	ACh	1,5	10	15	3
<b>ACh-d4</b>	150 > 91	ACh	1,5	28	12	3
<b>NE b</b>	151,75 > 106,94	1	0,69	15	20	0,5
<b>DA-d4 (IS)</b>	157,83 > 94,8	1,2,3	1,44	10	20	0,5
<b>DA</b>	153,93 > 90,57	1	1,46	10	20	0,5

A	B	C	D	E	F	G
<b>L-DOPA</b>	198,1 > 152,1	1	1,48	15	15	0,5
<b>3MT b</b>	150,7 > 90,96	1	3,09	35	20	0,5
<b>AC</b>	149,61 > 121,91	1	3,36	25	25	0,5
<b>DOMA c</b>	182,86 > 136,85	2	1,62	20	14	2
<b>VMA c</b>	197 > 136,9	2	3,61	20	20	2
<b>DOPAC c</b>	166,99 > 122,82	2	3,72	18	22	2
<b>5SCDA</b>	273,1 > 166,9	3	1,73	20	20	0,5
<b>5SCD</b>	317 > 154,86	3	2,01	24	30	0,5
<b>5HT</b>	177 > 160	4	0,97	10	5	3
<b>5HT-d4 (IS)</b>	181 > 164	4	0,97	10	5	3
<b>5HIAA</b>	192 > 146	4	1,5	25	20	3
<b>Trp</b>	205 > 188	4	2,1	15	10	3
a RT, retention time; D, dwell time; CV, cone voltage; CE, collision energy; CpV, capillary voltage.						
b Parent mass after loss of water.						

	A	B	C	D	E	F	G
	c Detected in negative mode.						

## Data analysis and representation

- 41 Samples with a concentration between the limit of detection (LOD) and limit of quantification (LOQ) or bigger than LOQ are considered acceptable; samples with a concentration lower than LOD are considered as the LOD value.
- 42 Catechol oxidation is measured using the formula  $AC + 5SCDA + 5SCDA - PB / DA + 5SCD + 5SCD - PB / L - DOPA$
- 43 DA synthesis is measured using the formula  $DA + NE + DOMA + VMA + 3MT + DOPAC / L - DOPA$
- 44 DA degradation is measured using the formula  $DOPAC + 3MT / DA$
- 45 Data is normalised by the protein concentration (determined by BCA) and presented as the percentage of the wt concentration or ratio.