

Oct 03, 2019

# Determining Chlorophyll Concentration using CuSO<sub>4</sub> Magnesium-Copper Exchange Titration

DOI

[dx.doi.org/10.17504/protocols.io.7xrhpm6](https://dx.doi.org/10.17504/protocols.io.7xrhpm6)



Victor vmr Rodriguez<sup>1</sup>

<sup>1</sup>Independent Researcher



Victor vmr Rodriguez

Independent Researcher

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.7xrhpm6](https://dx.doi.org/10.17504/protocols.io.7xrhpm6)

**Protocol Citation:** Victor vmr Rodriguez 2019. Determining Chlorophyll Concentration using CuSO<sub>4</sub> Magnesium-Copper Exchange Titration. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.7xrhpm6>

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development

We are still developing and optimizing this protocol

**Created:** October 03, 2019

**Last Modified:** October 03, 2019

**Protocol Integer ID:** 28369

**Keywords:** chlorophyll, botany, plant biology, copper II sulfate

## Abstract

This protocol is designed to be able to extract and analyze the concentration of chlorophyll within a sample on a molecular level. The procedures of this protocol require using copper II sulfate and hydrochloric acid as a means of titration to determine the approximate level a chlorophyll within a given sample.

## Guidelines

For proper extraction in titration of Chlorophyll concentration you must be able to measure and transfer liquids within a hundred micro-liters ensure that samples are separated and free of contaminants.

## Materials

### MATERIALS

☒ Magnesium sulfate heptahydrate

☒ Hydrochloric Acid

☒ Copper (II) sulfate pentahydrate **Bio Basic Inc. Catalog #CDB0063.SIZE.500g**

☒ Acetone **Nacalai Tesque Inc. Catalog #00310-95**

☒ Distilled Water **Thermo Fisher Catalog #15230196**

### STEP MATERIALS

☒ Acetone **Nacalai Tesque Inc. Catalog #00310-95**

☒ Magnesium sulfate heptahydrate

☒ Magnesium sulfate heptahydrate

☒ Acetone **Nacalai Tesque Inc. Catalog #00310-95**

☒ Hydrochloric Acid

☒ Copper (II) sulfate pentahydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #209198**

☒ Acetone **Nacalai Tesque Inc. Catalog #00310-95**

## Protocol materials

- ✕ Acetone **Nacalai Tesque Inc. Catalog #00310-95**
- ✕ Copper (II) sulfate pentahydrate **Bio Basic Inc. Catalog #CDB0063.SIZE.500g**
- ✕ Acetone **Nacalai Tesque Inc. Catalog #00310-95**
- ✕ Acetone **Nacalai Tesque Inc. Catalog #00310-95**
- ✕ Magnesium sulfate heptahydrate
- ✕ Magnesium sulfate heptahydrate
- ✕ Hydrochloric Acid
- ✕ Copper (II) sulfate pentahydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #209198**
- ✕ Acetone **Nacalai Tesque Inc. Catalog #00310-95**
- ✕ Hydrochloric Acid
- ✕ Distilled Water **Thermo Fisher Catalog #15230196**
- ✕ Magnesium sulfate heptahydrate
- ✕ Magnesium sulfate heptahydrate
- ✕ Magnesium sulfate heptahydrate
- ✕ Acetone **Nacalai Tesque Inc. Catalog #00310-95**
- ✕ Acetone **Nacalai Tesque Inc. Catalog #00310-95**
- ✕ Copper (II) sulfate pentahydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #209198**
- ✕ Hydrochloric Acid
- ✕ Acetone **Nacalai Tesque Inc. Catalog #00310-95**

## Safety warnings

- ⚠ This protocol requires the use of strong corrosive acids, corrosive and flammable solvents, and require the extraction of pigments that may stain clothing. Proper lab coat, eye protection, gloves in ventilation are required to conduct is chlorophyll extraction in titration protocol. Also, care must be taken to ensure that all materials used are disposed properly, as many of the chemicals may be hazardous to health and environment.



## Before start

In order to perform this chlorophyll extraction protocol you will need the following materials and chemicals:

the materials listed are based on one single sample, in must be multiplied based on the number of samples you would like to test

Two 20 ml (minimum) test tube

Ten 5 ml glass vials

A minimum of ten 200 micro-liters PCR tubes ( although more may be needed based on the accuracy necessary for the procedures performed)

A 100 to 1000 micro-liter adjustable pipette

30 disposable 1000 micro-liter pipette tips (Number varies based on need and mistakes)

1200 g Centrifuge

One 20 micron filtration filter paper

Mortar and Pestle

Transfer Pipettes (as needed for contamination prevention)

10 ml Graduated Cylinder

Test tube stirrers

Approximate protocol time: 6 hours total - 3 hours Preparation - 3 Hours Sample Sit time

100ul per PCR Tube Copper II sulfate Solution 60mg/ml Distilled H<sub>2</sub>O

As needed Distilled water


9 ml Acetone











1g Magnesium sulfate

2g of sample



## Preparation of Sample and Extraction







- 1 Add  2 g of sample **plant matter** or other test article to a mortar 1m

.
- 2 Add  1 g of  Magnesium sulfate heptahydrate crystals to the mortar. 1m
- 3 Grind the sample item and  Magnesium sulfate heptahydrate together until the leaves are ground up, the magnesium sulfate is completely dissolved into the plant liquid and a liquid starts to form. 5m
- 4 Transfer as much of the solid and liquid slurry to a 20ml test tube. 1m
- 5 Add  1 mL of  Acetone Nacalai Tesque Inc. Catalog #00310-95 to the mortar to wash up and collect as much of the sample as possible into the 20ml test tube. 2m
- 6 Add  8 mL  Acetone Nacalai Tesque Inc. Catalog #00310-95 to the 20ml test tube. This is the **raw mix**. 1m
- 7 Use a glass stirrer to stir the **raw mix**. Take care to press the plant slurry down to the base of the test tube, then agitate the slurry again to mix. Repeat this many times for  00:02:00 minutes. Ensure that the acetone is fully incorporated into the slurry. 2m
- 8 Let the **raw mix** stand for  00:05:00 minutes. 5m
- 9 Use a glass stirrer to stir the **raw mix** again for  00:01:00 minute. 1m




- 10 Using a **20 micron retention filtration paper** to filter the raw mix into another test tube. Ensure to use a new filter paper for each separate sample. This is the **refined mix**. 2m



## Dilution Phase

- 11 Add  1 mL  Hydrochloric Acid 1m
- to the **refined mix** and mix with clean glass stirrer for  00:01:00 minute to create the **Acid Activated Mix**.
- 12 Let the **Acid Activated Mix** stand for  00:02:00 minutes. 2m
- 13 Determine the resolution needed for the concentration. There resolution is the interval of dilution. For example, a 10% resolution will result in a dilution interval of 100%, 90%, 80% etc. The smaller the resolution the narrower the margin of error is. For this protocol, the resolution used will be 10%. 1m
- 14 Label ten 5ml vials in descending order from 100% to 10%. 1m
- 15 Add **Acid Activated mix** to each 5ml vial at volume based on the formula  $V_{ol} = 1000ul * P_{er}$  with  $P_{er}$  representing the percentage written on each vial. Use the percentage or decimal, not whole number. **Ex. 10% = .10** 5m
- 16 Add  Acetone Nacalai Tesque Inc. Catalog #00310-95 5m
- to each 5ml vial at volume based on the formula  $V_{ol} = 1000ul * (1 - P_{er})$  with  $P_{er}$  representing the percentage written on each vial. Use the percentage or decimal, not whole number. **Ex. 10% = .10**
- 17 Gently stir each vial by rotating them in a circular motion for  00:00:30 seconds each. 5m





- 18 Label The proper amount of PCR tubes corresponding to the percentages listed on the vials. 1m
- 19 Using an adjustable pipette, dispense  100  $\mu\text{L}$  of **Acid Activated Mix** from each vial to the corresponding PCR Tube. Ensure to change the pipette tip for every transfer in order to ensure contamination prevention. 5m

## Copper Transfer Phase

- 20 Prepare a **60mg $\text{CuSO}_4$  /ml $\text{H}_2\text{O}$**  concentration solution of  Copper (II) sulfate pentahydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #209198** to be referred to as the **copper (II) sulfate Solution**. 5m
- 21 Add  100  $\mu\text{L}$  **copper (II) sulfate Solution** to each PCR tube. 5m

## Test Phase

- 22 Close all PCR tubes and Centrifuge all tubes at  1200 x g for 2 minutes. 2m
- 23 In a PCR rack, let the samples stand for 6 hours 6h
- 24 Centrifuge all tubes at  1200 x g for 2 minutes. 2m

## Results Phase

- 25 Observe which PCR is the lowest percentage concentration not to have any trace of precipitated **Copper (II) Sulfate**. (A blue crystalline precipitate at the bottom. Take note of the percentage number to be used as the variable **P<sub>-</sub>**). 1m
- 26 Assign the highest percentage value vial to have visible **copper (II) Sulfate** precipitate the variable **P<sub>+</sub>**. 1m
- 27 Use this formula to calculate the mean value for the concentration in g/ml.



10m

$$\frac{(6 \times 10^{-4} * 893.51)/(P_- * 159.609) + (((6 \times 10^{-4} * 893.51)/(P_+ * 159.609) - (6 \times 10^{-4} * 893.51)/(P_- * 159.609)))/2}{2} = \text{Con}$$

- 28 Use this formula to calculate the margin of error of the concentration.

10m

$$((6 \times 10^{-4} * 893.51)/(P_+ * 159.609) - (6 \times 10^{-4} * 893.51)/(P_- * 159.609))/2 = \text{Moe}$$

- 29 Write the final results as **Con ± Moe**

1m