



Apr 22, 2025

Differential memory enrichment of cytotoxic CD4 T cells in Parkinson's disease patients reactive to α -synuclein - FACS of CELR2 expressing or non-expressing cells

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g7mn48gwz/v1>

Antoine Freuchet¹, Emil Johansson¹, cecilia¹, Alessandro Sette¹

¹La Jolla Institute for Immunology



Emil Johansson

La Jolla Institute for Immunology

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7mn48gwz/v1>

Protocol Citation: Antoine Freuchet, Emil Johansson, cecilia , Alessandro Sette 2025. Differential memory enrichment of cytotoxic CD4 T cells in Parkinson's disease patients reactive to α -synuclein - FACS of CELR2 expressing or non-expressing cells. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.q26g7mn48gwz/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: February 10, 2025

Last Modified: April 22, 2025

Protocol Integer ID: 119907

Keywords: cd4 tem cell, parkinson, facs of celr2, expressing celsr2, differential memory enrichment, subsequent scrnaseq analysis, celr2, celsr2, expressing cell, cells the purpose

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-000375

Abstract

The purpose of this protocol is to sort CD4 Tem cells either expressing or not expressing CELSR2, to allow subsequent scRNAseq analysis.

Differential memory enrichment of cytotoxic CD4 T cells in Parkinson's disease patients reactive to α -synuclein - FACS sorting

1 Setup

Thaw vial(s) in waterbath for ~1min30sec

Transfer into warm 10mL HR5 + 20uL of Benzonase. Centrifuge 400g, 5min, RT

Resuspend in 10mL of PBS-BSA 0.04% and count. Centrifuge 400g, 5min, RT

Resuspend at 20×10^6 cells/mL in PBS-BSA 1% (chilled). Add 1 ml

2 Staining

1. Centrifuge 400g, 5min, RT.

2. incubate cells with Fc blocking solution (90 μ L PBS and 10 μ L FC block) for 10 min at 4C.

3. Centrifuge 400g, 5min, RT.

4. Resuspend cells in staining solution. incubate for 30 in at 4C in the dark

	Antibody target	Fluorochrome	Amount per well (100 μ L per well)
	CD3	AF700	2.5
	CD45RA	PE-Cy7	2.5
	CD8a	BV650	1.25
	CCR7	PerCP-Cy5.5	1.25
	CD4	BV711	1
	CELSR2	Purified	0.25
	HASHTAG	-	1
	BV Buffer	-	10
	PBS-BSA 1%	-	80.25

1. Centrifuge 400g, 5min, RT.

2. Wash three times in PBS-BSA 1%.

3. Resuspend in 2 mL secondary antibody mix (1992 μ L PBS-BSA 1% and 8 μ L PE-F(ab')₂-Goat anti-Rabbit IgG (H+L) (A10542) per sample) and incubate for 30 min at 4C.

4. Wash three times in PBS-BSA 1%.

5. Resuspend in 2 mL of PBS-BSA 1%. Add DAPI prior to sorting to stain dead cells

6. Sort in 200 000 cells into PBS-FBS 20% pre-coated tubes