




Jan 15, 2024

Differentiation of Astrocytes from Human iPSC-derived NPCs

 [PLOS One](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.x54v9pb4mg3e/v1>



Balazs Szeky¹, Anita Feher¹, Melinda Zana¹, Radek Kucera^{2,3}, Jan Lochman^{4,5}, Andras Dinnyes^{1,6}

¹Biotalentum Ltd, 2100, Godollo, Hungary;

²Department of Pharmacology and Toxicology, Faculty of Medicine in Pilsen, Charles University, Alej Svobody 1655/76, 32300 Pilsen, Czech Republic;

³Department of Immunochemistry Diagnostics, University Hospital Pilsen, 30955 Pilsen, Czech Republic
⁵Laboratory of Neurobiology and Pathological Physiology, Institute of Animal Physiology and Genetics, Czech Academy of Sciences, Brno, Czech Republic;

⁴Department of Biochemistry, Faculty of Science, Masaryk University, Kotlářská 2, 61137 Brno, Czech Republic;

⁵Laboratory of Neurobiology and Pathological Physiology, Institute of Animal Physiology and Genetics, Czech Academy of Sciences, Brno, Czech Republic;

⁶Department of Physiology and Animal Health, Institute of Physiology and Animal Nutrition, Hungarian University of Agriculture and Life Sciences, 2100, Godollo, Hungary



Balazs Szeky

Biotalentum Ltd

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x54v9pb4mg3e/v1>

Protocol Citation: Balazs Szeky, Anita Feher, Melinda Zana, Radek Kucera, Jan Lochman, Andras Dinnyes 2024. Differentiation of Astrocytes from Human iPSC-derived NPCs . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.x54v9pb4mg3e/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 09, 2024

Last Modified: January 15, 2024

Protocol Integer ID: 93170

Keywords: Astrocyte, hiPSC, Neural progenitor cell, functional astrocytes with high purity, functional astrocyte, vivo astrocyte phenotype, astrocyte, differentiation of astrocyte, distinct subtypes of astrocyte, astrocyte phenotype, multifunctional glial cells of the central nervous system, neural progenitor cell, glial cell, multifunctional glial cell, neurodegenerative disorder, mechanisms of different CNS pathology, induced pluripotent stem cell, pluripotent stem cell, exacerbation of neurodegenerative disorder, central nervous system, different CNS pathology, essential roles in the metabolic support, neuron, drug metabolism

Funders Acknowledgements:

Gyorsítás

Grant ID: 2020-1.1.5-GYORSÍTÁSÁV-2021-00016

NDCi

Grant ID: 2019-2.1.7-ERA-NET-2022-00031

Ministry of Health, Czech Republic

Grant ID: NU20-09-00437

Abstract

Astrocytes are multifunctional glial cells of the central nervous system (CNS). They play essential roles in the metabolic support of neurons, synaptic transmission, maintaining blood-brain-barrier and immunological homeostasis. On the other hand, distinct subtypes of astrocytes play central roles in drug metabolism, neuroinflammation and in the exacerbation of neurodegenerative disorders. Understanding astrocyte-mediated mechanisms of different CNS pathologies requires a reliable in vitro model system recapitulating in vivo astrocyte phenotypes and functionalities. To provide such models, we developed an efficient, time- and cost-effective method for the generation of functional astrocytes with high purity from human induced pluripotent stem cell (iPSC)-derived neural progenitor cells (NPCs).

Materials

See the appropriate paragraphs of the method for product informations.

Safety warnings

 Biosafety level 1. Human cells should be cultured under aseptic conditions in a BSL2 laminar flow.

Ethics statement

Our research fully adhered to international ethical standards and guidelines for the establishment and maintenance of human iPSC lines, and iPSC-derived cell lines.

Preparation of Coated Plates and Solutions

1 Preparation of Coated Plates:

Neural progenitor cells can be grown on poly-L-ornithine/Laminin (POL) coated plates, whereas astroglial cells and astroglial progenitors must be seeded into Matrigel coated plates. Reagents required for preparing POL and matrigel plate coatings are described below.

1.1 Coating of Plates With POL

2h

	A	B	C	D
	Reagent	Manufacturer	Cat#	Dilution
	Poly-L-ornithine	Merck	P4957	1:5 (in PBS)
	Laminin	Sigma Aldrich	L2020	1 ug/cm ²

I. Prepare the desired plate format for POL coating

II. Prepare 5-times diluted poly-L-ornithine solution in PBS. Pipette the proper volume into the desired number of wells, incubate at  37 °C for  01:00:00 .

III. Prepare  1 µg /cm² laminin solution  On ice .

IV. Aspirate the Poly-L-ornithine from the plates. Wash the poly-L-ornithine-coated wells with PBS, then pipette the laminin solutions into the wells using chilled serological tips.

V. Seal the coated plates with parafilm, and store them at  4 °C overnight.

Incubate the POL-coated plates at  37 °C for at least  01:00:00 before use.

1.2 Coating of Plates With Matrigel

30m

hESC qualified Matrigel (Corning, #CLS354277) must be kept on ice to prevent condensation into gel phase. Matrigel stock solutions can be used only when the



solution becomes fully liquid. Matrigel must be diluted in ice-cold DMEM/F12 (Gibco, #11330032) medium.



I. Prepare the desired plate format for Matrigel coating

A	B
Plate format	Coating volume, mL/well
6-well plate	1.5 mL
24-well plate	0.5 mL
96-well plate	0.1 mL


II. Prepare Matrigel solution and DMEM/F12  On ice . Add 25 ul of matrigel into every 1.5 mL of DMEM/F12.

III. Pipette the prepared Matrigel working solution into the wells of the plate(s) using chilled serological tips. Incubate the plates at  37 °C for at least  00:30:00 before use. For extended storage, seal the plates with parafilm and store at  4 °C . Plates can be stored up to 14 days.

2 Preparation of Cell Culture Mediums

Neural maintenance medium (NMM) must be prepared for NPCs. NMM is completed with  10 ng/mL epidermal growth factor (EGF) and  10 ng/mL basic fibroblast growth factor (bFGF).

Astrocyte induction medium (AIM) is used for astroglial progenitor cell (APC) enrichment from NPCs.

Astrocyte maturation medium (AMM) is used for astrocyte maturation from APCs. AMM is completed with  20 ng/mL ciliary neurotrophic factor (CNTF).

2.1 Preparing NMM

Mix the components of NMM (listed in the table below) in a bottletop filter using filtered pipette tips and serological pipettes. Apply vacuum tube to the filter, turn on vacuum pump and sterile filter the solution. store NMM at  4 °C .



	A	B	C	D
	Component	Manufacturer	Stock concentration	Final Amount/Dilution
	DMEM/F12 with GlutaMax	Gibco		50% (%vol)
	Neurobasal Medium	Gibco		50% (%vol)
	B27	Gibco	50X	1X
	N2	Gibco	100X	1X
	GlutaMax	Gibco	100X	0.5X
	Non-Essential Amino Acid Solution	Merck	100X	1X
	Pen/Strep	Gibco	100X	1X
	EGF	Thermo Fisher	0.1 mg/mL	10 ng/mL
	bFGF	Thermo Fisher	0.1 mg/mL	10 ng/mL

2.2 Preparing AIM

Mix the components of AIM (listed in the table below) in a bottle-top filter using filtered pipette tips and serological pipettes. Apply vacuum tube to the filter, turn on vacuum pump and filter the solution. Store AIM at  4 °C .

	A	B	C	D
	Component	Manufacturer	Stock Solution	Final Amount/Dilution
	Astrocyte growth Medium	ScienCell		
	Foetal Bovine Serum (FBS)	ScienCell		2% (vol%)

	A	B	C	D
	Pen/Strep	Gibco	100X	1X

2.3 Preparing AMM

Mix the components of AMM (listed in the table below) in a bottletop filter using filtered pipette tips and serological pipettes. Apply vacuum tube to the filter, turn on vacuum pump and filter the solution. Store AMM at  4 °C .

	A	B	C	D
	Component	Manufacturer	Stock Solution	Final Amount/Dilution
	Astrocyte Growth Medium	ScienCell		
	CNTF	Peprtech	0.01 mg/mL	20 ng/mL
	Pen/Strep	Gibco	100X	1X

Neural Progenitor Cell Expansion

2w

3 Thawing NPCs

3.1 Prepare POL-coated plates and NMM before thawing NPCs. Prepare  5 mL NMM into a 15 mL Falcon tube.

3.2 Thaw one vial of frozen NPCs under sterile cell culture hood in sterile gloves. Pipette the cell suspension gently into the falcon tube containing the 5 mL NMM.

3.3 Centrifuge the cell suspension:  1000 rpm, 00:03:00 . Discard the supernatant.

3m



- 3.4 Resuspend the pellet in NMM, complete with EGF and bFGF. Pipette the suspension into the wells of the POL-coated plate
- 3.5 Put the plate into a cell culture incubator, ( 37 °C , 5% CO₂). Change medium next day to complete NMM.
- 3.6 Change medium every other day onwards. Passage NPCs upon reaching confluence.

4 Passaging NPCs

- 4.1 Prepare coated POL plate, complete NMM and Accutase at  Room temperature for passaging of NPCs.
- 4.2 Wash the cells with PBS, ( 2 mL /well for 6-well plates). Avoid scratching the cell monolayer.
- 4.3 Add Accutase, ( 1 mL /well for 6-well plates). Incubate at  37 °C for  00:03:00 . After the incubation, the rollup of the cells can be observed under the microscope. 3m
- 4.4 Add NMM ( 2 mL /well for 6-well plates) to stop the accutase digestion. Dissociate the monolayer by pipetting. Transfer the cell suspension into the Falcon tube.
- 4.5 Centrifuge the cell suspension:  1000 rpm, 00:03:00 . Discard the supernatant, except  200 µL . Resuspend the cells. 3m
- 4.6 Add completed NMM . Perform cell counting. Dilute the cell suspension with complete NMM to obtain 75.000-100.000 cells/cm² (depending on further experiments) for each well of the POL-coated plate, (1 million cells/well for 6-well plates).



4.7 Pipette the cell suspension into the POL-coated plates, ( 2 mL /well for 6-well plates). Distribute the cells by gentle shaking. Put the plates into a sterile cell culture incubator, ( 37 °C , 5% CO₂).

4.8 Change the medium next day to complete NMM. Change medium every other day onwards until the next passage.

5 Plating for Astrocyte Induction

5.1 By the second passaging of the NPCs, prepare Matrigel-coated plates, AIM and Accutase.

5.2 Passage and plate cells as described in steps 4.1-4.6.

5.3 Perform cell counting. Dilute the cell suspension with complete NMM to obtain 25.000-75.000 cells/cm² for each well of the Matrigel-coated plate.

5.4 Pipette the suspension into the matrigel coated plate (2mL/well for 6-well plates). Distribute the cells by gentle shaking. Put the plate into a sterile cell culture incubator, ( 37 °C , 5% CO₂).

Astroglial Progenitor Cell Induction

6 **Starting astrocyte induction at Day 0.** Next day after the second NPC passaging, replace completed NMM intro AIM, (2mL/well for 6-well plates). Change the medium every other day **at least for 21 days**. Passage the APCs upon reaching confluence.

7 Passaging APCs

7.1 Prepare Matrigel-coated plate, AIM and Accutase at  Room temperature for passaging of APCs.



- 7.2 Wash the cells with PBS, ( 2 mL /well for 6-well plates). Aspirate the PBS. Avoid scratching the cell monolayer.
- 7.3 Add Accutase, ( 1 mL /well for 6-well plates). Incubate at  37 °C for  00:03:00 -  00:07:00 until the cells begin to roll up.
- 7.4 Add AIM ( 2 mL /well for 6-well plates) to stop the accutase digestion. Dissociate the monolayer and add the 1mL cell suspension into a 15 mL Falcon tube.
- 7.5 Wash with  1 mL AIM and collect the rest of the cell suspension into the falcon tube.
- 7.6 Centrifuge the cell suspension:  1000 rpm,  00:03:00 . Discard the supernatant, except  200 µL . Resuspend the cells using 200 uL pipettes with filtered tips. 3m
- 7.7 Perform cell counting. Dilute the cell suspension with AIM to obtain 25.000-75.000 cells/cm² for each well of the Matrigel-coated plate, (0.5 million cells/well for 6-well plates).
- 7.8 Pipette the cell suspension into the Matrigel-coated plates. Distribute the cells by gentle shaking. Put the plates into a sterile cell culture incubator, ( 37 °C , 5% CO₂).
- 7.9 Change the medium next day to AIM. Change medium every other day onwards until the next passage.

Astrocyte Maturation

8 Initiating Astrocyte Maturation from APCs.

- 8.1 On the 21st day of differentiation, prepare Matrigel-coated plates, AIM, AMM and Accutase for the passaging of APCs.
- 8.2 Passage APCs as described in steps 7.1-7.9. Change medium to AMM completed with  20 ng/mL CNTF (2 mL/well for 6-well plates) next day, and every other day onwards **at least for 42 days**.



Protocol references

Tcw J, Wang M, Pimenova AA, Bowles KR, Hartley BJ, Lacin E, et al. An Efficient Platform for Astrocyte Differentiation from Human Induced Pluripotent Stem Cells. *Stem Cell Reports*. 2017;9(2):600-14.

Perriot S, Canales M, Mathias A, Du Pasquier R. Differentiation of functional astrocytes from human-induced pluripotent stem cells in chemically defined media. *STAR Protoc*. 2021;2(4):100902.

Chandrasekaran A, Avci HX, Ochalek A, Rösingh LN, Molnár K, László L, et al. Comparison of 2D and 3D neural induction methods for the generation of neural progenitor cells from human induced pluripotent stem cells. *Stem Cell Res*. 2017;25:139-51.