

Jan 10, 2025

Differentiation of bone marrow-derived macrophages

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.81wgbr2e3lpk/v1>

Susanne Herbst^{1,2}, Patrick Lewis^{1,2}

¹RVC;

²The Michael J. Fox Foundation for Parkinson's Research (MJFF) and the Aligning Science Across Parkinson's (ASAP) Initiative



Susanne Herbst

RVC

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgbr2e3lpk/v1>

Protocol Citation: Susanne Herbst, Patrick Lewis 2025. Differentiation of bone marrow-derived macrophages. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.81wgbr2e3lpk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 09, 2025

Last Modified: January 10, 2025

Protocol Integer ID: 117981

Keywords: ASAPCRN, BMDM, macrophages, differentiation, derived macrophage, macrophages this protocol, differentiation of bone marrow, bone marrow, using gm, differentiation factor

Abstract

This protocol describes the differentiation of bone marrow-derived macrophages using GM-CSF as a differentiation factor.

Materials

Plastic consumables:

- 20 ml syringe and 27G needles
- Falcon tubes
- 100 um cell strainer (eg. 130-098-463, Miltenyi Biotec)

Red blood cell lysis buffer, eg 00-4333-57, eBioscience

Media preparation:

- DMEM (eg. 31966047, Gibco, ThermoFisher) +
- 10 % v/v heat-inactivated FCS
 - 50 ng/ml murine GM-CSF (eg. 576306, BioLegend)

PBS/EDTA:

- PBS containing 2 mM EDTA

Troubleshooting



Differentiation

- 1 Collect femurs and clean bones of flesh.
- 2 Cut bones open on both ends and flush bone marrow into 50 ml Falcon tube with DMEM using a syringe.
- 3 Spin down cells at 300 x g for 5 min and resuspend pellet in red blood cell lysis buffer. Incubate for 10 min.
- 4 Add an equal amount of DMEM and filter cell suspension through a 100 um cell strainer.
- 5 Spin down cells at 300 x g for 5 min and resuspend cell pellet in DMEM containing 10 % FCS and 50 ng/ml murine GM-CSF.
- 6 Dispense cell suspension into 100mm Petri dishes containing a final volume of 10-12 ml of DMEM/ 10 % FCS/ 50 ng/ml murine GM-CSF.
- 7 Incubate cells at 37 °C, 5% CO₂.
- 8 At day 3 of differentiation: Add 5 ml of DMEM containing 10 % FCS and 50 ng/ml murine GM-CSF
- 9 At day 7 of differentiation: Harvest cells

Note

If there are many floating cells, cells can be split and differentiated further and initial cell densities can be decreased in the future.

Cell harvest

- 10 Remove medium and wash plates once in PBS.
- 11 Add 3 ml of PBS/EDTA to each Petri dish and place dishes in fridge for 10 to 15 min.



- 12 Harvest cells by flushing the surface of the Petri dish. Pool cells according to genotype and mouse ID.
- 13 Spin down cells at 300 x g for 5 min.
- 14 Resuspend cells in DMEM containing 10 % FCS and seed for experiments or freeze.

Freezing and thawing of BMDMs

- 15 To freeze cells, supplement DMEM/10% FCS with 10 % DMSO. Freeze over night at -80 °C and transfer to liquid nitrogen the following day.
- 16 To thaw cells, spin down cells to remove DMSO and resuspend pellet in DMEM/10% FCS containing 50 ng/ml murine GM-CSF. Plate cells accordingly and let recover for 24-48 hrs. On the day of the experiment, the GM-CSF can be omitted from the medium.