



Jun 17, 2025

Differentiation of hiPSC-derived microglia-like cells

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gpdrp5gzp/v1>

Alexis Penverne¹, James Evans²

¹Gandhi Lab; ²University College London, University of London

Alexis Penverne: University College London, University of London;



Alexis Penverne

Gandhi Lab, University College London, University of London

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gpdrp5gzp/v1>

Protocol Citation: Alexis Penverne, James Evans 2025. Differentiation of hiPSC-derived microglia-like cells. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.dm6gpdrp5gzp/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: April 03, 2025

Last Modified: June 17, 2025

Protocol Integer ID: 126085

Keywords: like cells from hipsc, differentiation of hipsc, like cells protocol for generation, hipsc, like cells protocol, derived microglia, microglia

Abstract

Protocol for generation of microglia-like cells from hiPSCs, based on protocol from Garcia-Reitboeck, Pablo et al., 2018.



Day 0 - Embryoid Body generation

8m

- 1 hiPSCs are maintained on Geltrex (ThermoFisher) in mTeSR (STEMCELL), and passaged using [M] 0.5 Mass Percent EDTA.
- 2 Once hiPSCs are healthy and ~80% confluent, incubate with 1 mL Accutase (STEMCELL) for 00:05:00 at 37 °C Dissociate using 1 mL of PBS and resuspend in a Falcon tube in 5 mL of PBS.
- 3 Centrifuge at 300xg for 00:03:00 and aspirate supernatant. Resuspend pellet in 1 mL of **Embryoid Medium** (1 mL of media per well of hiPSC initially dissociated).

Note

How to prepare **Embryoid Medium**:

- mTeSR (STEMCELL)
- [M] 10 micromolar (μ M) Rock Inhibitors (Tocris)
- 50 ng/mL BMP-4 (Peprotech)
- 50 ng/mL VEGF (Peprotech)
- 20ng/mL SCF (Peprotech)

- 4 Count cells and dilute in Embryoid medium to have suspension at 100,000 cells/mL.
- 5 Add 100 μ L of Embryoid Medium to each well of a low-attachment U-bottom 96-well plate.
- 6 Using a multi-channel pipette, add 100ul of cell suspension in each well of the 96-well plate.

**Note**

A full-well of embryoid bodies is preferable to ensure sufficient quantity of microglia-like cells downstream, therefore ~1 million hiPSCs are needed per induction.

7 Centrifuge the plate at 800rpm for 00:03:00 and transfer to 37 °C 5% CO2 incubator.

3m

8 After 48:00:00 , supplement each well with 50 μ L of fresh Embryoid media.

Day 4 - Embryoid Body Suspension and Myeloid differentiation

9 On day 4 after the start of the induction, collect every single Embryoid Body from each well and transfer in a single 15mL Falcon using a P1000 pipette and wide-bore tips.

10 Let EBs settle naturally at the bottom of the tube and aspirate Embryoid media, and Add 5mL of **Myeloid Differentiation media**.

Note

How to prepare **Myeloid Differentiation Medium**:

- X-Vivo 15 (Lonza)
- 1 % volume Glutamax (ThermoFisher)
- 0.5 % volume Penicilin-Streptomycin (ThermoFisher)
- 50 micromolar (μ M) 2-Mercaptoethanol (ThermoFisher)
- 100ng/mL MSCF (Peprotech)
- 25ng/mL IL3 (Peprotech)

11 Add 10 mL of Myeloid differentiation media to a T75 flask, and transfer the 5 mL with the EBs inside the flask for a total of 15 mL of media.



- 12 For the following 2 weeks (until day 18), change media with fresh Myeloid Differentiation media. Change by removing ~  10 mL of media from the flask and replacing it with fresh media. Be careful to move the plate gently: the EBs need to stick to the flask as this will improve the microglia-like cells yield.

Day 18 onward - Microglia Harvest

- 13 Starting on day 18 microglia-like cells can be harvested weekly by taking 10-12mL of media from T75 flask inside a 15mL Falcon. Replace with fresh myeloid differentiation media in flask.
- 14 Centrifuge harvested cells in media at 300xg for  00:03:00 . Resuspend pellet in  1 mL of **Maturation Medium**.

Note

How to prepare **Maturation Medium**:

- X-Vivo 15 (Lonza)
-  1 % volume Glutamax (ThermoFisher)
-  0.5 % volume Penicilin-Streptomycin (ThermoFisher)
- 100ng/mL MSCF (Peprotech)

- 15 Filter resuspended cells through a 40-micron filter, and count. Dilute to a concentration of ~300,000 cells/mL in Maturation media.
- 16 Plate microglia-like cells directly onto plasticware (well-plates, imaging chambers...) in maturation media at desired concentration.

Note

The first harvest tends to give a low yield of cells. Typically, from one induction it is possible to obtain roughly 5-6 weekly harvests before the yields decrease.

- 17 Perform experiments on microglia-like cells after 7 days from harvest.



Protocol references

Garcia-Reitboeck, Pablo et al. "Human Induced Pluripotent Stem Cell-Derived Microglia-Like Cells Harboring TREM2 Missense Mutations Show Specific Deficits in Phagocytosis." *Cell reports* vol. 24,9 (2018): 2300-2311. doi:10.1016/j.celrep.2018.07.094