



Dec 12, 2025

Differentiation of Induced Pluripotent Stem Cells (iPSCs) into Microglia-like cells (iMG)

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov11rw2vr2/v1>

Elisa Perciballi¹, Micheal John Dolan², Martine Therrien¹, Beth Stevens³

¹University of California, Davis; ²Smurfit Institute of Genetics School, Trinity College Dublin;

³HHMI, Broad Institute, Boston Children's Hospital



Martine Therrien

University of California, Davis

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov11rw2vr2/v1>

Protocol Citation: Elisa Perciballi, Micheal John Dolan, Martine Therrien, Beth Stevens 2025. Differentiation of Induced Pluripotent Stem Cells (iPSCs) into Microglia-like cells (iMG). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.ewov11rw2vr2/v1>

**Manuscript citation:**

Dolan, M.-J., Therrien, M., Jereb, S., Kamath, T., Gazestani, V., Atkeson, T., Marsh, S.E., Goeva, A., Lojek, N.M., Murphy, S., et al. (2023). Exposure of iPSC-derived human microglia to brain substrates enables the generation and manipulation of diverse transcriptional states in vitro. *Nat. Immunol.* 24, 1382–1390. <https://doi.org/10.1038/s41590-023-01558-2>.

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 06, 2025

Last Modified: December 12, 2025

Protocol Integer ID: 234321

Keywords: stem cells, microglia, derived microglia, role of microglia, microglia, microglial state, induced pluripotent stem cell, differentiation of induced pluripotent stem cell, immune alterations during development, immune alteration, human cell, response to specific brain challenge, use of human cell, lentiviral transduction, amenable to lentiviral transduction, specific brain challenge

Funders Acknowledgements:

Alzheimer's Association

Cure Alzheimer's Fund

Howard Hughes Medical Institute

Stanley Center for Psychiatric Research

Abstract

iPSC-derived microglia are an important tool for studying the role of microglia and neuro-immune alterations during development, adulthood, and aging, both in health and disease. The use of human cells enables the identification of molecular mechanisms driven by an individual's genetic susceptibility that would not be possible with animal models. Here, we described a step-by-step protocol used in our previous manuscript (Dolan, Therrien, et al 2023). This protocol is used to induce a microglial state that arises in response to specific brain challenges and is amenable to lentiviral transduction.

Materials

✂ ACCUTASE™ 100 mL STEMCELL Technologies Inc. Catalog #7920 Equipment

- Incubator for normoxic conditions (20% O₂, 5% CO₂, 37°C)
- Humidified tri-gas incubator for hypoxic conditions (5% O₂, 5% CO₂, 37°C)
- Centrifuge
- Cell counting equipment

iPSCs maintenance list of reagents

A	B
Reagent	Company; Cat #
PBS 1X	Thermo Fisher; 14190250
Matrigel	VWR; 47743-722
Essential E8 medium	Thermo Fisher; A1517001
Cryostor	Stem Cell Technologies; 100-1061
Accutase	Stem Cell Technologies; 07920

Microglia differentiation media composition and reagents

A	B	C
iHPC Base Growth Medium (first 10 days)		
Component	Final concentration	Company; Cat #
IMDM	50%	Thermo Fisher; 31980097
F12	50%	Thermo Fisher; 11765062

A	B	C
ITSG-X	2% v/v	Thermo Fisher; 51500-056
L-ascorbic acid 2-Phosphate	64 µg/ml	Sigma-Aldrich; A8960-5G
monothioglycerol	400 µM	Sigma-Aldrich; M1753-100ML
PVA	10 µg/ml	Sigma-Aldrich; P8136-250G
GlutaMax (1X)	1X	Thermo Fisher; 35050079
Chemically defined lipid concentrate (1X)	1X	Thermo Fisher; 11905031
Non-essential amino acids (NEAA) (1X)	1X	Thermo Fisher; 11140050

iHPC base media: All reagents are combined and filtered in .22µM bottle. Base media is used within 2 weeks. L-ascorbic acid 2-Phosphate and PVA solutions are used within a week.

A	B	C
Day 0 growth factors		
Component	Final concentration	Company; Cat #
FGF2	50 ng/ml	Thermo Fisher; PHG0261
BMP4	50 ng/ml	Thermo Fisher; PHC9534
Activin-A	12.5 ng/ml	Thermo Fisher; PHC9564
Rock Inhibitor (RI)	1 µM	Selleckchem; s1049
Lithium Chloride (LiCl)	2 mM	Sigma-Aldrich; L7026

Day 0 iHPC media: iHPC base media + growth factors listed. Growth factors are aliquoted in single use aliquots and store at -80C or -20C until expiration date following manufacturer instructions.



	A	B	C
	Day 2 growth factors		
	Component	Final concentration	Company; Cat #
	FGF2	50 ng/ml	Thermo Fisher; PHG0261
	VEGF	50 ng/ml	Peprotech; 100-20-50UG

Day 2 iHPC media: iHPC base media + growth factors listed. Growth factors are aliquoted in single use aliquots and store at -80C or -20C until expiration date following manufacturer instructions.

	A	B	C
	Day 4 growth factors		
	Component	Final concentration	Company; Cat #
	FGF2	50 ng/ml	Thermo Fisher; PHG0261
	VEGF	50 ng/ml	Peprotech; 100-20-50UG
	TPO	50 ng/ml	Thermo Fisher; 300-18-50UG
	SCF	10 ng/ml	Thermo Fisher; PHC2115
	IL-6	50 ng/ml	Thermo Fisher; 200-06
	IL-3	10 ng/ml	Thermo Fisher; 200-03-10UG

Day 4-10 iHPC media: iHPC base media + growth factors listed. Growth factors are aliquoted in single use aliquots and store at -80C or -20C until expiration date following manufacturer instructions.

	A	B	C
	iMG differentiation basal medium		

	A	B	C
	Component	Final concentration	Company; Cat #
	DMEM	50%	Thermo Fisher; 10569044
	F12	50%	Thermo Fisher; 11765062
	ITS-G	2x v/v	Thermo Fisher; 41400045
	B-27	2x v/v	Thermo Fisher; 17504044
	N2	0.5x	Thermo Fisher; 17502048
	Monothioglycerol	200 μ M	Sigma-Aldrich; M1753-100ML
	GlutaMax (1X)	1X	Thermo Fisher; 35050079
	Non-essential amino acids (NEAA) (1X)	1X	Thermo Fisher; 11140050
	Insulin	5 μ g/ml	Sigma-Aldrich; I2643

Mg base media: All reagents are combined and filtered in .22uM bottle. Base media is used within 2 weeks. B27 and N2 are frozen in single used aliquots.

	A	B	C
	Day 10 growth factors		
	Component	Final concentration	Company; Cat #
	M-CSF	25 ng/ml	Thermo Fisher; 300-25-100
	IL-34	10 ng/ml	Thermo Fisher; 200-34-250UG
	TGF β	50 ng/ml	Thermo Fisher; 100-21-100UG



Day 10 Mg media: Mg base media + growth factors listed. Growth factors are aliquoted in single use aliquots and store at -80C or -20C until expiration date following manufacturer instructions.

	A	B	C
	Day 30 growth factors		
	Component	Final concentration	Company; Cat #
	CD200	100 ng/ml	Thermo Fisher; 103872-010
	CX3CL1	100 ng/ml	Thermo Fisher; 300-31
	M-CSF	25 ng/ml	Thermo Fisher; 300-25-100
	IL-34	10 ng/ml	Thermo Fisher; 200-34-250UG
	TGFβ	50 ng/ml	Thermo Fisher; 100-21-100UG

Day 10 Mg media: Mg base media + growth factors listed. Growth factors are aliquoted in single use aliquots and store at -80C or -20C until expiration date following manufacturer instructions.

Materials for microglia cell growth

	A	B
	Product name	Company; Cat #
	6-well Ultra-Low Attachment plate	Corning Life Sciences; 3471
	Cell lifter	Life Sciences; 76036-004



Protocol materials

- ✕ Essential 8 Life Technologies Catalog #A1517001
- ✕ TeSR™-E8™ Kit for hESC/hiPSC Maintenance 1 Kit **STEMCELL Technologies Inc. Catalog #5990**
- ✕ Corning® Matrigel® Basement Membrane Matrix, Phenol Red-free, LDEV-free, 10 mL **Corning Catalog #356237**
- ✕ CryoStor CS10-100ml **STEMCELL Technologies Inc. Catalog #07930**
- ✕ ACCUTASE™ 100 mL **STEMCELL Technologies Inc. Catalog #7920**
- ✕ DPBS, no calcium, no magnesium **Thermo Fisher Scientific Catalog #14190144**
- ✕ Y-27632 (ROCK inhibitor) (Selleck Chemicals Catalog # S1049) **Selleckchem Catalog #S1049**
- ✕ Costar® 6-well Clear Flat Bottom Ultra-Low Attachment Multiple Well Plates **Corning Catalog #3471**
- ✕ Fisherbrand™ Cell Lifters **Thermo Fisher Scientific Catalog #08-100-240**

iPSC maintenance

1 Maintenance of stem cells

Successful differentiation is achieved with high quality stem cells. It is recommended to split iPSCs at least 2 times after thawing before starting differentiation and check frequently at their morphology to avoid spontaneous differentiations that may interfere with mesodermal induction. Make sure iPSCs are growing in

 **Essential 8 Life Technologies Catalog #A1517001** or

 **TeSR™-E8™ Kit for hESC/hiPSC Maintenance 1 Kit STEMCELL Technologies Inc. Catalog #5990**

otherwise change accordingly.

Generation of HPC

47m

2 Day -1: Generation of embryoid bodies

2.1 Detach cells

15m

- Aspirate the culture media and wash with 1X PBS.
- Dissociate iPSCs by incubating with

 **ACCUTASE™ 100 mL STEMCELL Technologies Inc. Catalog #7920**

 **00:10:00** (time may vary depending on different cell lines) at  **37 °C** , 5%

CO₂. Lightly tap the side of the plate to dislodge the cells.

- Collect the entire plate volume and transfer the cell suspension to a 15 ml conical tube.
- Wash the plate again with

 **DPBS, no calcium, no magnesium Thermo Fisher Scientific Catalog #14190144**

to collect cells left behind and transfer to the same 15 ml tube.

- Centrifuge the cells at  **300 x g, 00:05:00** . Aspirate supernatant to waste and resuspend the pellet to triturate to create a single cell suspension using

 **Essential 8 Life Technologies Catalog #A1517001** with **[M] 10 micromolar (μM)**

 **Y-27632 (ROCK inhibitor) (Selleck Chemicals Catalog # S1049) Selleckchem Catalog #S1049**



2.2 Count the cells and adjust cell density to 5×10^5 per well in a

 Costar® 6-well Clear Flat Bottom Ultra-Low Attachment Multiple Well Plates **Corning Catalog #3471**

in

 Essential 8 **Life Technologies Catalog #A1517001**

+

 Y-27632 (ROCK inhibitor) (Selleck Chemicals Catalog #S1049) **Selleckchem Catalog #S1049**

(2 ml in a 6-well plate).

- Count cells using Trypan blue. Viability should be at least 90%.
- Plate evenly on growth surface and begin normoxic culture conditions (20% O₂, 5% CO₂, 37°C).

Note

Because of the different proliferation rates among cell lines, it is recommended to plate different concentrations (2×10^5 and 8×10^5) and start differentiation with the one with 50-100 embryoid bodies (EBs).

3 Day 0: Change media

5m

- Carefully collect the entire plate medium with EBs in suspension by using a 5 ml pipette to not dislodge them, and transfer to a 15 ml tube.
- Centrifuge EBs at  100 x g, 00:05:00 . Aspirate supernatant to waste and resuspend EBs in iHPC Base Growth medium supplemented with  2 mL Day 0 growth factors in a new low-attachment plate well.
- Plate evenly on growth surface and begin hypoxic culture conditions by incubating cells in a humidified tri-gas incubator with 5% O₂, 5% CO₂,  37 °C .

4 Day 2: Making HPC media and change media

5m

- Carefully collect the entire plate medium with EBs in suspension by using a 5 ml pipette to not dislodge them, and transfer to a 15 ml tube.
- Centrifuge EBs at  100 x g, 00:05:00 . Aspirate supernatant to waste and scroll the tune on the grill of the hood to dislodge the pellet. Resuspend EBs in  2 mL iHPC Base Growth medium supplemented with Day 2 growth factors in a new low-attachment plate well.



- Plate evenly on the growth surface and incubate under hypoxic culture conditions (5% O₂, 5% CO₂, 37 °C).

5 Day 4: Change media

5m

- Carefully collect the entire plate medium with EBs in suspension by using a 5 ml pipette to not dislodge them, and transfer to a 15 ml tube..
- Centrifuge EBs at 100 x g, 00:05:00 . Aspirate supernatant to waste and scroll the tube on the grill of the hood to dislodge the pellet. Resuspend EBs in 2 mL iHPC Base Growth medium supplemented with Day 4 growth factors in a new low-attachment plate well.
- Plate evenly on growth surface and incubate under normoxic culture conditions (20% O₂, 5% CO₂, 37 °C).

6 Day 6: Feed cells

Add 1 mL of iHPC Base Growth medium supplemented with Day 4 growth factors.

7 Day 8: Feed cells

Add 1 mL of iHPC Base Growth medium supplemented with Day 4 growth factors.

8 Day 10: Collect HPC

7m

- Carefully collect the entire plate medium with EBs in suspension and transfer to a 15 ml tube.
- Let the EBs precipitate by gravity (around 00:02:00).
- Collect the supernatant containing hematopoietic progenitor cells (hpc) in single cells in a tube and filter through a 40 µm screen (to avoid any residual EBs).
- Centrifuge the single cell suspension at 300 x g, 00:05:00 Aspirate supernatant to waste and resuspend the pellet in Mg differentiation basal medium supplemented with Day 10 growth factors.
- Count the cells using Trypan blue

**Note**

Quality control: To assess that differentiation is successful, on Day 10, it is recommended to run a quality control to check viability, and the expression of CD45 and CD43 markers via flow cytometry. Viability should be > 70%, CD43 positive cells > 95% and CD45 positive cells > 50%

Note

Cells can be frozen at this time in

 CryoStor CS10-100ml **STEMCELL Technologies Inc. Catalog #07930** . We recommend freezing at least 300 000 cells/aliquot so they can be thawed in 1 well of a 6-well plate

8.1 Cell replating to continue differentiation

- Resuspend HPC in Mg differentiation basal media with Day 10 growth factors
- Adjust cell density to 3×10^5 per well in a

 Corning® Matrigel® Basement Membrane Matrix, Phenol Red-free, LDEV-free, 10 mL **Corning Catalog #356237**

coated 6-well plate.

- Incubate under normoxic culture conditions (20% O₂, 5% CO₂,  37 °C).

9 Every 2 days: Feed cells

Add 1 ml of Mg differentiation basal medium supplemented with Day 10 growth factors

10 Day 22: Collect, count and replate cells**10.1** Collect and count cells

- Collect the entire plate medium and transfer to a 15 ml tube (do not aspirate, as many cells might still be in suspension).
- Collect early iMGL by adding

 DPBS, no calcium, no magnesium **Thermo Fisher Scientific Catalog #14190144**

to the well and let sit for  00:05:00 at room temperature.

10m



- Gently detach cells with a

 Fisherbrand™ Cell Lifters **Thermo Fisher Scientific Catalog #08-100-240** .

Wash the well with

 DPBS, no calcium, no magnesium **Thermo Fisher Scientific Catalog #14190144**

and transfer to a new 15 ml tube.

- Centrifuge both tubes (one containing plate medium and one PBS) at

 300 x g, 00:05:00

- From the tube with PBS, aspirate supernatant to waste. From the tube with conditioned Mg medium, collect the supernatant, and transfer to a new tube, leaving 1 ml to cover the pellet. Resuspend the pellet with the remaining 1 ml and transfer this cell suspension into the PBS tube.
- Combine and gently triturate the two pellets together.
- Count using Trypan blue

10.2 Replate cells

- Adjust cell density to 3×10^5 per well in a

 Corning® Matrigel® Basement Membrane Matrix, Phenol Red-free, LDEV-free, 10 mL **Corning Catalog #356237**

coated 6-well plate.

- Plate cells in 50% old conditioned medium + 50% fresh Mg differentiation basal medium supplemented with Day 10 growth factors.

Note

For high-proliferating cell lines, it might be necessary to split cells at Day 20 instead.

Microglia generation and maturation

5m

- 11 Starts from either freshly collected iHPC or from frozen iHPC vials. If cells were frozen, thaw rapidly at  37 °C water bath, resuspend cells in warm DMEM-F12, centrifuge at

 300 x g, 00:05:00 and count to plate at the appropriate density.

5m

- 12 Every 2 days: Feed cells

Add 1 ml of Mg differentiation basal medium supplemented with Day 10 growth factors

13 Day 30: Collect, count and replate cells

Note

Cells at Day 30 must be plated in the final experimental vessel.

13.1 Collect and count cells

- Collect the entire plate medium and transfer to a 15 ml tube (do not aspirate, as many cells might still be in suspension).
- Collect early iMGL by adding
 - ⊗ DPBS, no calcium, no magnesium **Thermo Fisher Scientific Catalog #14190144**to the well and let sit for  00:05:00 at room temperature.
- Gently detach cells with a
 - ⊗ Fisherbrand™ Cell Lifters **Thermo Fisher Scientific Catalog #08-100-240**Wash the well with
 - ⊗ DPBS, no calcium, no magnesium **Thermo Fisher Scientific Catalog #14190144**and transfer to a new 15 ml tube.
- Centrifuge both tubes (one containing plate medium and one PBS) at
 - ⦿ 300 x g, 00:05:00
- From the tube with PBS, aspirate supernatant to waste. From the tube with conditioned Mg medium, collect the supernatant, and transfer to a new tube, leaving 1 ml to cover the pellet. Resuspend the pellet with the remaining 1 ml and transfer this cell suspension into the PBS tube.
- Combine and gently triturate the two pellets together.
- Count using Trypan blue.

13.2 Replate cells

- Adjust cell density to 3×10^5 per well in a
 - ⊗ Corning® Matrigel® Basement Membrane Matrix, Phenol Red-free, LDEV-free, 10 mL **Corning Catalog #356237**coated 6-well plate.
- Plate cells in 50% old conditioned medium + 50% fresh Mg differentiation basal medium supplemented with Day 30 growth factors.

14 Every 2 days: Feed cells

Add 1 ml of Mg differentiation basal medium supplemented with Day 30 growth factors.



15 Day 40: Collect cells for desired experiments.