




Jun 28, 2023

Differentiation of mature neurons from mouse NPCs

DOI

dx.doi.org/10.17504/protocols.io.6qpvr4oo2gmk/v1



Sonia Cruciani¹, Eva Maria Novoa¹

¹Center for Genomic Regulation (CRG)



Sonia Cruciani

CRG

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvr4oo2gmk/v1>

Protocol Citation: Sonia Cruciani, Eva Maria Novoa 2023. Differentiation of mature neurons from mouse NPCs. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.6qpvr4oo2gmk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: September 21, 2022

Last Modified: June 28, 2023

Protocol Integer ID: 70334

Keywords: differentiation of mature neuron, mature neuron, mouse npc, neuron, differentiation

Funders Acknowledgements:

MICINN Spanish Grant Funding

Grant ID: PID2021-128193NB-I00

Abstract

This is a step by step protocol used to differentiate mouse NPCs to mature neurons.

Materials

- Mouse NPCs: vial of mNPCs derived from mES mettl3#6Cre from TEF
- Neurobasal (Thermo scientific, cat. no. 21103-049)
- DMEM/F12 (Thermo scientific, cat. no. 11320-033)
- N2 supplement (Gibco, cat. no. 17502-048)
- Phosphate-buffered saline (PBS; Life technologies, cat. no. 14190-169)
- EDTA 0.5M pH 8
- B27 minus vitamin A (Gibco, cat. no. 11500446)
- Non-essential amino acids (Thermo scientific, cat no.11140-035)
- 2-mercaptethanol (Thermo scientific, cat no.31350-010)
- Glutamax (life technologies, cat no. 35050-061)
- P/S (Thermo scientific, cat no. 15140-122)
- Laminin (Sigma cat no. L2020)
- Poli-L-Ornithine solution (Sigma cat no. A004-C)
- EGF (Gibco, cat no. PMG3043)
- FGF2 (aka bFGF) (Stem cell technologies, cat no. 78003)
- DAPT (Selleck,cat no. S2215)
- BDNF (Sigma,cat no. B3795)
- FBS
- DMSO

•N2B27 medium:

Component	Volume for 50ml
Neurobasal	24 mL
DMEM/F12	24 mL
P/S	0.5 mL
MEM Non-Essential Amino Acids Solution	0.5 mL
GlutaMAX™- 100X	0.5 mL
B-27 without vitamin A	0.5 mL
N2	0.25 mL
β-mercaptoethanol, 1000X. Sigma M-7522	50 uL

This medium could be stored 1 month at 4 C



●**NPC medium:**

N2B27 + EGF 10 ng/ml + FGF2 10 ng/ml

●**differentiation medium:**

N2B27 + BDNF 40 ng/ml + DAPT 10 uM

Troubleshooting

Coating

- 1 -Dilute Poly-L-Ornithine (PLO) 1:6 in H₂O
- 2 -Fully coat the wells with diluted PLO (i.e. 500 ul per well in a 12WP)
- 3 -Leave in the incubator for minimum 2h (also overnight works)
- 4 -Aspirate the PLO and wash three times with sterile H₂O
- 5 -Dilute Laminin 1:500 in sterile H₂O
- 6 -Fully coat the wells with diluted Laminin (i.e. 500 ul per well in a 12WP)
- 7 -Leave in the incubator for minimum 2h (also overnight works)

Plating - Day -1

- 8 -count the desired number of cells: Plate ~100.000 cells on a coverslip, ~150.000 cells in a 24WP well, ~250.000 cells in a 12WP well.
- 9 -seed NPCs for the final differentiation in NPC medium

Day 0

- 10 -aspirate all NPC medium
- 11 -immediately add differentiation medium not directly on the bottom of the well but making it go through the walls to avoid disturbing the cells.



Day 2-14

- 12 -Change medium (remove half volume and add half volume of new medium) avoiding the cells' contact with air every 2-3 days. If at diff day 2-3 you notice overgrowth of glial cells, add AraC 1uM to the media in order to stop proliferation.
Cells can be kept in colture for >14 days, but they are already MAP2+ and TUJ1+ at day 7.

Tamoxifen treatment

- 13 -Add 2.5 uM Tamoxifen to the medium, and the corresponding volume of metOH as control. Incubate for at least 6 days.