

Feb 26, 2019

Differentiation of NPC into cortical neurons

 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.x9bfr2n

Celeste M M. Karch¹, Rita Martinez¹, Jacob Marsh¹


¹Washington University in St Louis

Neurodegeneration Method Development Community
Tech. support email: ndcn-help@chanzuckerberg.com



Celeste M M. Karch

Washington University in St Louis



OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.x9bfr2n

Protocol Citation: Celeste M M. Karch, Rita Martinez, Jacob Marsh 2019. Differentiation of NPC into cortical neurons.
protocols.io <https://dx.doi.org/10.17504/protocols.io.x9bfr2n>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 17, 2019

Last Modified: February 26, 2019

Protocol Integer ID: 20483



Attachments



IPSC CORTICAL

DIFFER...

179KB

Guidelines

This protocol is part of the IPSC CORTICAL DIFFERENTIATION collection.

This method should be performed using sterile technique.

Materials

Please refer to the attached full manuscript for required materials.

Safety warnings




Please refer to the SDS (Safety Data Sheet) for information about hazards, and to obtain advice on safety precautions.

Before start

IMPORTANT: To generate cortical neurons, NPCs must be at passage 4 or lower at the time of plating for terminal differentiation. NPCs that are beyond passage 4 may exhibit higher densities of astrocyte contamination or inefficient neuronal differentiation.



- 1 Ensure cells are plated at an appropriate density (e.g.: 150K/well in 12-well plate; 75-50K/well for 48 well plates; 30K/well for 8-well chamber slides) in NIM on pre-coated PLO/laminin plates.

- 2 After  24:00:00 , replace with cortical neuron differentiation medium (Neurobasal medium, 1x B27, 20 ng/mL BDNF, 20 ng/mL GDNF, 0.5mM cAMP, 1% Glutamax, 1% penicillin/streptomycin)

- 3 Feed cells every 2-3 days for 30 days.