

Oct 31, 2025

Digital PCR Analysis for Mitochondrial DNA Integrity and Copy Number (Absolute Q MAP16 Format)

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gpm1w5gzp/v1>

Jenny Ghelfi¹, Soraya Hezzaz¹, Sylvie Delcambre¹, Anne Grünewald^{1,2}

¹Luxembourg Centre for Systems Biology (LCSB), University of Luxembourg, Luxembourg;

²Institute of Neurogenetics, University of Lübeck, Germany



jenny.ghelfi

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gpm1w5gzp/v1>

Protocol Citation: Jenny Ghelfi, Soraya Hezzaz, Sylvie Delcambre, Anne Grünewald 2025. Digital PCR Analysis for Mitochondrial DNA Integrity and Copy Number (Absolute Q MAP16 Format). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.dm6gpm1w5gzp/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 29, 2025

Last Modified: November 01, 2025

Protocol Integer ID: 231061

Keywords: dPCR, mitochondrial deletion, mitochondrial transcription, copy number, digital pcr analysis for mitochondrial dna integrity, digital pcr analysis, mitochondrial dna integrity, digital pcr, mitochondrial dna, setup of digital pcr, mtdna

Funders Acknowledgements:

Michael J. Fox Foundation

Grant ID: 15744

Michael J. Fox Foundation

Grant ID: 025802

Michael J. Fox Foundation

Grant ID: 021147

Abstract

This protocol describes the preparation and setup of digital PCR assays to quantify mitochondrial DNA (mtDNA) integrity, deletion load, and copy number using the Absolute Q MAP16 format.

Materials

	A	B	C	D
	Item	Supplier	Catalog / Article No.	Notes
	QuantStudio™ Absolute Q™ MAP16 Plate Kit and Master Mix	Thermo Fisher Scientific	Cat. #A53301	Includes MAP16 plates and 5× Master Mix.
	DLOOP-JUN/QSY Probe	Thermo Fisher Scientific	Custom, Article No. CCU002NR	Target: D-Loop region.
	ND4-ABY/QSY Probe	Thermo Fisher Scientific	Custom, Article No. CCU002NR	Target: ND4 gene.
	B2M-FAM/MGB Probe	Thermo Fisher Scientific	Assay ID APCFGXE, Cat. #43320	Target: B2M (nuclear reference).
	ND1-VIC/MGB Probe	Thermo Fisher Scientific	Assay ID APAANDG, Cat. #44485	Target: ND1 gene.
	Water, nuclease-free	Thermo Fisher Scientific	Cat. #R0582	For all dilutions and reactions.
	TE Buffer, 1× (100 mL)	Promega	Cat. #V6231	Used for probe dilution.
	Ethanol (70%)	VWR International (Avantor)	930.031.006	For surface cleaning.

	A	B	C	D
	Equipment	Brand	Model / SKU	Link / Notes
	QuantStudio™ Absolute Q™ Digital PCR System	Thermo Fisher Scientific	—	Used for MAP16 plate runs.
	Mini Microcentrifuge	Corning	Model 6770	Corning LSE Mini Microcentrifuge
	Vortex Mixer	Thermo Scientific	Cat. #88882011	Basic Vortex Mixer



Protocol materials

70% ethanol **VWR International (Avantor)**

Absolute Q™ DNA Digital PCR Master Mix (5X) **Thermo Fisher Scientific Catalog #A52490**

Water, nuclease-free **Thermo Fisher Catalog #R0582**

Safety warnings

- Wear appropriate personal protective equipment (PPE): laboratory coat, gloves, and safety glasses at all times.
- Clean all work surfaces and equipment with 70% ethanol before and after use.
- Dispose of biological waste, tips, and consumables according to your institution's biosafety and chemical waste disposal procedures.
- Follow the manufacturer's safety data sheets (SDS) for all reagents and probes.

Before start

Samples should be extracted using a clean, preferably column-based method to minimize imaging artifacts during digital PCR. It is recommended to test a small range of concentrations before scaling up to a full plate.

Probe Preparation

10m

1 Check the format of probes:

- Use 20× probes directly
- If supplied as 60×, dilute to 20× in 1× TE buffer (20 µL probe + 40 µL TE)

Quencher guideline:

A maximum of two MGB-quenched assays should be combined per reaction to ensure optimal fluorescence separation. QSY quenchers are compatible without limitation.

Digital PCR Setup (Absolute Q MAP16 Format)

40m

2 Plate Preparation

2.1 DNA template: Use DNA from CSF exosomes, PBMCs, or whole blood

- Avoid impure extracts to prevent imaging artifacts

2.2 MAP16 plate capacity:

Up to 16 samples

Unused columns may be sealed and stored for future runs

2.3 Prepare the master mix for each reaction:

- 2 µL

Absolute Q™ DNA Digital PCR Master Mix (5X) Thermo Fisher Scientific Catalog #A52490

, vortex well before pipetting

- 0.5 µL of each 20× probe mix (or 0.166 µL if using 60× probes)
- Add Water, nuclease-free Thermo Fisher Catalog #R0582 , variable volume (depending on DNA input later, final volume of master mix and DNA = 10 µL)
- Prepare one extra reaction to account for pipetting loss
- Vortex

2.4 Combine DNA and master mix:

- For CSF exosomes: 5 µL DNA + 5 µL master mix
- For PBMC or whole blood DNA: 1.1 µL µL DNA + 8.9 µL master mix

2.5 Vortex and spin down to collect contents

2.6 Load the MAP16 plate:

- Add  9 μ L sample (DNA and master mix) to each well
- Add  15 μ L isolation buffer
- Pipette at a 45° angle to avoid damaging the bottom of wells
- Apply gasket strips after each loaded column
- Seal unused columns

Running the Plate on the Absolute Q

1h 20m

3 Run preparation

3.1 Thermal Cycling Parameters

- Pre-heat: 96 °C for 10 min
- 40 cycles of:
 - 96 °C for 5 s
 - 60 °C for 30 s

3.2 Instrument Preparation and Run Start

- Clean the instrument surface with  70% ethanol **VWR International (Avantor)** and allow to dry completely.
- Load the prepared MAP16 plate into the Absolute Q instrument
- Verify plate orientation matches the software layout
- Start the run

Expected Results and Data Analysis

20m

4 Expected Results and Analysis

4.1 Expected Results

- Clean partition images
- Clean separation of positive and negative partitions
- Multiplex fluorophores:
 - VIC (ND1), FAM (B2M), ABY (ND4), JUN (D-loop).

4.2 Analysis

Quantitative results should be analyzed using the Absolute Q software or equivalent digital PCR analysis platform.

Calculate the following ratios for interpretation of mitochondrial DNA status:

20m



- ND4 / ND1 → indicator of *mtDNA deletion load*
- D-Loop / ND1 → indicator of *mtDNA integrity and replication*
- ND1 / B2M → indicator of *mtDNA copy number per cell*