

Jun 03, 2020

O Digoxigenin Labeled In Vitro Transcription and Ribogreen Quantification

DOI

dx.doi.org/10.17504/protocols.io.bg4vjyw6

Allen Institute for Brain Science¹

¹Allen Institute

BICCN / BICAN

Allen Institute for Brain S...



Allen Institute

Allen Institute

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account





DOI: https://dx.doi.org/10.17504/protocols.io.bg4vjyw6

Protocol Citation: Allen Institute for Brain Science 2020. Digoxigenin Labeled In Vitro Transcription and Ribogreen Quantification. **protocols.io** https://dx.doi.org/10.17504/protocols.io.bg4vjyw6

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: June 03, 2020

Last Modified: March 17, 2021

Protocol Integer ID: 37749

Keywords: MB0002, IVT, in vitro transcription, DIG-labeled RNA, rna from polymerase chain reaction, digoxigenin, vitro transcription, rna, ribogreen quantification this protocol, ribogreen quantification, transcription, using ribogreen, polymerase chain reaction, pcr product, generating dig,

Abstract

This protocol is used for in vitro transcription (IVT) reaction generating DIG-labeled anti-sense (AS) RNA from polymerase chain reaction (PCR) product for in situ hybridization (ISH), as well as fluorescent quantification of products using Ribogreen immediately following. This protocol is written for 96 reactions but can also be used for single reactions or in multiples of 96 by adjusting the amount of materials and equipment appropriately.

Note: Research reported in this publication was supported by the National Institute Of Mental Health of the National Institutes of Health under Award Number U19MH114830. The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.

Attachments



MB0002_Digoxigenin_L.

<u>--</u> 47KB

Troubleshooting

