

Feb 27, 2025

DIO-SPOTlight Image Processing and Analysis v.250227

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov1dr3kvr2/v1>

Matt Oliver^{1,2}

¹Calakos Lab; ²Duke University



Matt Oliver

Duke University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov1dr3kvr2/v1>

Protocol Citation: Matt Oliver 2025. DIO-SPOTlight Image Processing and Analysis v.250227. protocols.io
<https://dx.doi.org/10.17504/protocols.io.ewov1dr3kvr2/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 27, 2025



Last Modified: February 27, 2025

Protocol Integer ID: 123545

Keywords: ASAPCRN, DIO-SPOTlight image processing/analysis, spotlight immunohistochemistry, spotlight image processing, image acquisition protocol, spotlight, image processing, abstract this protocol, dio, processing, protocol, image, partner protocol, imaging

Funders Acknowledgements:

ASAP

Grant ID: ASAP-020607

Abstract

Abstract

This protocol contains the image processing and analysis steps that follow the image acquisition protocol.

Partner protocols: [DIO-SPOTlight Immunohistochemistry v.250227](#) & [DIO-SPOTlight Confocal Imaging Acquisition v.250227](#)

1 Image Analysis

- 2 Open Fiji (ImageJ) perform batch conversion from IMS to TIFF format
 - a. **Process > Batch**
 - i. Batch convert window opens. Select input and output files. Output format: TIFF; Interpolation: none; Scale factor: 1.00. Press convert.
- 3 Background subtraction and ROI cellular masking.
 - a. Open image and duplicate
 - b. Select duplicated image and perform a 15-pixel Gaussian blur.
 - i. **Process > Filters > Gaussian Blur**
 1. Gaussian Blur window opens. Type 15 in the "Sigma (Radius)" box. Select OK.
 - c. Subtract blurred image from the original.
 - i. **Process > Image Calculator**
 1. Image Calculator window opens. For Image 1, select the original. For Image 2, select the blurred duplicate. For Operation, select subtract. Check box "Create in new window." Select OK.
 - d. Select the new background subtracted image and perform a 2-pixel Gaussian blur.
 - i. **Process > Filters > Gaussian Blur**
 1. Gaussian Blur window opens. Type 2 in the "Sigma (Radius)" box. Select OK and save the image for mask making.
 - e. Locate the desired plane containing a high number of the cells of interest and duplicate that plane for the original and background subtracted image.
 - f. Select the background subtracted image and create binary masks of the anti-ChAT and nuclear stained channels (channels 1 and 4):
 - i. **Image > Adjust > Threshold**
- 4 ROI selection of binary images
 - a. Select the wand tool and select a binary cholinergic neuron and add to the ROI manager. For the same cell switch to the binary nuclear stained channel and add the nucleus of the cell to the ROI manager.
 - b. Segregate the cytoplasmic and nuclear signal using the XOR function in the ROI manager.
 - i. **ROI manager > More > XOR**
 - c. Repeat for the other cells of interest
 - d. Save the coordinates.
- 5 Select the original raw image corresponding to the masked cells. Measure the fluorescent signal of the selected ROI and save or directly transfer to a spread sheet. Label as ChAT.

- 6 Create a binary mask of neighboring non-cholinergic neurons.
 - a. Select the RFP and GFP channels (channels 2 and 3) from the background subtracted image.
 - b. Create an additive image of both RFP and GFP channels to limit biasing ROI selection toward cells expressing low levels of either fluorescent protein.
 - i. **Process > Image Calculator**
 1. Image Calculator window opens. For Image 1, select RFP channel. For Image 2, select GFP channel. For Operation, select Add. Check box "Create in new window." Select OK.
 - ii. Create a threshold binary mask for the new additive image the same as above for the ChAT and nuclear stained channels.
- 7 ROI selection of binary images for neighboring neuron selection.
 - a. Add ROI to ROI manager and exclude nucleus the same as described above.
 - b. Include background measurements to ROI manager from neurons where Cre recombination failed.
 - c. Save neighboring neuron ROI coordinates.
- 8 If cells did not include a nucleus in the analyzed plan no XOR function was performed. In some instances, manual ROIs were drawn.
- 9 Select the original raw image corresponding to the masked cells. Measure the fluorescent signal of the selected ROI and save or directly transfer to a spread sheet. Label as neighboring neurons
- 10 Cytoplasmic measurements only were used for analysis.