

Mar 10, 2025

DIO-SPOTlight Live 2-Photon Imaging v.250227

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l6dj35vqe/v1>

Matt Oliver^{1,2}, Victoria Hall^{1,2}

¹Calakos Lab; ²Duke University



Matt Oliver

Duke University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l6dj35vqe/v1>

Protocol Citation: Matt Oliver, Victoria Hall 2025. DIO-SPOTlight Live 2-Photon Imaging v.250227 . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bp2l6dj35vqe/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 27, 2025

Last Modified: March 10, 2025

Protocol Integer ID: 123548

Keywords: ASAPCRN, DIO-SPOTlight 2-photon, photon imaging, live tissue preparation, spotlight, tissue, imaging

Funders Acknowledgements:

ASAP

Grant ID: ASAP-020607

Abstract

This protocol contains the steps taken for live tissue preparation and two-photon imaging.

Prepare acute brain slices

- 1 Anesthetize mouse and perform intracardial perfusion with ice cold N-methyl-D-glucamine (NMDG)-based artificial cerebrospinal fluid (ACSF).
 - a. ACSF (90 mM NMDG, 2.5 mM KCl, 1.2 mM NaH₂PO₄, 35 mM NaHCO₃, 20 mM HEPES, 25 mM glucose, 5 mM sodium ascorbate, 2 mM thiourea, 3 mM sodium pyruvate, 10 mM MgSO₄ heptahydrate, 0.5 mM CaCl₂ dihydrate, and saturated with 95% O₂ and 5% CO₂)
- 2 Quickly remove the brain and slice into 250 µm sagittal slices using a vibratome.
- 3 Immediately transfer slices to a recovery chamber containing NMDG ACSF for 10 to 15 min at 32°C.
- 4 Transfer heat recovered slices to a holding chamber containing ACSF and let rest for approximately 1 hour
 - a. ACSF (124 mM NaCl, 2.5 mM KCl, 1 mM MgCl₂, 2 mM CaCl₂, 26 mM NaHCO₃, 1.2 mM NaH₂PO₄, and 10 mM glucose, saturated with 95% O₂ and 5% CO₂, pH 7.4, 300 mOsm/l)

Two-photon imaging

- 5 Place the slices into a perfusion chamber and superfuse with oxygenated ACSF.
- 6 Image slices using a Bruker Investigator Multiphoton Imaging System or equivalent.
- 7 Set the emitted light wave-length to 920 nm on the Coherent Chameleon I Ti:Sapphire femtosecond laser.
 - a. Keep laser power between 35 and 50 mW
- 8 Set up image stack using galvo-galvo scanning with frame rate of 1 Hz with 3-4 µm depth increments at either 1X or 4X optical zoom.
- 9 Collect fluorescence signal using a 16X W objective and Quad GaAsP photomultiplier tubes (PMTs).
- 10 Split the SPOTlight red and green channels using a 565lpxr dichroic mirror followed by a 525/70 bandpass PMT filter.

