

May 30, 2025

DIO-SPOTlight mice tissue processing

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gpqno1lzp/v1>

Matt Oliver¹

¹Duke University



Matt Oliver

Duke University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gpqno1lzp/v1>

Protocol Citation: Matt Oliver 2025. DIO-SPOTlight mice tissue processing. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.dm6gpqno1lzp/v1>

Manuscript citation:

Matthew L Oliver, Zachary F Caffall, Callie B Eatman, Timothy D Faw, Nicole Calakos
(2025) **DIO-SPOTlight Transgenic Mouse to Functionally Monitor Protein Synthesis Regulated by the Integrated Stress Response**
*eLife*14:RP104457 <https://doi.org/10.7554/eLife.104457.1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 30, 2025

Last Modified: May 30, 2025

Protocol Integer ID: 219161

Keywords: ASAPCRN, Perfusion , Tissue processing, spotlight mice tissue, perfusing mice, mice tissue, tissue slicing, tissue, cryostat, tissue processing, mice

Funders Acknowledgements:

ASAP

Grant ID: ASAP-020607

Abstract

This protocol contains a relatively standard description of transcardially perfusing mice followed by tissue slicing using a cryostat.

DIO-SPOTlight mice tissue processing

- 1 Anesthetize mice using the drop jar method containing 3-5% aerosolized isoflurane.
- 2 Perform a toe pinch to ensure the mouse is fully anesthetized prior to pinning the limbs down.
- 3 Use scissors to open the thoracic cavity and perform a transcardial perfusion by inserting a needle in to the right ventricle and cut the top of the left lobe, immediately perfuse using ~25 ml ice cold PBS followed by ~30 ml ice cold 4% paraformaldehyde (PFA).
- 4 Cut off the head and remove the skin behind the ears to behind the eyes.
- 5 Open the skull and remove the brain.
- 6 Place the brain in 4% PFA + 20% sucrose in 1x PBS for 24 hours at 4C.
- 7 The following day transfer the brains to 30% sucrose in 1x PBS for an additional 24 hours at 4C.
- 8 Embed the tissue in optimal cutting temperature compound (OCT) and freeze.
- 9 The OCT embedded brains are then sliced cryostat into 40 um sagittal sections and store in 1x PBS with 0.01% sodium azide at 4C until IHC is performed.
- 10 If spinal cord dissections are performed following the removal of the brain use the following landmarks to dissect the spine: cervical, base of skull to caudal T1 lamina; thoracic, 0.6 cm segment centered on T9 lamina; and lumbar, identified based on L1-L6 nerve roots.
- 10.1 Repeat steps 6 - 9 with the exception that the spinal cord sections were directly mounted onto microscope slides.

Acknowledgements

Spinal cord dissections performed by Timothy Faw.