

Jun 25, 2025

Version 3

Direct Detect D-HMQC V.3



In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.kqdg32qo7v25/v3>

Alexander L. Paterson¹

¹National Magnetic Resonance Facility at Madison (NMRFAM), University of Wisconsin-Madison, Madison, WI, United States



NMRFAM Facility

UW-Madison

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg32qo7v25/v3>

Protocol Citation: Alexander L. Paterson 2025. Direct Detect D-HMQC. [protocols.io](https://doi.org/10.17504/protocols.io.kqdg32qo7v25/v3)

<https://dx.doi.org/10.17504/protocols.io.kqdg32qo7v25/v3> Version created by [NMRFAM Facility](#)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 25, 2025

Last Modified: June 25, 2025

Protocol Integer ID: 221049

Keywords: Materials Quadrupolar Direct-Detect D-HMQC, space correlations between quadrupolar nuclei, quadrupolar nuclei, nuclei pair, robust choice for nuclei, quantum coherence, spin nuclei, hmqc purpose, examples of nuclei pair, hmqc benefits from alternate implementation, nuclei, hmqc, hmqc benefit, dipolar heteronuclear

Funders Acknowledgements:

National Science Foundation

Grant ID: 1946970

Abstract

Purpose

To observe through-space correlations between quadrupolar nuclei and spin-1/2 nuclei.

Scope

The dipolar heteronuclear multiple-quantum coherence (D-HMQC) experiment reveals through-space proximities between two nuclei, facilitating spectral assignments. This implementation is a robust choice for nuclei that are significantly averaged at moderate MAS rates, e.g., 20 kHz. Examples of nuclei pairs which benefit from this experiment are $^{27}\text{Al}\{^1\text{H}\}$, $^{23}\text{Na}\{^{31}\text{P}\}$, and $^{11}\text{B}\{^{29}\text{Si}\}$. The indirect-detect, or proton-detect, version of D-HMQC benefits from alternate implementations.

Guidelines

As this implementation of D-HMQC places the dipolar recoupling pulses on the indirect channel, the quadrupolar nucleus should be the one that is detected.

This pulse sequence should not be used for observing the correlations between two quadrupolar nuclei.

This protocol uses the SR4_1^2 recoupling scheme.



Materials

Definitions:

1. D-HMQC: Dipolar Heteronuclear Multiple-Quantum Coherence
2. MAS: Magic angle spinning
3. ν_r : Magic angle spinning speed
4. ν_{rf} : Radiofrequency power
5. X: Direct nucleus
6. Y: Indirect nucleus
7. CT: Central transition

Appendix:

While the $SR4_1^2$ recoupling sequence is robust against minor rf imperfections and MAS instability, in dilute systems this sequence can suffer from t_1 noise. Reducing the **d1** to a very short value (e.g., $0.2 \times T_1$) can mitigate this noise. If doing so, ensure that the overall probe duty cycle is respected.

Before start

Operator should be familiar with acquisition of 2D NMR spectra.

Operator must be familiar with both the MAS spinning limitations of the probe and the RF power limitations of the probe, particularly on the indirect channel.

Operator must be familiar with the duty cycle of the probe.

Expected amount of time SOP will use: 24 hours.



Procedure

1 Identify an appropriate MAS rate, ν_r , for the sample under investigation.

1.1 Note that the Y rf power, ν_{rf} , will need to be set such that $\nu_{rf} = 2\nu_r$.



Safety information

Ensure that this is within the safe limits of the probe.

Note

A MAS speed which returns a rotor period which is an integer number of microseconds can aid in rotor synchronization. E.g., $\nu_r = 13333.3$ Hz, 25000 Hz.

2 Optimize a CT-selective 90° pulse length and power for the X nucleus.

Note

Overly long pulse lengths can cause issues with rotor synchronization.

3 Optimize a Y nucleus 90° pulse width and power that $\nu_{rf} = 2\nu_r$.

Note

For example, an MAS rate of 25 kHz would require $\nu_{rf} = 50$ kHz, *i.e.*, a 5 μ s pulse.

4 Load the **dp_d-hmqc_indirect_rec.nmrfam** pulse sequence and set the following parameters:

4.1 **cnst31**: MAS rate in Hz



4.2 **p3**: X channel CT-selective 90° pulse length

4.3 **plw21**: X channel CT-selective 90° pulse power

4.4 **p1**: Y channel 90° pulse length.

Note

This does not have to be the same as in Step 3.

4.5 **plw3**: Y channel 90° pulse power.

Note

This does not have to be the same as in Step 3.

4.6 **plw23**: Y channel 90° pulse power for SR4_1^2 recoupling, optimized in Step 3.

Note

This must be the same as in Step 3.

4.7 **p16**: Initial SR4_1^2 recoupling time, in μs . 700 μs is a reasonable starting point.

4.8 **d1**: $1.3 \times T_1$.

5 Acquire an initial 1D spectrum. Repeat with increased **ns** until the signal-to-noise is high enough to allow for optimization.

6 Optimize the value of **plw23** to maximize the signal intensity. This will likely be within 1 or 2 W of the previously optimized value.

7 Optimize the value of **p16** to maximize the signal intensity.



**Safety information**

The total time used for SR4₁² recoupling is reported in **p17** in μ s. Ensure that this time is short enough to respect the probe duty cycle and maximum sustained high-power pulse limit.

8 Copy the parameters to a new dataset. Convert this dataset to a 2D experiment.

9 Set the following 2D parameters:

9.1 **FnMODE**: States-TPPI

9.2 **F1 TD**: 64, or as appropriate to allow for adequate F1 resolution.

9.3 **F1 SWH**: MAS rate in Hz.

Note

If this spectral window is too narrow, **F1 SWH** can be set to an integer multiple of the MAS rate in Hz.

10 Acquire the 2D spectrum.



Protocol references

Tricot, G.; Trébosc, J.; Pourpoint, F.; Gauvin, R.; Delevoye, L. The D-HMQC MAS-NMR Technique. Annual Reports on NMR Spectroscopy, 2014; pp 145-184. DOI: [10.1016/B978-0-12-800185-1.00004-8](https://doi.org/10.1016/B978-0-12-800185-1.00004-8)

Example Data at 1.1 GHz: usnan.org/ark:/83454/e1df5b48ac-c838-42c2-9a2c-3e7fc6e5bae8

Example Data at 900 MHz: usnan.org/ark:/83454/e178eceb8c-aae4-407d-982f-27f95ca1b4e7

Example Data at 600 MHz: usnan.org/ark:/83454/e1141a31fc-13f2-4a57-9976-8c718eb4e37c

Protocol

NAME

Quadrupole CT-Selective Pulse Width Optimization

CREATED BY
NMRFAM Facility

Preview

Protocol

NAME

Quadrupole Non-Selective Pulse Width Optimization

CREATED BY
NMRFAM Facility

Preview



Protocol

NAME

Saturation Recovery for Quadrupolar Nuclei

CREATED BY

NMRFAM Facility

Preview