

Aug 19, 2020

## Direct ELISA for investigating the binding of Protein-G to immunoglobulins.

DOI

[dx.doi.org/10.17504/protocols.io.bjxrkpm6](https://dx.doi.org/10.17504/protocols.io.bjxrkpm6)

Angel A Justiz-Vaillant<sup>1</sup>, Monica F. Smikle<sup>2</sup>

<sup>1</sup>University of the West Indies St. Augustine; <sup>2</sup>University of the West Indies. Mona Campus

University of the West In...

angel.vaillant@sta.uwi.e...



Angel A Justiz-Vaillant

University of the West Indies St. Augustine

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bjxrkpm6>

**Protocol Citation:** Angel A Justiz-Vaillant, Monica F. Smikle 2020. Direct ELISA for investigating the binding of Protein-G to immunoglobulins.. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bjxrkpm6>

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 19, 2020

**Last Modified:** August 19, 2020

**Protocol Integer ID:** 40657

**Keywords:** streptococcal protein, mammalian immunoglobulin, binding bacterial protein, immunoglobulin, bacterial protein, fc region of many mammalian immunoglobulin, many mammalian immunoglobulin, linked immunosorbent assay, reactivity with immunoglobulin, binding of protein, immunosorbent assay, binding protein, protein, direct elisa

## Abstract

Streptococcal protein G is an immunoglobulin-binding protein that interacts with the Fc region of many mammalian immunoglobulins [1].

## References

1. Vaillant AJ, McFarlane-Anderson N, Wisdom B, Mohammed W, Vuma S, et al. (2013) Immunoglobulin-binding Bacterial Proteins (IBP) Conjugates and their Reactivity with Immunoglobulin in Enzyme-Linked Immunosorbent Assays (ELISA). J Anal Bioanal Tech 4: 175. doi:10.4172/2155-9872.1000175

## Materials

### MATERIALS



Nunc®; 96-Well Polystyrene Round Bottom Microwell Plates, V 96 well plate, Non-Treated, clear, without lid, Sterile **Thermo Fisher Catalog #260210**



Streptococcal protein G by Sigma Aldrich

## Troubleshooting

- 1 This ELISA is used to study the interaction of Streptococcal protein-G (SpG) with diverse immunoglobulins.
- 2 The 96 well microtitre plate is coated overnight at 4°C with 1 µg/µl per well of purified immunoglobulins or 50 µl of any animal sera in carbonate-bicarbonate buffer pH 9.6.
- 3 Then plate is treated with bovine serum albumin solution and washed 4X with PBS-Tween.
- 4 Then 50 µl of peroxidase-labeled-protein-G conjugate diluted 1:5000 in PBS-non-fat milk is added to each well and incubated for 1.30h at RT. After that the plate is washed 4X with PBS-Tween.
- 5 Pipette 50 µl of 3,3',5,5' - tetramethylbenzidine (TMB; Sigma-Aldrich) to each well.
- 6 The reaction is stopped with 50 µl of 3M H<sub>2</sub>SO<sub>4</sub> solution.
- 7 The plate is visually assessed for the development of colour and read in a microplate reader at 450 nm.
- 8 A cut-off point should be calculated as the mean of the optical density of negative controls x 2 SD.